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
213221 Anti-Jk^a for card method monoclonal, human IgM 5 ml

213223 Anti-Jk^b for card method monoclonal, human IgM 5 ml

730-15-5306 Version 006 / 20.05.2021

CE 0483

Antitoxin GmbH - Industriestraße 88 - 69245 Bammalatt

FOR IN VITRO DIAGNOSTIC USE ONLY

INTENDED USE

Monoclonal agglutinating Anti-Jk^a and -Jk^b reagents are produced from cell culture supernatants of heterohybridoma-cell lines. The cells are secreting an antibody of IgM-type that reacts specific with the corresponding antigen. The antibody is human protein. The reagents are used to in-vitro determine qualitative whether human red blood cells possess or lack the corresponding blood group antigens Jk^a and Jk^b. The reagents are intended to be used by qualified and technical personnel only.

PRINCIPLE OF PROCEDURE

The test method used with these reagents is based on the principle of column agglutination. Normal human erythrocytes, possessing one of these antigens, will agglutinate in the presence of the specific antibody directed toward the antigen.

REAGENTS

The listed reagents contain antibodies of the following cell clones:

- Anti-Jk^a for card method monoclonal, human IgM clone: MS-15
- Anti-Jk^b for card method monoclonal, human IgM clone: MS-8

The reagents contain <0.1% (w/v) sodium azide as preservative. Additionally the reagents are prepared of active antibody, sodium chloride, macromolecules and bovine albumin, which has been tested and certified by the US Veterinary service inspectors.

WARNING

These reagents are prepared from supernatants of cell cultures. As biological products it should be looked upon as potentially infectious because of never complete exclusion of danger through excipients of disease. The reagents contain sodium azide that may be toxic and may react with lead or copper to form highly explosive salts.

On disposal, flush with large quantities of water.

For the reasons mentioned above the reagents should be handled with proper care.

STORAGE REQUIREMENT

Store opened and unopened products at 2 to 8°C. May be at room temperature while in use.

In principle, store and use the reagents only to declared expiry date.

REMARKS

1. With each testing positive and negative controls should be performed.
2. Inappropriate storage impairs efficacy of the reagents.
3. Weak turbidity of the reagents does not affect its reactivity. Bacteria and chemical contamination of the product should be avoided. If a visible change is detected, the reagent should no longer be used, this sign may indicate a microbiological contamination.
4. Strength of positive reactions also depends on age of used blood.
5. Use the appropriate card centrifuge. The use of another card-specific centrifuge (each card centrifuge has its specified unchangeable g-force) may lead to false results due to the changed g-force.
6. The test method identified below is for manual testing only. When using automated or semi-automated instruments, follow the procedures that are contained in the operator's manual provided by the device manufacturer. Laboratories must follow approved validation procedures to demonstrate compatibility of this product on automated systems.
7. For usage of these reagents all effective national laws, directives and guidelines has to be observed, in Germany especially the „Richtlinien zur Gewinnung von Blut und Blutbestandteilen und zur Anwendung von Blutprodukten (Hämotherapie)“¹ in the current version.
8. The information on the use of the Grifols test cards in the relevant insert must be observed.

SAMPLE PREPARATION

1. Blood samole should be collected by approved medical procedure.
2. Blood samples to be tested should be used as soon as possible after blood collection to reduce the risk of false-positive and false-negative results due to improper storage or contamination of the reagents. If a delay in testing occurs, samples should be stored at 2 to 8°C. Blood drawn into EDTA should be tested within 7 days and samples treated with sodium citrate within 14 days after collection. Blood bag / Donor Blood can be tested until the expiry date.

REAGENT PREPARATION

There is no preparation of the reagents required. Take and use reagent directly from the vials.

PROCEDURE

Not provided material additionally needed

Card Method

1. Cards: Grifols „DG Gel Neutral“ or Grifols „DG Gel Coombs“
2. Microliter pipette
3. Tubes
4. Grifols Card Centrifuge
5. DG Gel Sol (card specific diluent)

Test procedure

Grifols Card Method (Neutral- or Coombs card)

1. Prepare 0,8 % suspension of red blood cells in DG Gel Sol (card specific diluent).
2. Add 50 µL of appropriate cell suspension to each marked micro tube.
3. Add 25 µL of appropriate reagent to each micro tube.
4. Centrifuge the card in appropriate card centrifuge with the for this centrifuge unchangeable g-force.
5. Check macroscopically for agglutination within 30 minutes.
6. Document the results.

INTERPRETATION OF RESULTS

Positive results (+): visible agglutination of erythrocytes is a positive result and indicates the presence of the corresponding antigen.

Negative results (-): No visible agglutination of erythrocytes is a negative result and indicates the absence of the corresponding antigen.

Read and interpret the results of the card method according to the Grifols card instruction for use.

LIMITATIONS OF THE PROCEDURE

1. Inaccuracy at compliance with instructions written under section “Procedures” and “Interpretation of results” may lead to incorrect results.
2. No valid conclusion concerning the test result can be reached, if controls with uncertain or false results occur.
3. Enzyme treated erythrocytes or addition of bovine albumin and/or other solutions containing protein may cause unspecific reactions.
4. Hemolyzed, turbid, contaminated or clotted blood samples must not be used in this test.
5. Due to variability of antigen expression on human red blood cells, reactivity of the reagents, mentioned above, against certain phenotypes may give weaker reactivity compared to control cells.
6. No one specific antiserum or technique can be guaranteed to detect all rare, weak or variant antigens.²
7. Red blood cells coated with alloantibodies or autoantibodies of the same or similar specificity as the reagent used for the test (i.e. cells that are positive in the direct antiglobulin-test (DAT)) are not suitable for this test procedure.
8. Red blood cells with a positive direct Coombs-test may cause false-positive reactions in card method.
9. It is described in the literature, that samples from patients treated with anti-CD38 monoclonal antibodies can cause false positive results in the Coombs-test.⁵
10. Pay attention to all statements to limitations in inserts of Grifols card.
11. Incubation at room temperature may increase potentially weak reactions on „DG Gel Neutral“.
12. For donors, a small pellet of non-agglutinated red blood cells in the bottom of the microtube may be observed in positive reactions (image of double population, DP). In patients, DPs could also be caused by prior transfusion. Additional information on patient history could be necessary for resolution.

PERFORMANCE CHARACTERISTICS

A performance evaluation for the products was conducted.

Different samples material was used and compared with a reference method.

When testing with the „DG Gel Neutral“ card Blood bag with the additive solutions SAG-M and PAGGS-M, as well as EDTA- und Citratblood were used.

When testing with the „DG Gel Coombs“ card Blood bag with the additive solution PAGGS-M, as well as EDTA- und Citratblood were used.

Anti-Jk^a agreement of the results and performance data.

Method	Positive	Negative	Specificity	Specificity
DG Gel Neutral	103/103	43/43	100%	100%
DG Gel Coombs	54/54	20/20	100%	100%

Anti-Jk^b agreement of the results and performance data:

Method	Positive	Negative	Specificity	Specificity
DG Gel Neutral	106/106	36/36	100%	100%
DG Gel Coombs	54/54	20/20	100%	100%

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
PRESENTATION

- 213221 Anti-Jk^a for card method monoclonal, human IgM 5 ml
- 213223 Anti-Jk^b for card method monoclonal, human IgM 5 ml

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Antitoxin GmbH - Industriestraße 88 - 69245 Bammental