NUR ZUR IN-VITRO-DIAGNOSTIK

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FOR IN VITRO DIAGNOSTIC USE ONLY

INTENDED USE

Polyclonal Coombs-reactive Anti-Fy^a and Fy^b reagents are produced from human plasma that contains antibody of IgG-type, which reacts specific with the corresponding blood group antigen. The reagents are used to in-vitro determine qualitative whether red blood cells possess or lack the corresponding blood group antigen Fy^a and Fy^b.

The reagents are intended to be used by qualified and technical personnel only.

PRINCIPLE OF PROCEDURE

The test method used with these reagents is based on the principle of column agglutination. Normal human erythrocytes, possessing the appropriate antigen, will be recognized and coated by the corresponding specific antibody and then the cells will be agglutinated by a secondary antibody that reacts with human IgG-molecules.

REAGENTS

The listed reagents are available in following formulation:

Anti-Fy^a for card method Coombs-reactive, polyclonal, human
Anti-Fy^b for card method Coombs-reactive, polyclonal, human

The reagents contain <0.1% (w/v) sodium azide as preservative. Beside the parts active antibody and human serum, the reagents contain sodium chloride, macromolecules and bovine albumin.

WARNING

The reagents are prepared from human plasma. The raw materials for these reagents are tested for HBsAg, HIV 1/2- and HCV-antibodies and found to be negative. Nevertheless as biological product it should be looked upon as potentially infectious because of never complete exclusion of danger through excitants of disease. The reagents contain sodium azide, that may be toxic and may react with lead or copper to form highly explosive salts. On disposal, flush with large quantities of water. For the reasons mentioned above these reagents should be handled with proper care.

STORAGE REQUIREMENT

Store opened and unopened products at 2 to 8°C. May be at room temperature (15 to 30°C) while in use. In principle, store and use the reagent to declared expiry date only.

REMARKS

- With each testing positive and negative controls should be performed.
- Inappropriate storage impairs efficacy of the reagents.
- Weak turbidity of the reagents does not affect its reactivity.
- Strength of positive reactions also depends on age of used blood.
- The use of another card-specific centrifuge (each card centrifuge has its specified unchangeable g-force) may lead to false results due to the changed g-force.
- Blood samples to be tested should be used as soon as possible. If a delay in testing occurs, samples should be stored at 2 to 8°C. Blood drawn into EDTA should be tested within 7 days and samples treated with sodium citrate within 14 days after collection.
- The test method identified below is for manual testing only. When using automated or semi-automated instruments, follow the procedures that are contained in the operator's manual provided by the device manufacturer. Laboratories must follow approved validation procedures.
- For usage of these reagents all effective national laws, directives and guidelines have to be observed, in Germany especially the „Richtlinien zur Gewinnung von Blut und Blutbestandteilen und zur Anwendung von Blutprodukten Hämotherapie“*1 in its actual form.
- The information on the use of the Grifols test cards in the relevant insert must be observed.

REAGENT PREPARATION

There is no preparation of the reagents required. Take and use the reagents directly from the vials.

PROCEDURE

Not included but necessary materials:

Card Method

- Cards: Grifols „DG Gel Coombs“
- Microtiter pipette
- Tubes
- Grifols Card centrifuge
- Grifols Card Incubator
- DG Gel Sol (card specific diluent)

Test procedure**Card Method**

- Prepare 0,8 % suspensions of red blood cells in DG Gel Sol (card specific diluent).
- Add 50 µL of appropriate cell suspension to each marked micro tube.
- Add 25 µL of the appropriate reagent to each micro tube.
- Incubate card at 37°C for 15 minutes in a Grifols card incubator.
- Centrifuge the card in a Grifols card centrifuge with the for this centrifuge unchangeable g-force.
- Check macroscopically for agglutination within 30 minutes.
- Document the results.

INTERPRETATION OF RESULTS

Positive results (+): visible agglutination of erythrocytes is a positive result and indicates the presence of the corresponding antigen.

Negative results (-): No visible agglutination of erythrocytes is a negative result and indicates the absence of the corresponding antigen.

Read and interpret the results of the card method according to the Grifols card instruction for use.

LIMITATIONS OF THE PROCEDURE

- Inaccuracy at compliance with instructions written under section "Procedures" and "Interpretation of results" may lead to incorrect results.
- No valid conclusion concerning the test result can be reached, if controls with uncertain or false results occur.
- Enzyme treated erythrocytes or addition of bovine albumin and/or other solutions containing protein may cause unspecific reactions.
- Hemolyzed, turbid, contaminated or clotted blood samples must not be used in this test.
- Due to variability of antigen expression on human red blood cells, reactivity of the reagents, mentioned above, against certain phenotypes may give weaker reactivity compared to control cells.
- Red blood cells coated with alloantibodies or autoantibodies of the same or similar specificity as that reagent used for the test (i.e. cells that are positive in the direct antiglobulin-test (DAT)) are not suitable for this test procedure.
- Red blood cells with a positive direct Coombs-test may cause false-positive reactions.
- Pay attention to all statements to limitations in inserts of Grifols card.

PERFORMANCE CHARACTERISTICS

A performance evaluation was conducted. Different samples (donor, patient, panel blood) was used and compared with other products.

Product	Positive	Negative	False Positive	False Negative
Anti-Fy ^a	129	57	1	0
Anti-Fy ^b	142	43	2	0

The calculated values found for:

Product	Sensitivity	Specificity
Anti-Fy ^a	100%	98,28%
Anti-Fy ^b	100%	95,56%

LITERATURE

- Peter D. Issit, David J. Anstee Applied Blood Group Serology, 4. edition, Montgomery Scientific Publications 1998
- Christian Mueller-Eckhardt, Volker Kiefel Transfusionsmedizin, 3. Auflage, Springer-Verlag 2004.

PRESSENTATION

- 213187 Anti-Fy^a for card method Coombs-reactive, polyclonal, human 5ml
213188 Anti-Fy^b for card method Coombs-reactive, polyclonal, human 5ml

730-15-5105 Version 005 / 17. June 2019

	In vitro diagnostics device
	Batch Code
	Use by
	Temperature limitation
	Consult Instructions for use
	Catalogue Number
	Cards
	Manufacturer

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Antitoxin GmbH - Industriestraße 88 - 69245 Bammental