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INTENDED USE

Mollison and Polley discovered in 1964 that decrease in reducing ionic strength through a low ionic strength solution (LISS) enhances the antigen-antibody response. The present Liss solution is a modification of the low-salt medium for antibody diagnostics described by Löw, B. and Messeter, L. 1974. LISS solution modified is used in blood group serology to increase sensitivity / enhancement of agglutinations, as well as to shorten the reaction time between antibodies and erythrocytes.

**Excluded** is the use in connection with serological blood typing of:

- ABO system [A (ABO1), B (ABO2), AB (ABO3)];
- Rhesus system [RH1 (D), RHW1, RH2 (C), RH3 (E), RH4 (c), RH5 (e)];
- Kell system [Kel1 (K)];
- Kidd system [JK1 (Jka), JK2 (Jkb)];
- Duffy system [FY1 (Fya), FY2 (Fyb)].

The reagent is intended to be used by qualified and technical personnel only.

PRINCIPLE OF PROCEDURE

The test method use with this solution is based on the principle of agglutination. The low ionic strength solution in combination with a small proportion of bovine albumin solution and the lower ion concentration than in a physiological NaCl solution causes fluid inflow into the erythrocytes, thereby forming a larger surface area, making various antigens more accessible to the corresponding antibodies (increasing sensitivity), increasing the amount of AK outside the cell, however, the specificity is maintained. Normal human erythrocytes carrying the corresponding antigen are recognized and agglutinated by a corresponding antibody.

REAGENTS

The listed reagent is offered in the following form:

LISS solution *modified*

The reagent contains Thiomerosal as preservative. The solution also contains sodium chloride, macromolecules and bovine albumin, which has been tested and certified by the US Veterinary service inspectors.

WARNING

LISS solution is prepared from biological material. As biological product it should be looked upon as potentially infectious because of never complete exclusion of danger through excitants of disease. The solution contains thiomerosal that may be toxic. On disposal, flush with large quantities of water. For the reasons mentioned above the reagent should be handled with proper care.

STORAGE REQUIREMENT

Store opened and unopened products at +2 to +8°C. May be at room temperature while in use. In principle, store and use the reagent to declared expiry date only.

REMARKS

1. With each testing positive and negative controls should be performed.
2. Inappropriate storage impairs efficacy of the reagent
3. Weak turbidity of the reagent does not affect its reactivity. Bacteria and chemical contamination of the reagent should be avoided. If a visible change is detected, the reagent should no longer be used, this sign may indicate a microbiological contamination.
4. Strength of positive reactions also depends on age of used blood.
5. Centrifugation highly different from appointed relative centrifugal force may lead to false results.
6. The test method identified below is for manual testing only. When using automated or semi-automated instruments, follow the procedures that are contained in the operator's manual provided by the device manufacturer. Laboratories must follow approved validation procedures.
7. For usage of this reagent all effective national laws, directives and guidelines have to be observed in its current valid version, in Germany especially the „Richtlinien zur Gewinnung von Blut und Blutbestandteilen und zur Anwendung von Blutprodukten (Hämotherapie)“<sup>1</sup>.

SAMPLE PREPARATION

1. Blood sample should be collected by approved medical procedure.
2. Blood samples to be tested should be used as soon as possible after blood collection to reduce the risk of false-positive and false-negative results due to improper storage or contamination of the reagent. If a delay in testing occurs, samples should be stored at +2 to +8°C. Blood drawn into EDTA should be tested within 7 days and samples treated with sodium citrate within 14 days after collection. Blood bag / Donor Blood can be tested until the expiry date.

REAGENT PREPARATION

There is no preparation of the reagent required. Take and use the reagent directly from the vials.

PROCEDURE

Not included but necessary materials:

- Tube Centrifugation Method:
1. tubes (10 x 75 mm or 12 x 75 mm)
  2. microliter pipette
  3. timer
  4. incubator
  5. centrifuge
  6. isotonic saline (0,85 - 0,9% sodium chloride)
  7. Anti-Human-Globulin-Reagent (Coombs-Serum / AHG- Serum)

Test procedure

Tube Centrifugation Method

1. Prepare 2% to 5% suspensions of red blood cells in isotonic saline (red blood cells may be washed 1-3 times with isotonic saline).
2. In a labeled test tube add 50µl (alternative one drop = approximately 50µl) of appropriate cell suspension
3. Add 100µl (alternative two drops = approximately 100µl) serum/plasma of a patient/proband or bloodgrouping reagent.
4. Subsequently add 100µl (alternative two drops = approximately 100µl) LISS solution to the tube.
5. Mix Erythrocytes / Reagentmixture well by slightly shaking.
6. Incubate tube in an incubator at +37 °C for 10 min.
7. Centrifugation of the tube for 1 minute at 1.000 rpm (approximately 180-270 x g).
8. Gently shake the red cells completely from the bottom of the tube and check macroscopically for agglutination within 3 minutes.
9. Document the result.
10. Wash red cells 3 times with (cold) isotonic saline.
11. Subsequently add 100 µL Anti-Human-Globulin-Reagent (Coombs-Serum / AHG-Serum) to the tube, release the cells from the bottom of the tube by slightly shaking and mix with the Coombs-Serum / AHG-Serum.
12. Centrifuge tube for 1 min at 1.000 rpm (approximately 180-270 x g).
13. Gently shake the red cells from the bottom of the tube and check macroscopically for agglutination within 3 minutes.
14. Document the result

INTERPRETATION OF RESULTS

"Slightly shaking" at Tube Centrifugation Method:

- Positive results (+): visible agglutination of erythrocytes is a positive result and indicates the presence of the corresponding antigen.
- Negative results (-): No visible agglutination of erythrocytes is a negative result and indicates the absence of the corresponding antigen.

LIMITATIONS OF THE PROCEDURE

1. Inaccuracy at compliance with instructions written under section "Procedures" and "Interpretation of results" may lead to incorrect results.
2. No valid conclusion concerning the test result can be reached, if controls with uncertain or false results occur.
3. Hemolyzed, turbid, contaminated or clotted blood samples must not be used in this test.
4. No one specific antiserum or technique can be guaranteed to detect all rare, weak or variant of antigens.<sup>2</sup>
5. Due to variability of antigen expression, reactivity of the reagent against certain phenotypes may give weaker reactivity compared to control cells.
6. Red blood cells coated with alloantibodies or autoantibodies (i.e., cells that are positive in the direct antiglobulin-test (DAT)) are unsuitable for testing.
7. As described in the literature, samples from patients treated with anti-CD38 monoclonal antibodies can cause false positive results in the Coombs test.<sup>5</sup>

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 REF	Product Code	 LOT	Lot
	Store from - to		Expiration Date
 IVD	In-Vitro Diagnostic		EU CE symbol
	Manufacturer as to 98/79/EG		Observe instruction for use
 UDI	Unique Device Identification	Distributor	

REF

BG1651 LISS solution 10 ml

730-13-2008 Version 008 / 25.02.2025

CE

 ANTITOXIN GmbH Industriestraße 88 69245 Bammental Germany

 +49 (0) 6223/ 8661-0  +49 (0) 6223/ 8661-13  [gara@antitoxin-gmbh.de](mailto:gara@antitoxin-gmbh.de)

sifin

sifin diagnostics gmbh  
Berliner Allee 317-321  
13088 Berlin / Germany  
☎ +49 30 700 144-0