NUR ZUR IN-VITRO-DIAGNOSTIK

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FOR IN VITRO DIAGNOSTIC USE ONLY

INTENDED USE

Anti-Js^b reagent contains monoclonal antibodies from IgG-Type, specifically directed against the corresponding antigen.

The reagent is used to in-vitro determine qualitative whether red blood cells possess or lack the corresponding blood group antigen Js^b.

The reagent is intended to be used by qualified and technical personnel only.

PRINCIPLE OF PROCEDURE

The test method used with this reagent is based on the principle of column agglutination. Normal human erythrocytes, possessing the appropriate antigen, will be recognized and coated by the corresponding specific antibody and then the cells will be agglutinated by a secondary antibody that reacts with human IgG-molecules.

REAGENTS

The listed reagent contains antibodies of the following clone:

Anti-Js^b for card method monoclonal, IgG clone: PRIM-1

The reagent contain <0.1% (w/v) sodium azide as preservative. Beside the parts active antibody the reagent contain sodium chloride, macromolecules and bovine albumin, which has been tested and certified as disease free by the US Veterinary service inspectors.

WARNING

This reagent is prepared from supernatants of cell cultures.

Nevertheless as biological product it should be looked upon as potentially infectious because of never complete exclusion of danger through excitants of disease.

The reagent contains sodium azide, that may be toxic and may react with lead or copper to form highly explosive salts.

On disposal, flush with large quantities of water.

For the reasons mentioned above the reagent should be handled with proper care.

STORAGE REQUIREMENT

Store opened and unopened products at 2 to 8°C. May be at room temperature while in use. In principle, store and use the reagent to declared expiry date only.

REMARKS

1. With each testing positive and negative controls should be performed.
2. Inappropriate storage impairs efficacy of the reagent.
3. Weak turbidity of the reagent does not affect its reactivity. Bacteria and chemical contamination of the product should be avoided. If a visible change is detected, the reagent should no longer be used, this sign may indicate a microbiological contamination.
4. Strength of positive reactions also depends on age of used blood.
5. Use the Grifols card centrifuge. The use of another card-specific centrifuge (each card centrifuge has its specified unchangeable g-force) may lead to false results due to the changed g-force.
6. The test method identified below is for manual testing only. When using automated or semi-automated instruments, follow the procedures that are contained in the operator's manual provided by the device manufacturer. Laboratories must follow approved validation procedures.
7. For usage of this reagent all effective national laws, directives and guidelines have to be observed, in Germany especially the „Richtlinien zur Gewinnung von Blut und Blutbestandteilen und zur Anwendung von Blutprodukten Hämotherapie“¹ in its actual form.
8. The information on the use of the Grifols test cards in the relevant insert must be observed.

SAMPLE PREPARATION

1. Blood sample should be collected by approved medical procedure.
2. Blood samples to be tested should be used as soon as possible after blood collection to reduce the risk of false-positive and false-negative results due to improper storage or contamination of the sample. If a delay in testing occurs, samples should be stored at 2 to 8°C. Blood drawn into EDTA should be tested within 7 days and samples treated with sodium citrate within 14 days after collection. Blood bag / Donor Blood can be tested until the expiry date.

REAGENT PREPARATION

There is no preparation of the reagent required. Take and use the reagent directly from the vials.

PROCEDURE

Not included but necessary materials:

Card Method

1. Cards: Grifols „DG Gel Coombs“
2. Microliter pipette
3. Tubes
4. Isotonic Saline (0.85 - 0,9% sodium chloride)
5. Centrifuge
6. Grifols Card centrifuge
7. Grifols Card Incubator

Test procedure

Card Method

1. prepare 0,8 % suspensions of red blood cells **in isotonic saline**.
2. Add 50 µL of appropriate cell suspension to each marked micro tube.
3. Add 25 µL of the reagent to each micro tube.
4. Incubate card at 37°C for 15 minutes in a Grifols incubator.
5. Centrifuge the card in a Grifols card centrifuge with the for this centrifuge unchangeable g-force.
6. Check macroscopically for agglutination within 30 minutes.
7. Document the results.

INTERPRETATION OF RESULTS

Positive results (+): visible agglutination of erythrocytes is a positive result and indicates the presence of the corresponding antigen.

Negative results (-): No visible agglutination of erythrocytes is a negative result and indicates the absence of the corresponding antigen.

Read and interpret the results of the card method according to the Grifols card instruction for use.

LIMITATIONS OF THE PROCEDURE

1. Inaccuracy at compliance with instructions written under section “Procedures” and “Interpretation of results” may lead to incorrect results.
2. No valid conclusion concerning the test result can be reached, if controls with uncertain or false results occur.
3. Enzyme treated erythrocytes or addition of bovine albumin and/or other solutions containing protein may cause unspecific reactions.
4. Hemolyzed, turbid, contaminated or clotted blood samples must not be used in this test.
5. Suspensions of red blood cells prepared in the specific dilution medium (DG Sol) lead to false negative reactions. The use of isotonic saline solution is absolutely necessary.
6. Due to variability of antigen expression on human red blood cells, reactivity of the reagent, mentioned above, against certain phenotypes may give weaker reactivity compared to control cells.
7. No one specific antiserum or technique can be guaranteed to detect all rare, weak or variant antigens.²
8. Red blood cells coated with alloantibodies or autoantibodies of the same or similar specificity as that reagent used for the test (i.e. cells that are positive in the direct antiglobulin-test (DAT)) are not suitable for this test procedure.
9. Red blood cells with a positive direct Coombs-test may cause false-positive reactions in card method.
10. It is described in the literature, that samples from patients treated with anti-CD38 monoclonal antibodies can cause false positive results in the Coombs test.⁵
11. Pay attention to all statements to limitations in inserts of Grifols card.
12. Observe that the contents of the reaction chamber(s) are combined. Assure that the reagents remain in the reaction chamber during incubation. There should be no mixing of reactants with reagents in the column prior to centrifugation. Mixing of reactants with reagents in the column can lead to weakened to negative results.

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
PRESENTATION

213238 Anti-Js^b for card method monoclonal, IgG clone: PRIM-1 2 ml

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Antitoxin GmbH - Industriestrasse 88 - 69245 Bammental