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213181 Anti-M for card method monoclonal, mouse clone: LM110/140 5ml

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INTENDED USE

Monoclonal agglutinating Anti-M reagent is produced from cell culture supernatants of hybridoma-cell lines, which secreting antibodies from IgG-Type that reacts specifically with the corresponding antigen. The reagent is used to in-vitro-determine qualitative whether human red blood cells possess or lack the corresponding blood group antigen M. The reagent is intended to be used by qualified and technical personnel only.

PRINCIPLE OF PROCEDURE

The test method used with this reagent is based on the principle of column agglutination. Normal human erythrocytes, possessing the corresponding antigen, will agglutinate in the presence of the specific antibody directed toward the antigen.

REAGENTS

The listed reagent contain antibodies of the following clone:

Anti-M for card method monoclonal, mouse clone: LM110/140

The reagent contain <0.1% (w/v) sodium azide as preservative. Beside the parts active antibody the reagent contain sodium chloride, macromolecules and bovine albumin, which has been tested and certified as disease free by the US Veterinary service inspectors.

WARNING

This reagent was prepared from supernatants of cell cultures. Nevertheless, as biological product, it should be looked upon as potentially infectious because of never complete exclusion of danger through excitants of disease. The reagent contain sodium azide that may be toxic and may react with lead or copper to form highly explosive salts. On disposal, flush with large quantities of water. For the reasons mentioned above the reagent should be handled with proper care.

STORAGE REQUIREMENT

Store opened and unopened products at 2 to 8°C. May be at room temperature while in use. In principle, store and use the reagent to declared expiry date only.

REMARKS

1. With each testing positive and negative controls should be performed.
2. Inappropriate storage impairs efficacy of the reagent.
3. Weak turbidity of the reagent does not affect its reactivity. Bacterial and chemical contamination of the products should be avoided. If a visible change is detected, the reagents should no longer be used this sign may indicate a microbiological contamination.
4. Strength of positive reactions also depends on age of used blood.
5. Use the Grifols card centrifuge. The use of another card-specific centrifuge (each card centrifuge has its specified unchangeable g-force) may lead to false results due to the changed g-force.
6. The test method identified below is for manual testing only. When using automated or semi-automated instruments, follow the procedures that are contained in the operator's manual provided by the device manufacturer. Laboratories must follow approved validation procedures.
7. For usage of this reagent all effective national laws, directives and guidelines have to be observed, in Germany especially the „Richtlinien zur Gewinnung von Blut und Blutbestandteilen und zur Anwendung von Blut-produkten Hämotherapie“⁴ in its actual form.
8. The information on the use of the Grifols test cards in the relevant insert must be observed.

SAMPLE PREPARATION

1. Blood sample should be collected by approved medical procedure.
2. Blood samples to be tested should be used as soon as possible after blood collection to reduce the risk of false-positive and false-negative results due to improper storage or contamination of the sample. If a delay in testing occurs, samples should be stored at 2 to 8°C. Blood drawn into EDTA should be tested within 7 days and samples treated with sodium citrate within 14 days after collection. Blood bag / Donor Blood can be tested until the expiry date.

REAGENT PREPARATION

There is no preparation of the reagents required. Take and use the reagents directly from the vials.

PROCEDURE

Not included but necessary materials:

Card Method

1. Cards: Grifols „DG Gel Neutral“ or „DG Gel Coombs“
2. Microliter pipette
3. Tubes
4. Grifols Card Centrifuge
5. Grifols Card Incubator
6. DG Gel Sol (card specific diluent)

Test procedure

Card Method for Neutral- as well as Coombs-Card

1. Prepare 0,8 % suspensions of red blood cells in DG Gel Sol (card specific diluent).
2. Add 50 µL of appropriate cell suspension to each labeled micro tube.
3. Add 25 µL of appropriate reagent to each micro tube.
4. „DG Gel Neutral“ - without incubation
“DG Gel Coombs“ - incubate card at 37 °C for 15 minutes
5. Centrifugation of card in appropriate centrifuge with the, for the respective centrifuge, unchangeable g-force.
6. Check macroscopically for agglutination within 30 minutes.
7. Document the result.

INTERPRETATION OF RESULTS

Positive results (+): visible agglutination of erythrocytes is a positive result and indicates the presence of the corresponding antigen.
Negative results (-): No visible agglutination of erythrocytes is a negative result and indicates the absence of the corresponding antigen.

Read and interpret the results of the card method according to the Grifols card instruction for use.

LIMITATIONS OF THE PROCEDURE

1. Inaccuracy at compliance with instructions written under section “Procedures” and “Interpretation of results” may lead to incorrect results.
2. No valid conclusion concerning the test result can be reached, if controls with uncertain or false results occur.
3. Enzyme treated erythrocytes or addition of bovine albumin and/or other solutions containing protein may cause unspecific reactions.
4. Hemolyzed, turbid, contaminated or clotted blood samples must not be used in this test.
5. Due to variability of antigen expression on human red blood cells, reactivity of the reagent mentioned above, against certain phenotypes, may give weaker reactivity compared to control cells.
6. No one specific antiserum or technique can be guaranteed to detect all rare, weak or variants of antigens.²
7. Red blood cells coated with alloantibodies or autoantibodies of the same or similar specificity as that reagent used for the test (i.e. cells that are positive in the direct antiglobulin-test (DAT)) are not suitable for this test procedure.
8. Red blood cells with a positive direct Coombs-test may cause false-positive reactions.
9. It is described in the literature, that samples from patients treated with anti-CD38 monoclonal antibodies can cause false positive results in the Coombs test.⁵
10. Pay attention to all statements to limitations in inserts of Grifols card.
11. In rare cases, Anti-M for card method may cause weak positive reactions with M-negative cells. In case of weak positive reactions, it should be confirmed with further tests.

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Antitoxin GmbH - Industriestrasse 88 - 69245 Bammental - Nemecko