Phosphate Buffered Saline, Concentrate

FOR IN VITRO DIAGNOSTIC USE ONLY

INTENDED USE

Phosphate Buffered Saline, Concentrate is a 10-fold concentrated phosphate buffer solution. After dilution with demineralised water to a ready-to-use solution, it is suitable for use as a buffer solution or washing solution in in vitro diagnostic procedures that require a phosphate buffered solution. The product is intended to be used by qualified and technical personnel only.

MAIN COMPONENTS

Sodium chloride (NaCl)	87.8 g/L
Di-Sodium hydrogenphosphate-Dihydrate ($\text{Na}_2\text{HPO}_4 \cdot 2\text{H}_2\text{O}$)	27.8 g/L
Potassium chloride (KCl)	4.2 g/L
Sodium fluoride (NaF)	2.5 g/L
Sodium dihydrogen phosphate-Monohydrate ($\text{NaH}_2\text{PO}_4 \cdot \text{H}_2\text{O}$)	2.9 g/L

SPECIFICATIONS

Colour:	colourless
Form:	liquid
pH value of 1:10 dilution:	7.3 – 7.5
Osmolality of 1:10 dilution:	335 – 345 mOsmol/kg

BUFFER PREPARATION

Transfer the complete volume (1 L) into a 10 L container and add 9 L of demineralised water. Ensure that the buffer is mixed.

ADDITIONAL REQUIRED MATERIAL NOT SUPPLIED

- 10 litre container for diluting the concentrate.
- Measuring device suitable for 9 L of demineralised water.
- demineralised water
- 0.1 M NaOH or 70% isopropanol for cleaning

REMARKS

1. Regularly empty and clean the container used for the diluted solution to prevent microbial contamination. Use 0.1 M NaOH or 70% isopropanol for cleaning, followed by rinsing with demineralised water. Do not use bleach.
2. Replace the container used for the diluted solution after a maximum of one year of use.
3. When withdrawing the diluted solution, use clean equipment only and avoid contaminating the remaining stock.
4. It is recommended to use a drain tap on the canister to avoid contamination during withdrawal.
5. Do not reuse dispensed solution. Discard any portion that has been removed from the canister.
6. Salt crystals may form at low temperatures. If crystals are visible, warm the solution until fully dissolved before use.
7. If turbidity or particulate matter is visible, discard the solution.

STORAGE REQUIREMENTS

Store at +10 to +30 °C in a dry place, protected from sunlight. Do not freeze.

Unopened concentrate containers are stable until the declared expiry date (YYYY-MM-DD).

Once diluted, the ready-to-use solution must be used within 3 months of preparation. Record the preparation date on the container.

DISPOSAL

Dispose of all solutions, cleaning agents, and containers in accordance with applicable local regulations.

INCIDENTS RELATED TO THE DEVICE

Any serious incident that has occurred in relation to the device must be reported to the manufacturer and the competent authority of the Member State where the user and/or patient is established.

GLOSSARY OF SYMBOLS

The following symbols are used in the labelling of the product:

 REF	Product code	 LOT	Batch
	Store from - to		Expiry Date
 IVD	In-Vitro Diagnostic		EU CE symbol
	Manufacturer according to (EU) 2017/746		Consult instruction for use
 UDI	Unique Device Identification		

REF

05.800-1K.VI Phosphate Buffered Saline, Concentrate 6x1L

 ANTITOXIN GmbH Industriestraße 88 69245 Bammental Germany
 +49 (0) 6223/ 8661-0
 +49 (0) 6223/ 8661-13
 gara@antitoxin-gmbh.de

 05.800- / Version R002 / 2026-05-26

Marking of changes:

Changes to the previous version are shaded grey.

Phosphate Buffered Saline, Concentrate

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ANTITOXIN GmbH Industriestraße 88 69245 Bammental Tyskland