Phosphate Buffered Saline, Isotonic

FOR IN VITRO DIAGNOSTIC USE ONLY

INTENDED USE

This phosphate buffered saline is a ready-to-use solution intended for use by qualified and technical personnel only as a buffer solution or washing solution suitable for in vitro diagnostic procedures that require a phosphate buffered solution.

MAIN COMPONENTS

Sodium chloride (NaCl)	8.78 g/L
Di-Sodium hydrogenphosphate-Dihydrate ($\text{Na}_2\text{HPO}_4 \cdot 2\text{H}_2\text{O}$)	2.78 g/L
Potassium chloride (KCl)	0.42 g/L
Sodium fluoride (NaF)	0.25 g/L
Sodium dihydrogen phosphate-Monohydrate ($\text{NaH}_2\text{PO}_4 \cdot \text{H}_2\text{O}$)	0.29 g/L

SPECIFICATIONS

Colour:	colourless
Form:	liquid
pH value:	7.3 – 7.5
Osmolality:	335 – 345 mOsmol/kg

BUFFER PREPARATION

The buffer is ready to use.

ADDITIONAL REQUIRED MATERIAL NOT SUPPLIED

No additional material is required.

REMARKS

1. Do not reuse dispensed solution. Discard any portion that has been removed from the canister.
2. When withdrawing buffer, use clean equipment only and avoid contaminating the remaining stock.
3. It is recommended to use a drain tap on the canister to avoid contamination during withdrawal.
4. If turbidity or particulate matter is visible, discard the solution.

STORAGE REQUIREMENTS

Store at +10 to +30 °C in a dry place, protected from sunlight. Do not freeze.

Once in use, the solution must be used within 3 months of first use. Record the date of first use on the container.

Do not use the solution beyond the declared expiry date (YYYY-MM-DD).

DISPOSAL

Dispose of all solutions and containers in accordance with applicable local regulations.

INCIDENTS RELATED TO THE DEVICE

Any serious incident that has occurred in relation to the device must be reported to the manufacturer and the competent authority of the Member State where the user and/or patient is established.

GLOSSARY OF SYMBOLS

The following symbols are used in the labelling of the product:

 REF	Product code	 LOT	Batch
	Store from - to		Expiry Date
 IVD	In-Vitro Diagnostic		EU CE symbol
	Manufacturer according to (EU) 2017/746		Consult instruction for use
 UDI	Unique Device Identification		

REF

05.801-10K	Phosphate Buffered Saline, Isotonic	1x 10L
05.801-5K.II	Phosphate Buffered Saline, Isotonic	2x 5L

CE

 ANTITOXIN GmbH Industriestraße 88 69245 Bammental Germany
 +49 (0) 6223/ 8661-0
 +49 (0) 6223/ 8661-13
 gara@antitoxin-gmbh.de

 05.801- / Version R002 / 2026-05-26

Marking of changes:

Changes to the previous version are shaded grey.

Phosphate Buffered Saline, Isotonic

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ANTITOXIN GmbH Industriestraße 88 69245 Bammental Tyskland