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For the indirect Coombs Test
FOR IN VITRO DIAGNOSTIC USE ONLY

INTENDED USE

The reagent is used in vitro to qualitatively determine whether human red blood cells possess or lack the corresponding blood group antigen W_r^a . The reagent is intended to be used by qualified and technical personnel only to perform immunohematology screening tests as part of the practice of transfusion medicine in the general population. The test method used with this reagent is based on the principle of agglutination, performed manually on human blood samples.

INDICATION / CONTRA-INDICATION

The Anti- W_r^a monoclonal Coombs-reactive blood grouping reagent is used to test patient or donor red cells for the presence of the W_r^a antigen. Typing of donor cells facilitates the selection of suitable antigen-negative units for transfusion to patients with corresponding antibody. Cell typing also serves as final verification process for the identification of Anti- W_r^a in patient or donor sera. The product was validated using samples collected in Europe from patients of unknown ethnic background.

Contraindication: The reagent is not validated for neonatal blood samples and must not be used with such specimens.

The approximate frequencies of W_r^a antigen¹:

Phenotype	All populations
W_r^a+	<0.01%
W_r^b+	100%

REAGENTS

The reagent is derived from cell culture supernatant of a hybridoma cell line secreting a specific antibody of IgG-type, which reacts specifically with the corresponding antigen. The monoclonal blood group reagent is prepared from the BGU1-WR clone. The reagent contains < 0.1% (w/v) sodium azide as preservative. Additionally, the reagent contains sodium chloride, macromolecules and bovine albumin (BSA). The BSA used in this formulation is derived from US-sourced animals from USDA and APHIS-approved facilities. It complies with Regulations (EC) No 1069/2009 and (EU) No 142/2011 for use in in vitro diagnostic reagents and has been certified free from both Vesicular Stomatitis Virus and Bluetongue Virus.

WARNING

This reagent is prepared from cell culture supernatant. As a biological product, this reagent should be regarded as potentially infectious, since the risk of pathogen transmission can never be completely excluded. The reagent contains sodium azide, a toxic substance that may react with lead or copper to form highly explosive salts. Upon disposal, flush with large quantities of water. For the reasons mentioned above, the reagent should be handled with proper care.

STORAGE REQUIREMENT

Store opened and unopened products at +2 to +8 °C. It may be kept at room temperature while in use. In-use stability testing demonstrated that 30 cycles of 2-hour storage at room temperature did not compromise the qualitative test results until the indicated expiry date. Use only until the declared expiry date, which is given in the format year-month-day (YYYY-MM-DD).

REMARKS

- The strength of positive reactions also depends on the age of the used blood.
- With each testing positive and negative controls should be performed.
- The Coombs /AHG serum is to be tested for reactivity/functionality using positive and negative controls.
- Inappropriate storage impairs efficacy of the reagent.
- Centrifugation outside the specified speed range may lead to false results.
- The test method described below is intended for manual testing only and must be carried out according to the instructions for use.
In case of
a) changes in technology or deviations from the instructions for use;
b) use of automated or semi-automated systems;
the laboratories must follow the procedures specified in the operator's manual provided by the device manufacturer and perform validation according to recognized procedures.
- The reagent must be used in compliance with all applicable national laws, directives and guidelines currently in force, in particular, for Germany, the "Richtlinien zur Gewinnung von Blut und Blutbestandteilen und zur Anwendung von Blutprodukten (Hämotherapie)".²
- Slight turbidity does not affect the reactivity of the reagent. Protect the reagent from bacterial and chemical contamination. Do not use the reagent if a visible change, such as increased turbidity or a colour change, is observed, as this may indicate microbiological contamination.
- Do not use damaged (e.g. leaking or broken bottles as well as broken dropper pipettes) or unlabeled products.
- For all materials required for use but not supplied with this reagent, the respective instructions for use and maintenance requirements must be followed.

SAMPLE PREPARATION

- Blood samples should be obtained using an appropriate phlebotomy technique. Sample drawn into tubes containing EDTA, sodium citrate, CPD-A, ACD or blood bags containing PAGGS-M may be used.
- Blood samples to be tested should be used immediately after collection to reduce the risk of false-positive and false-negative results due to improper storage or sample contamination. Samples that cannot be tested immediately should be stored at +2 to +8 °C. Blood drawn into EDTA must be tested within 7 days after collection. Samples treated with sodium citrate, CPD-A or ACD must be tested within 14 days after collection. Donor blood collected in blood bags containing PAGGS-M can be tested until the expiry date.
- Do not freeze the samples.

REAGENT PREPARATION

No preparation of the reagent is required.
Use the reagent directly from the vials.

TEST PROCEDURE

Tube Centrifugation Method:

Materials required but not supplied:

- Tubes (10 x 75 mm or 12 x 75 mm)
- Microliter pipette
- Timer
- Incubator
- Centrifuge
- Isotonic saline (0.85 – 0.9 % sodium chloride)
- AHG serum blend ANTITOXIN REF: 03.140-10

Workflow:

- Prepare 2 to 5 % suspensions of red blood cells in isotonic saline. Red blood cells may be washed 1 - 3 times with isotonic saline.
- Add 100 µl (alternative: one drop = approximately 50 µl) of appropriate reagent to each tube.
- Add 100 µl (alternative: one drop = approximately 50 µl) of appropriate cell suspension to each tube.
- Mix well by gentle shaking.
- Incubate tube in an incubator at +37 °C for 30 min.
- Wash red cells 3 times with (cold) isotonic saline.
- Add 100 µl of anti-human globulin reagent (Coombs serum/AHG serum) to the tube. Gently shake the tube to release the cells from the bottom and mix the cells with the serum.
- Centrifuge tube for 1 minute at 800-1000 x g.
- Gently shake the red blood cells and check macroscopically for agglutination within 3 minutes.
- Document the result.

INTERPRETATION OF RESULTS

- Positive results (+): If agglutination of the erythrocytes occurs within the accepted limitations of the procedure, the test result is to be interpreted as positive and indicates the presence of the corresponding antigen.
- Negative results (-): If no agglutination of the erythrocytes occurs within the accepted limitations of the procedure, the test result is to be interpreted as negative and the corresponding antigen is not detectable.

The reading and interpretation of the results after "gentle shaking" using the Tube Centrifugation Method:

Negative	No detectable agglutinates, homogeneous red coloration of the suspension.
Positive	One complete agglutinate.
	No complete agglutination; individual agglutinates visible.
	Red coloration of the suspension with small/ miniature agglutinates.

LIMITATIONS OF THE PROCEDURE

- Failure to comply with the instructions provided in the sections "TEST PROCEDURE" and "INTERPRETATION OF RESULTS" may lead to incorrect results.
- If the controls yield invalid or incorrect results, do not interpret the test results and repeat the test.
- The use of enzyme-treated erythrocytes and the addition of LISS and/or BSA have not been validated and may cause non-specific reactions.
- Commercially available test cells may contain stabilization solutions that differ from the anticoagulants validated for this reagent and may lead to incorrect results.
- Hemolyzed, turbid, contaminated or clotted blood samples must not be used in this test.
- A red blood cell suspension that deviates from the specified concentration may lead to false positive or false negative results.
- The addition of volumes that deviate from the volumes specified in the method may lead to altered reaction behaviour.
- Due to variability in antigen expression on human red blood cells, the reactivity of the reagent against certain phenotypes may be weaker than that observed in control cells.
- No single antiserum or technique can guarantee the detection of all rare, weak or variant antigens.³
- Therapeutic monoclonal antibodies (e.g. those targeting CD38) may interfere with serological testing.⁴
- Red blood cells coated with alloantibodies or autoantibodies (i.e., cells that are positive in the direct antiglobulin-test (DAT)) are unsuitable. They may cause false-positive reactions, even in the absence of the reagent.
- Loss of reactivity of the anti-human globulin serum (Coombs serum / AHG serum i.e. the second antibody that recognises human IgG molecules) due to an incorrect test procedure such as insufficient washing with a possibly too warm physiological saline solution after incubation (see point 6 of the test procedure).
- Adding an insufficient volume of anti-human globulin serum (Coombs serum / AHG serum) that differs from the specified volume can lead to a weaker or negative reaction.
- This reagent has been validated with the AHG serum blend ANTITOXIN REF: 03.140-10 mentioned in the "TEST PROCEDURE" section. Using a different AHG serum may lead to false negative results.

INCIDENTS RELATED TO THE DEVICE

Any serious incident that has occurred in relation to the device must be reported to the manufacturer and the competent authority of the Member State in which the user and / or patient is established.

PERFORMANCE CHARACTERISTICS

A performance evaluation for the device was carried out.
The required blood samples were tested and compared with other reference methods / devices.

device Method	Anti-Wr(a) monoclonal IgG			
	Positive samples	Sensitivity	Negative samples	Specificity
Tube Centrifugation Method	8 / 8	100%	102 / 102	100%

Diagnostic Sensitivity: The probability that the device gives a positive result in the presence of the target marker.
Diagnostic Specificity: The probability that the device gives a negative result in the absence of the target marker.

Anti-Wr(a) monoclonal IgG is equivalent and does not differ in quality from comparable reagents available on the market.

DIFFERENCES BETWEEN BATCHES

Validation between three batches throughout the entire shelf life demonstrated no differences in performance.

INTERFERENCE STUDY

The interference studies showed no impairment of the qualitative tests when the following potentially interfering substances were used at the concentrations listed below:
Heparin 720 U/dl, Albumin 15000 mg/dl, Triglycerides 1500 mg/dl, Bilirubin 40 mg/dl, Ethanol 620 mg/dl, Glucose 1000 mg/dl.
For the anticoagulants and additive solutions used (EDTA, sodium citrate, CPD-A, ACD, PAGGS-M), three times the respective recommended concentration was tested.

SUMMARY OF SAFETY AND PERFORMANCE

The Summary of Safety and Performance of the device is available via ANTITOXIN (www.antitoxin-gmbh.de) and can be accessed via the EUDAMED database.

LITERATURE

1. Reid ME, Lomas-Francis C, Olsson ML: The Blood Group Antigen FactsBook, 3rd ed. Academic Press, 2012
2. Richtlinien zur Gewinnung von Blut und Blutbestandteilen und zur Anwendung von Blutprodukten (Hämotherapie)
3. CLSI, I/LA33-A Validation of Automated System for Immunhematological Testing Before Implementation; Approved Guideline CLINICAL AND LABORATORY STANDARDS INSTITUTE Dezember 2009
4. Jones, Andrew D et al. "Impact of new myeloma agents on the transfusion laboratory." Pathology vol. 53,3 (2021): 427-437.
5. Peter D. Issit, David J. Anstee Applied Blood Group Serology, 4. edition, Montgomery Scientific Publications 1998
6. Christian Mueller-Eckhardt, Volker Kiefel Transfusionsmedizin, 3. Auflage, Springer-Verlag 2004.
7. Körmöczy, G. et al. Immunhämatologische Untersuchungen bei Patienten. Österreichische Gesellschaft für Blutgruppenserologie, Transfusionsmedizin, Regenerative Medizin und Immungenetik (ÖGBT), Ausgabe Mai 2015.

GLOSSARY OF SYMBOLS

The following symbols may be used in the labelling of the device:

 REF	Product code	 LOT	Batch
	Store from - to		Expiry Date
 IVD	In-Vitro Diagnostic		EU CE symbol
	Manufacturer according to (EU) 2017/746		Consult instruction for use
 UDI	Unique Device Identification		Distributor

69260 Anti-Wr(a) monoclonal IgG 2 ml

CE 0483

 ANTITOXIN GmbH Industriestraße 88 69245 Bammental Germany

 +49 (0) 6223/ 8661-0
 +49 (0) 6223/ 8661-13
 gara@antitoxin-gmbh.de

 01.255- / version R003 / 2026-06-18

Marking of changes:
Changes to the previous version are shaded grey.

In case of any inconsistencies between different language versions of the IFU, the English version shall be considered the governing text.

ANTITOXIN GmbH Industriestraße 88 69245 Bammental Germany

BAG Diagnostics GmbH
Amtsgerichtsstr. 1-5
35423 Lich/Germany
A BAG Group company

T. +49 (0) 6404/925-100
F. +49 (0) 6404/925-460
M. info@bag-diagnostics.com
W. www.bag-diagnostics.com

Ordering: order@bag-diagnostics.com
T. +49 (0) 6404 / 925 – 450
F. +49 (0) 6404 / 925 - 460
M. order@bag-diagnostics.com

Customer Service:
T. +49 (0) 6404 / 925 - 125
F. +49 (0) 6404 / 925 - 421
M. service@bag-diagnostics.com