Anti-D blend monoclonal, human IgG + IgM

Deutsch

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Anti-D blend monoclonal, human IgG + IgM

English

For Slide-, Tube-, Card- and Microplate Method
FOR IN VITRO DIAGNOSTIC USE ONLY

INTENDED USE

Monoclonal Anti-D blend reagent is produced from cell culture supernatants of heterohybridoma - cell lines. The cells are secreting an antibody of IgG- or IgM-type which reacts specific with the corresponding blood group antigen. The antibody in each case is human protein. The reagent is used for In-Vitro-Diagnostic, to determine qualitative whether human red blood cells possess or lack the corresponding blood group antigen D. This reagent reacts negative with erythrocytes of category D^{VI} with the IgM-antibody and positive with the IgG-antibody in indirect Coombs-test.

The reagent is intended to be used by qualified and technical personnel only.

PRINCIPLE OF PROCEDURE

The test methods used with this reagent is based on the principle of agglutination. Normal human erythrocytes, possessing the corresponding antigen, will agglutinate in the presence of the specific antibody directed towards the antigen.

REAGENT

The listed reagent contains antibodies of the following cell clones:

Anti-D blend monoclonal, human IgG+IgM clones: MS-26, TH-28

This reagent contains <0.1% (w/v) sodium azide as preservative.

Additionally the reagent is prepared of active antibody, sodium chloride, macromolecules and bovine albumin, which has been tested and certified by the US Veterinary service inspectors.

WARNING

This reagent is prepared from supernatants of cell cultures. As biological product it should be looked upon as potentially infectious because of never complete exclusion of danger through excipients of disease. The reagent contains sodium azide that may be toxic and may react with lead or copper to form highly explosive salts.

On disposal, flush with large quantities of water.

For the reasons mentioned above, reagent should be handled with proper care.

STORAGE REQUIREMENT

Store opened and unopened product at +2 to +8°C. May be at room temperature while in use. In principle, store and use the reagent to indicated expiry date only.

REMARKS

1. With each testing positive and negative controls should be performed.
2. Inappropriate storage impairs efficacy of the reagent.
3. Weak turbidity of the reagent does not affect its reactivity. Bacteria and chemical contamination of the product should be avoided. If a visible change is detected, the reagent should no longer be used, this sign may indicate a microbiological contamination.
4. Strength of positive reactions also depends on age of used blood
5. Centrifugation outside the specified speed range may lead to false results.
With the card method, the use of another card-specific centrifuge (each card centrifuge has its specified unchangeable g-force) may lead to false results due to the changed g-force.
6. The test methods identified below are for manual testing only. When using automated or semi-automated instruments, follow the procedures that are contained in the operator's manual provided by the device manufacturer. Laboratories must follow approved validation procedures.
7. For usage of this reagent all effective national laws, directives and guidelines have to be observed, in Germany especially the „Richtlinien zur Gewinnung von Blut und Blutbestandteilen und zur Anwendung von Blutprodukten (Hämotherapie)“ in its current version.
8. This reagent reacts negative with erythrocytes of category D^{VI} with the IgM-antibody and positive with the IgG-antibody in indirect Coombs-test.
9. The information on the use of the test cards in the relevant insert must be observed.

SAMPLE PREPARATION

1. Blood sample should be collected by approved medical procedure.
2. Blood sample to be tested should be used as soon as possible after blood collection to reduce the risk of false-positive and false-negative results due to improper storage or contamination of the sample.
If a delay in testing occurs, sample should be stored at +2 to +8°C.
Blood drawn into EDTA should be tested within 7 days and sample treated with sodium citrate within 14 days after collection.
Blood bag / Donor Blood can be tested until the expiry date.

REAGENT PREPARATION

There is no preparation of the reagent required.
Take and use reagent directly from the vials.

PROCEDURE

Not provided material, additionally needed

Slide Method:	glass slide; Pasteur pipette; mixing stick; timer; tubes (10 x 75 mm or 12 x 75 mm);
Tube Centrifugation Method:	microliter pipette; centrifuge; timer; incubator; isotonic saline (0,85 - 0,9% sodium chloride); Anti-Human-Globulin-Serum (Coombs-Serum / AHG-Serum)
Card Method:	Cards: - BIO-RAD (DiaMed) „NaCl, Enzyme Test and cold Agglutinins“ - BIO-RAD (DiaMed) „LISS/Coombs“ - Grifols „DG Gel Neutral“ - Grifols „DG Gel Coombs“ tubes; microliter pipette; centrifuge; timer; corresponding card centrifuge; card incubator; card specific diluent; isotonic saline (0,85 - 0,9% sodium chloride)
Microplate Method:	microplate with 96 U-wells; microliter pipette; centrifuge; timer; microplate shaker; microplate centrifuge; isotonic saline (0,85 - 0,9% sodium chloride)

Test procedure

Slide Method

1. Use erythrocyte sediment only.
2. Place one drop (approximately 50 µL) of the reagent on a marked glass slide.
3. Using a Pasteur pipette add one drop erythrocyte sediment (approximately 50 µL) to the drop of reagent on the glass slide.
4. Mix the erythrocytes with reagent well with a stick and spread to a circle with a diameter of approximately 2 cm.
5. By slightly rotating the slide, check for agglutination within 1 minute (reaction starts within seconds).
6. Document the result.

Tube Centrifugation Method

1. Prepare a 2% to 5% suspension of red blood cells in isotonic saline (red blood cells may be washed 1-3 times with isotonic saline).
2. At first put 100 µL of appropriate reagent in a marked tube subsequently add 100 µL of appropriate cell suspension in the tube.
Alternative one drop = approximately 50 µL cell suspension can be added to one drop = approximately 50 µL test serum.
3. Mix Erythrocytes-/Reagent mixture well by slightly shaking.
4. Incubate tube at room temperature for 1-15 min.
5. Centrifugation of tube for 1 minute at 2.000 rpm (approximately 800-1.000 x g).
6. Gently shake the red cells completely from the bottom of the tube and check macroscopically for agglutination within 3 minutes.
7. Document the result

If the results are negative, continue with the indirect Coombs

8. Repeat step 1. to 3. with fresh material or continue after reading of result directly.
9. Incubate tube at +37 °C for 30 min.
10. Wash red cells 3 times with (cold) isotonic saline.
11. Add 100 µL Anti-Human-Globulin-serum (Coombs-Serum / AHG-Serum) to the tube, release the cells from the bottom of the tube by slightly shaking and mix with the Coombs-Serum / AHG-Serum.
12. Centrifugation of tube for 1 min at 1.000 rpm (approximately 180-270 x g).
13. Gently shake the red cells completely from the bottom of the tube and check macroscopically for agglutination within 3 minutes.
14. Document the result.

Card Method (manual method / valid for the cards:

- BIO-RAD (DiaMed) „NaCl, Enzyme Test and cold Agglutinins“
- Grifols „DG Gel Neutral“)

1. Prepare 0,8 % suspension of red blood cells in card specific diluent (red blood cells may be washed 1-3 times with isotonic saline).
2. Add 50 µL of appropriate cell suspension to a marked micro tube.
3. Add 25 µL of reagent to the micro tube.
4. Centrifuge the card in appropriate card centrifuge with the for this centrifuge unchangeable g-force after latest 30 minutes.
5. Check macroscopically micro tube for agglutination within 30 minutes.
6. Document the result.

If the results are negative, continue with the Coombs-card (indirect Coombs-test)

(manual method / valid for the cards:

- BIO-RAD (DiaMed) „LISS/Coombs“
- Grifols „DG Gel Coombs“)

1. Prepare 0,8 % suspension of red blood cells in card specific diluent (red blood cells may be washed 1-3 times with isotonic saline).
2. Add 50 µL of appropriate cell suspension to a marked micro tube.
3. Add 25 µL of reagent to the micro tube.
4. Incubate card at +37°C for 15 minutes.
5. Centrifuge the card in appropriate card centrifuge with the for this centrifuge unchangeable g-force.
6. Check macroscopically micro tube for agglutination within 30 minutes.
7. Document the result.

Microplate Method

1. Prepare 2-5% suspension of red blood cells in isotonic saline (red blood cells may be washed 1-3 times with isotonic saline).
2. Add 50 µL of reagent to a marked well.
3. Add 50 µL of appropriate cell suspension to the well.
4. Shake the microplate for 30 seconds on microplate shaker with medium speed.
5. Centrifuge the microplate in appropriate microplate centrifuge for 30 seconds at 400 x g.
6. Shake microplate for 30 seconds on a microplate shaker with medium speed.
7. Evaluated macroscopically for agglutination directly after shaking.
8. Document the result.
9. Tests with negative or doubtful results have to be incubated for 5 to 10 minutes at room temperature.
10. Repeat steps 5 to 8 after incubation at room temperature

INTERPRETATION OF RESULTS

"Slightly rotating / shaking"

at Slide Method / at Tube Centrifugation Method / Microplate Method:

Positive results (+): visible agglutination of erythrocytes is a positive result and indicates the presence of the corresponding antigen.

Negative results (-): No visible agglutination of erythrocytes is a negative result and indicates the absence of the corresponding antigen.

Read and interpret the results of the card according to the instruction for use of the corresponding card.

LIMITATIONS OF THE PROCEDURE

1. Inaccuracy at compliance with instructions written under section "Procedures" and "Interpretation of results" may lead to incorrect results.
2. No valid conclusion concerning the test result can be reached, if controls with uncertain or false results occur.
3. Enzyme treated erythrocytes or addition of bovine albumin and/or other solutions containing protein may cause unspecific reactions with this reagent.
4. Hemolyzed, turbid, contaminated or clotted samples should not be used.
5. With the slide method, unspecific reactions might appear due to drying of the reaction-formation or if the slide is heated
6. With the slide method most of weak D-antigens (Dweak) and categories are not recognized.
7. This reagent reacts negative with erythrocytes of category D^{VI} with the IgM-antibody and positive with the IgG-antibody in indirect Coombs-test.
8. Due to variability of antigen expression, reactivity of this reagent against certain phenotypes may give weaker reactivity compared to control cells.
9. No one specific antiserum or technique can be guaranteed to detect all rare, weak or variant antigens².
10. Red blood cells coated with alloantibodies or autoantibodies of the same or similar specificity as the reagent (i.e., cells that are positive in the direct antiglobulin-test (DAT)) are not suitable for this test procedure.
11. Red blood cells coated with antibodies (cells that are positive in the direct antiglobulin-test (DAT)) may give false-positive results in card method. These cells react positive without testserum too.
12. As described in the literature, samples from patients treated with anti-CD38 monoclonal antibodies can cause false positive results in the Coombs test.⁵
13. Pay attention to all statements to limitations in the instruction for use of the cards.

BAG
DIAGNOSTICS

BAG Diagnostics GmbH
Amtsgerichtsstr. 1-5
35423 Lich/Germany
A BAG Group company

T. +49 (0) 6404/925-100
F. +49 (0) 6404/925-460
M. info@bag-diagnostics.com
W. www.bag-diagnostics.com

Ordering:
T. +49 (0) 6404 / 925 - 450
F. +49 (0) 6404 / 925 - 460
M. order@bag-diagnostics.com

Customer Service:
T. +49 (0) 6404 / 925 - 125
F. +49 (0) 6404 / 925 - 421
M. service@bag-diagnostics.com

PERFORMANCE CHARACTERISTICS

In compliance with Common Technical Specifications (CTS Commission Decision of 03. February 2009) a performance evaluation was conducted. Different samples (donor, patient, neonates-, panel blood) was used and compared with other reference methods / products.

Product	Positive	Negative	False Positive	False Negative
Anti- D blend	2164	590	0	0 *

The calculated values found for:

Product	Sensitivity	Spezifitcity
Anti- D blend	100% *	100%

* for Anti-D blend monoclonal, human IgG+IgM, clones: MS-26; TH-28 in slide method 4 D^{weak} samples were tested false negative, so in this method sensitivity was 99,82% only.

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REF Product Code	LOT Batch
Store from - to	Expiration Date
IVD In-Vitro Diagnostic	CE EU CE symbol
Manufacturer according to 98/79/EG	Consult instructions for use
UDI Unique Device Identification	Distributor

REF

69210 Anti-D blend monoclonal, human IgG+IgM clones: MS-26, TH-28 10 ml
69211 Anti-D blend monoclonal, human IgG+IgM clones: MS-26, TH-28 10 x 10ml

730-13-2414 Version 014 / 15. January 2022

CE 0483

BAG Diagnostics GmbH
 Amtsgerichtsstr. 1-5
 35423 Lich/Germany
A BAG Group company

T. +49 (0) 6404/925-100
F. +49 (0) 6404/925-460
M. info@bag-diagnostics.com
W. www.bag-diagnostics.com

Ordering:
T. +49 (0) 6404 / 925 - 450
F. +49 (0) 6404 / 925 - 460
M. order@bag-diagnostics.com

Customer Service:
T. +49 (0) 6404 / 925 - 125
F. +49 (0) 6404 / 925 - 421
M. service@bag-diagnostics.com