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# INSTRUCTIONS FOR USE

## Anti-H (A<sub>2</sub>), monoclonal IgM (mouse) OrthoClone

**REF** 690091A

### INTENDED USE

Monoclonal agglutinating Anti-H (A<sub>2</sub>) reagent is produced from cell culture supernatants of hybridoma-cell line. The cells are secreting an antibody of IgM-type. The reagent is used to in-vitro determine qualitative whether red blood cells possess or lack free H-substance on human red blood cells. It is used to differentiate between A<sub>1</sub> and A<sub>2</sub> blood. There is no free H-substance on A<sub>1</sub> A<sub>1</sub>-blood, but with A<sub>1</sub>B-, A<sub>1</sub>O- and especially with BO-types there may still be free rests of H-substance and show weak reactions. Strong reactions are visible with A<sub>2</sub>-blood or weaker A-variants. As defined strongest reaction are seen with O-bloods. The reagent is intended to be used by qualified and technical personnel only.

### PRINCIPLE OF PROCEDURE

The test method used with this reagent is based on the principle of agglutination. Normal human erythrocytes which have free H-substance will agglutinate in the presence of the specific antibody.

### REAGENT

The listed reagent contains antibodies of the following clone:

Anti- H monoclonal, mouse IgM clone: 10934C11

The reagent contains <0.1% (w/v) sodium azide as preservative. Beside the part active antibody, the reagent contains sodium chloride, macromolecules and bovine albumin, which has been tested and certified by the US Veterinary service inspectors.

**CAUTION:** The reagent is prepared from supernatants of cell cultures. Regardless, as biological product it should be looked upon as potentially infectious because of never complete exclusion of danger through excitants of disease. The reagent contains sodium azide that may be toxic and may react with lead or copper to form highly explosive salts.

On disposal, flush with large quantities of water.

For the reasons mentioned above, the reagent should be handled with proper care.

### STORAGE REQUIREMENTS

Store at +2 to +8 °C (unopened / opened) or at room temperature while in use. Do not use reagent beyond its labelled expiration date.

### REMARKS

1. With each testing positive and negative controls should be performed.
2. Inappropriate storage impairs efficacy of the reagent.
3. Weak turbidity of the reagent does not affect its reactivity. Bacteria and chemical contamination of the product should be avoided. If a visible change is detected, the reagent should no longer be used, this sign may indicate a microbiological contamination.
4. Strength of positive reactions also depends on age of used blood
5. Centrifugation outside the specified speed range may lead to false results.
6. The test method identified below is for manual testing only. When using automated or semi-automated instruments, follow the procedures that are contained in the operator's manual provided by the device manufacturer. Laboratories must follow approved validation procedures.
7. For usage of the reagent all effective national laws, directives and guidelines in its current version have to be observed, in Germany especially the "Richtlinien zur Gewinnung von Blut und Blutbestandteilen und zur Anwendung von Blutprodukten (Hämotherapie)"<sup>1</sup>.

### SAMPLE PREPARATION

1. Blood sample should be collected by approved medical procedure.
2. Blood sample to be tested should be used as soon as possible after blood collection to reduce the risk of false-positive and false-negative results due to improper storage or contamination of the sample.  
If a delay in testing occurs, sample should be stored at +2 to +8°C.  
Blood drawn into EDTA should be tested within 7 days and sample treated with sodium citrate within 14 days after collection.  
Blood bag / Donor Blood can be tested until the expiry date.

### REAGENT PREPARATION

There is no specific preparation of the reagent required. Take and use reagent directly from the vial.

### PROCEDURE

Materials required but not provided:

Tube Centrifugation Method

1. Test tubes, 10 x 75 mm or 12 x 75 mm
2. Pipettes designed to deliver 50 µL/100 µL
3. Disposable pipette tips
4. Timer
5. Centrifuge
6. Isotonic saline (0.85 – 0.9% sodium chloride)

### Test procedure

#### Tube Centrifugation Method

1. Prepare 2% to 5% suspensions of red blood cells in isotonic saline (red blood cells may be washed 1-3 times with isotonic saline).
2. Add 100 µL (alternative: one drop, approximately 50 µL) of reagent to a marked tube.
3. Add 100 µL (alternative: one drop, approximately 50 µL) of appropriate cell suspension to the tube.
4. Mix well by shaking slightly.
5. Incubate tube at room temperature for 15 min.
6. Centrifugate tube for 1 minute at 2.000 rpm (approximately 800 – 1.000 g).
7. Gently shake the red cells from the bottom of the tube and check macroscopically for agglutination within 3 minutes.
8. Document the result.

### INTERPRETATION OF RESULTS

"Slightly shaking" at Tube Centrifugation Method.

**Positive result (+):** visible agglutination of erythrocytes indicates the presence of the corresponding antigen.

**Negative result (-):** no visible agglutination of erythrocytes indicates the absence of the corresponding antigen.


### LIMITATION OF THE PROCEDURE

1. The "Procedure" and "Interpretation of Results" sections must be followed closely to assure the accuracy of the test results.
2. No valid conclusion concerning the test result can be obtained, if controls with unclear or false results should occur.
3. Enzyme treated erythrocytes or addition of bovine albumin and/or other solutions containing protein may cause unspecific reactions.
4. Hemolyzed, turbid, contaminated or clotted samples should not be used.
5. Due to variability of free H-substance on human red blood cells, reactivity of the reagent against certain phenotypes may give weaker reactivity compared to control cells.
6. No one specific antiserum or technique can be guaranteed to detect all rare, weak or variant Antigens.<sup>2</sup>
7. Red blood cells coated with alloantibodies or autoantibodies of the same or similar specificity as that reagent used for the test (i.e., cells that are positive in the direct antiglobulin-test (DAT)) are not suitable for this test procedure.

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## SYMBOL - LEGENDE

 Store from - to	<b>REF</b>	Product Code	 Expiration Date	 Manufacturer according to 98/79/EU			
<b>LOT</b>	Lot	<b>CLON</b>	Clone(s)	<b>IVD</b>	In vitro diagnostic medical device	<b>CE</b>	EU CE-symbol
<b>UDI</b>	Unique Device Identification		Consult Instrucion for use				

730-22-2108 Version 008 / 15.8.2021

 Antitoxin GmbH, Industriestraße 88, 69245 Bammental, Germany

 +49 (0) 6223/ 8661-0  +49 (0) 6223/ 8661-13  [gara@antitoxin-gmbh.de](mailto:gara@antitoxin-gmbh.de)

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### Anti-H (A2), monoclonal IgM (mouse) OrthoClone

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Anti-H monoclonale, mouse IgM clone: 10934C11

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Antitoxin GmbH, Industriestraße 88, 69245 Bammental, Germania