Anti-K (KEL1) monoclonal, human IgM clone: MS-56

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Anti-K (KEL1) monoclonal, human IgM clone: MS-56

English

For Slide-, Tube-, Card- and Micoplate-Method
FOR IN VITRO DIAGNOSTIC USE ONLY

INTENDED USE

The reagent is used to in-vitro determine qualitative whether human red blood cells possess or lack the corresponding blood group antigen K.
The reagent is intended to be used by qualified and technical personnel only to perform immuno-hematology screening tests as part of the practice of transfusion medicine in the general population.
The test method used with this reagent is based on the principle of agglutination, performed not automated.
Normal human erythrocytes, possessing the corresponding antigen, will agglutinate in the presence of the specific antibody directed toward the antigen.

INDICATION / CONTRA-INDICATION

The Anti-K monoclonal blood grouping reagent is used to test patient or donor red cells for the presence of the K antigen. Typing of donor cells facilitates the selection of suitable antigen-negative units for transfusion to patients with this antibody. Cell typing also serves as final verification of the identification of Anti-K in patient or donor sera.
There is no contra-indication to perform the in-vitro-test on blood samples.
The product was validated with sample collected in Europe from patients of unknown ethic background.

The approximate frequencies of K antigen:

Phenotype	European	African
K+k-	0,2%	seldom
K+k+	8,8%	2%
K-k+	91%	98%

REAGENTS

The listed reagent contains antibodies of the following cell clone:

Anti-K (KEL1) monoclonal, human IgM clone: MS-56

The reagent is produced from cell culture supernatant of hetero-hybridoma-cell lines, which secreting antibodies from IgM-Type that reacts specifically with the corresponding antigen.
The antibodies are human protein.
The reagent contains <0.1% (w/v) sodium azide as preservative.
In addition to the active antibody component the test serum contains, sodium chloride, macromolecules, bovine albumin, which has been tested negative for Vesicular Stomatitis Virus and Blue-tongue.
The BSA was derived from US sourced animals from USDA and APHIS approved facilities, to be used for in vitro diagnostic reagents, in compliance with EC 1069/2009 and EU 142/2011.

WARNING

This reagent is prepared from supernatants of cell cultures. Nevertheless, as biological product, it should be looked upon as potentially infectious because of never complete exclusion of danger through excitants of disease. The reagent contains sodium azide that may be toxic and may react with lead or copper to form highly explosive salts.
On disposal, flush with large quantities of water.
For the reasons mentioned above the reagent should be handled with proper care.

STORAGE REQUIREMENT

Store opened and unopened products at +2 to +8°C. May be at room temperature while in use.
For the stimulation of an In Use Study, the sera were stored 30 times at room temperature for 2 hours and subsequently showed no differences in the qualitative tests until the expiry date.
In principle, store and use the reagent to declared expiry date only.
(Format for the expiry date: Year xxxx Month xx Day xx)

REMARKS

- With each testing positive and negative controls should be performed.
- Inappropriate storage impairs efficacy of the reagent.
- Weak turbidity of the reagent does not affect its reactivity. Bacteria and chemical contamination of the product should be avoided. If a visible change (increase in turbidity or colour change due to temperature) is detected, the reagent should no longer be used, this sign may indicate a microbiological contamination.
- Do not use leaking, unlabeled or broken vials.
- Strength of positive reactions also depends on age of used blood.
- Centrifugation outside the specified speed range may lead to false results.
The use of another card-specific centrifuge (each card centrifuge has its specified unchangeable g-force) may lead to false results due to the changed g-force.
- The test method identified below is for manual testing only and must be carried out according to the instructions for use.
 - In case of changes in technology / deviations from the Instructions for Use
 - use of automated or semi-automated systemsthe laboratories have to follow the procedures that are contained in the operator's manual provided by the device manufacturer and carry out validations according to recognized procedures.
- For usage of the reagent all effective national laws, directives and guidelines have to be observed in its actual form, in Germany especially the „Richtlinien zur Gewinnung von Blut und Blutbestandteilen und zur Anwendung von Blutprodukten (Hämotherapie)“
- The information on the use of the test cards (a plastic frame with 6 or 8 microtubes filled with a buffered medium) in the relevant insert must be observed.
- The information on the use of the different materials in the relevant insert must be observed.
- This reagent has been validated and approved with the Slide Method, Tube Centrifugation Method, the Card Method and the Microplate Method.

SAMPLE PREPARATION

- Blood sample should be obtained using an acceptable phlebotomy technique.
Sample drawn into tubes containing EDTA, Sodium Citrate, CPD-A, ACD or Blood bag containing PAGGS-M may be used.
- Blood samples to be tested should be used immediately after blood collection to reduce the risk of false-positive and false-negative results due to improper storage or contamination of the samples.
Samples that cannot be tested immediately should be stored at +2 to +8°C.
Blood drawn into EDTA should be tested within 7 days and samples treated with sodium citrate, CPD-A and ACD within 14 days after collection.
Blood bag / Donor Blood (with PAGGS-M) can be tested until the expiry date.
- Do not freeze samples.

REAGENT PREPARATION

There is no preparation of the reagent required.
Take and use the reagent directly from the vials.

ADDITIONAL REQUIRED MATERIAL NOT SUPPLIED AND THE ASSOCIATED PROCEDURAL INSTRUCTION

Slide Method

Materials:

- Glass slide
- Pastdeur pipette
- Mixing stick
- Timer

Procedure:

- Use erythrocyte sediment only.
- Place one drop (approximately 50 µL) of the reagent on a marked glass slide.
- Using a Pasteur pipette add one drop erythrocyte sediment (approximately 50 µL) to the drop of reagent on the glass slide.
- Mix the erythrocytes with reagent well with a stick and spread to a circle with a diameter of approximately 2 cm.
- By slightly rotating the slide, check for agglutination within 1 minute (reaction starts within seconds).
- Document the result.

Tube Method:

Materials:

- Tubes (10 x 75 mm or 12 x 75 mm)
- Microliter pipette
- Centrifuge
- Timer
- isotonic saline (0,85 - 0,9% sodium chlorid)

Procedure:

- Prepare 2 to 5 % suspensions of red blood cells in isotonic saline (red blood cells may be washed 1 - 3 times with isotonic saline).
- At first put 100 µL of appropriate reagent in a marked tube subsequently add 100 µL of appropriate cell suspension in the tube.
Alternative one drop = approximately 50 µL cell suspension can be added one drop = approximately 50 µL test serum.
- Mix Erythrocytes-/Reagent mixture well by slightly shaking.
- Incubate tube at room temperature for 15 min.
- Centrifugation of tube for 1 minute at 2.000 rpm (approximately 800-1.000 x g).
- Gently shake the red cells completely from the bottom of the tube and check macroscopically for agglutination within 3 minutes.
- Document the result.

Micoplate-Method:

Materials:

- Microplate with 96 U-wells
- Tubes (10 x 75 mm or 12 x 75 mm)
- Microliter pipette
- Centrifuge
- Timer
- Microplate shaker
- Microplate centrifuge
- Isotonic saline (0,85 - 0,9% sodium chlorid)

Procedure:

- Prepare 2 to 5 % suspensions of red blood cells in isotonic saline (red blood cells may be washed 1 - 3 times with isotonic saline).
- Add 50 µL of appropriate reagent into a marked well.
- Add 50 µL of appropriate cell suspension to the well.
- Shake the red cell/test serum mixture in the microtitre plate on a microtitre plate shaker for 30 seconds on medium speed.
- Centrifugation of the microplate in appropriate microplate centrifuge for 30 seconds at 400 x g.
- Shake microplate for 30 seconds on a microplate shaker with medium speed.
- Examine the test results macroscopically for agglutinaion immediately after shaking.
- Document the result.
- Tests with negative or doubtful results have to be incubated for 5 to 10 minutes at room temperature.
- Repeat steps 5 to 8 after incubation.

Card Method:

Materials:

Card Method manually Grifols:

- Card: Grifols „DG Gel Neutral“ REF 210343
- Tubes (10 x 75 mm or 12 x 75 mm)
- Microliter pipette
- Isotonic saline (0,85 - 0,9% sodium chlorid)
- Centrifuge
- Card specific diluent "DG Gel Sol" REF 210354
- Grifols Card Centrifuge "DG-Spin"

Procedure: Neutral Card

- Prepare 0,8 % suspensions of red blood cells in "DG Gel Sol" (card specific diluent). (red blood cells may be washed 1 - 3 times with isotonic saline).
- Add 50 µL of appropriate cell suspension to each labeled microtube.
- Add 25 µL of the reagent to each microtube.
- No Incubation- / Reaction time for the "DG Gel Neutral" card
- Centrifugation in Grifols card centrifuge with the, for the respective centrifuge, unchangeable g-force.
- Check macroscopically for agglutination within 30 minutes.
- Document the result.

Materials:

Card Method manually BIO-RAD (DiaMed):

- 1. Card: Bio-Rad "NaCl, Enzyme test and cold Agglutinins " REF 005014
- 2. Tubes (10 x 75 mm or 12 x 75 mm)
- 3. Microliter pipette
- 4. Isotonic saline (0,85 - 0,9% sodium chlorid)
- 5. Centrifuge
- 6. Card specific diluent: "ID-Diluent 2" REF 009280
- 7. Bio-Rad Card Centrifuge: ID-Centrifuge 24S

Procedure: „NaCl, Enzyme Test and cold Agglutinins card“

- 1. Prepare 0,8 % suspensions of red blood cells in "ID-Diluent 2" (card specific diluent). (red blood cells may be washed 1 - 3 times with isotonic saline).
- 2. Add 50 µL of appropriate cell suspension to each labeled microtube.
- 3. Add 25 µL of the reagent to each microtube.
- 4. No Incubation- / Reaction time for the „NaCl, Enzyme Test and cold Agglutinins" card.
- 5. Centrifugation in Bio-Rad card centrifuge with the, for the respective centrifuge, unchangeable g-force.
- 6. Check macroscopically for agglutination within 30 minutes.
- 7. Document the result.

Materials:

Card Method manually BioVue® System:

- 1. Card: Ortho BioVue® "Reverse Diluent Cassette" REF 707550
- 2. Tubes (10 x 75 mm or 12 x 75 mm)
- 3. Microliter pipette
- 4. Isotonic saline (0,85 - 0,9% sodium chlorid)
- 5. Centrifuge
- 6. Card specific diluent: Isotonic saline (0,85 - 0,9% NaCl)
- 7. Ortho Card centrifuge: „BioVue Workstation"

Procedure: Reverse Card

- 1. Prepare 3% to 5% suspensions of red blood cells in isotonic saline (red blood cells may be washed 1-3 times with isotonic saline)
- 2. Add 40 µL of the reagent to each labeled microtube.
- 3. Add 10 µL of appropriate cell suspension to each microtube.
- 4. No Incubation- / Reaction time for the "Ortho BioVue® Reverse Diluent Cassette".
- 5. Centrifugation in Ortho BioVue® card centrifugation with the, for the respective centrifuge, unchangeable g-force.
- 6. Check macroscopically for agglutination immediately after the end of centrifugation.
- 7. Document the result.

INTERPRETATION OF RESULTS

Positive results (+): If agglutination of the erythrocytes occurs within the accepted limitations of the test procedure, the test result is positive and indicates the presence of the corresponding antigen.

Negative results (-): If agglutination of the erythrocytes does not occur within the accepted limitations of the test procedure, the test result is negative and indicates the absence of the corresponding antigen.

The reading and interpretation of the results after "careful shaking" at the Slide, Tube Centrifuga-tion Method / Microplate Method:

Negative	No detectable agglutinates, homogeneous red coloration of the liquid.
Positive	One complete agglutinate.
	No complete agglutinate, some individual agglutinates
	Red colouration of the liquid containing only small / miniature agglutinates.

Read and interpret the results of the card methods in accordance with the respective card instruc-tions.

Negative	All erythrocytes are through the column. A strip of red cells at the bottom of the column and no visible agglutinated cells in the rest of the column.
Positive	Few small agglutinations in the lower half of the column and a strip of cells at the bottom of the column.
	Some small agglutinations in the lower half of the column.
	Small or medium-sized agglutinations distributed throughout the column.
	Medium-sized agglutinations in the upper half of the column.
	Strip / band of agglutinated erythrocytes in the upper area / at the upper end of the column.

LIMITATIONS OF THE PROCEDURE

- 1. Inaccuracy at compliance with instructions written under section "Procedures" and "Interpretation of results" may lead to incorrect results.
- 2. No valid conclusion concerning the test result can be reached, if controls with uncertain or false results occur.
- 3. Enzyme treated erythrocytes or addition of bovine albumin and / or other solutions containing protein may cause unspecific reactions.
- 4. Hemolyzed, turbid, contaminated or clotted blood samples must not be used in this test.
- 5. A red blood cell suspension with a concentration that deviates from the indicated concentration may lead to false positive or false negative results.
- 6. The use of a diluent other than the one indicated on the individual cards for the erythrocyte suspensions may lead to an altered reaction behaviour.
- 7. The addition of volumes that deviate from the volumes specified in the method may lead to altered reaction behavior.
- 8. With the slide method, unspecific reactions might appear due to drying of the reaction-formation or if the slide is heated.
- 9. Due to variability of antigen expression on human red blood cells, reactivity of the reagent mentioned above, against certain phenotypes, may give weaker reactivity compared to control cells.
- 10. No one specific antiserum or technique can be guaranteed to detect all rare, weak or variants of antigens.²
- 11. Red blood cells with a strong positive direct Coombs-test may give false positive results in rare cases.
- 12. Pay attention to all statements to limitations in the instruction for use of the cards.

INCIDENTS RELATED TO THE DEVICE

Any serious incident that has occurred in relation to the device must be reported to the manufactur-er and the competent authority of the Member State where the user and/or patient is established.

PERFORMANCE CHARACTERISTICS

A performance evaluation for the products was carried out in accordance with the Common Specification (CS Common Specification of 04. July 2022). The required samples were used and compared with other reference methods / products.

Technique \ Product	Anti-K (KEL1) monoclonal, human IgM clone: MS-56			
	Positive Bleeds n	Sensitivity	Negative Bleeds n	Specificity
Slide methode	7 / 7	100 %	123 / 123	100 %
Tube methode	15 / 15	100 %	110 / 110	100 %
Mikroplate methode	15 / 15	100 %	110 / 110	100 %
Grifols „DG Gel Neutral" manually	7 / 7	100 %	123 / 123	100 %
Bio-Rad "NaCl, Enzyme test and cold Agglutinins" manually	7 / 7	100 %	123 / 123	100 %

Diagnostic Sensitivity: The probability that the device gives a positive result in the presence of the target marker.

Diagnostic Specificity: The probability that the device gives a negative result in the absence of the target marker.

Anti-K (KEL1) monoclonal human IgM clone: MS-56 is equivalent and does not differ in quality from comparable reagents available on the market.

DIFFERENCES BETWEEN BATCHES

Validation between three batches over the entire shelf life showed no differences.

INTERFERENCE STUDY

The interference studies showed no impairment for the qualitative test when using the following interfering substances:
Heparin 720 U/dl, Albumin 15000 mg/dl, Triglycerides 1500 mg/dl, Bilirubin 40mg/dl, Ethanol 620 mg/dl, Glucose 1000 mg/dl.
For the anticoagulants (EDTA, Sodium Citrate, ACD, CPD-A, PAGGS-M) three times the recom-mended concentration was tested.

SUMMARY OF SAFETY AND PERFORMANCE

The Summary of Safety and Performance of this reagent is available via ANTITOXIN (www.antitoxin-gmbh.de) and can be accessed via the EUDAMED database.

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REF	Product code	LOT	Batch
	Store from - to		Experation Date
IVD	In-Vitro Diagnostic	CE	EU CE symbol
	Manufacture according to (EU) 2017/746		Consult instruction for use
UDI	Unique Device Identification		Distributor

REF

- 01.047-10 Anti-K (KEL1) monoclonal, human IgM clone: MS-56 10ml
- 01.047-10.X Anti-K (KEL1) monoclonal, human IgM clone: MS-56 10 x 10ml

CE 0483

 ANTITOXIN GmbH Industriestraße 88 69245 Bammental Germany

 +49 (0) 6223/ 8661-0

 +49 (0) 6223/ 8661-13

 gara@antitoxin-gmbh.de

 01.047- / Version R002 / 2025-04-22

Marking of changes:

Changes to the previous version are shaded grey.

BAG Diagnostics GmbH
Amtsgerichtsstr. 1-5
35423 Lich/Germany
A BAG Group company

T. +49 (0) 6404/925-100
F. +49 (0) 6404/925-460
M. info@bag-diagnostics.com
W. www.bag-diagnostics.com

Ordering: order@bag-diagnostics.com
T. +49 (0) 6404 / 925 – 450
F. +49 (0) 6404 / 925 - 460
M. order@bag-diagnostics.com

Customer Service:
T. +49 (0) 6404 / 925 - 125
F. +49 (0) 6404 / 925 - 421
M. service@bag-diagnostics.com