Anti-D monoclonal, human IgM

Deutsch

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Anti-D monoclonal, human IgM

English

For Slide-, Tube-, Card- and Microplate - Method
FOR IN VITRO DIAGNOSTIC USE ONLY

INTENDED USE

Monoclonal agglutinating Anti-D reagents are produced from cell culture supernatants of heterohybridoma - cell lines. The cells are secreting an antibody of IgM-type that reacts specific with the corresponding blood group antigen. The antibody in each case is human protein. The reagents are used for In-Vitro-Diagnostic, to determine qualitative whether human red blood cells possess or lack the corresponding blood group antigen D. These reagents react negative with erythrocytes of category D^{VI}. The reagents are intended to be used by qualified and technical personnel only.

PRINCIPLE OF PROCEDURE

The test methods used with these reagents are based on the principle of agglutination. Normal human erythrocytes, possessing the corresponding antigen, will agglutinate in the presence of the specific antibody directed towards the antigen.

REAGENT

The listed reagents contain antibody of the following cell clones:

Anti-D monoclonal, human IgM clone: RUM-1

Anti-D monoclonal, human IgM clone: P3x61

These reagents contain <0.1% (w/v) sodium azide as preservative.

Additionally the reagents are prepared of active antibody, sodium chloride, macromolecules and bovine albumin, which has been tested and certified by the US Veterinary service inspectors.

WARNING

These reagents were prepared from supernatants of cell cultures. As biological products they should be looked upon as potentially infectious because of never complete exclusion of danger through excipients of disease. The reagents contain sodium azide that may be toxic and may react with lead or copper to form highly explosive salts.

On disposal, flush with large quantities of water.
For the reasons mentioned above, reagents should be handled with proper care.

STORAGE REQUIREMENT

Store opened and unopened products at +2 to +8°C. May be at room temperature while in use. In principle, store and use the reagents to indicated expiry date only.

REMARKS

1. With each testing positive and negative controls should be performed.
2. Inappropriate storage impairs efficacy of the reagents.
3. Weak turbidity of the reagents mentioned above does not affect its reactivity. Bacteria and chemical contamination of the products should be avoided. If a visible change is detected, the reagent should no longer be used, this sign may indicate a microbiological contamination.
4. Strength of positive reactions also depends on age of used blood
5. Centrifugation outside the specified speed range may lead to false results.
With the card method, the use of another card-specific centrifuge (each card centrifuge has its specified unchangeable g-force) may lead to false results due to the changed g-force.
6. The test methods identified below are for manual testing only. When using automated or semi-automated instruments, follow the procedures that are contained in the operator's manual provided by the device manufacturer. Laboratories must follow approved validation procedures.
7. For usage of these reagents all effective national laws, directives and guidelines have to be observed, in Germany especially the „Richtlinien zur Gewinnung von Blut und Blutbestandteilen und zur Anwendung von Blutprodukten (Hämotherapie)“¹⁾ in its current version.
8. These reagents react negative with erythrocytes of category D^{VI}.
9. The information on the use of the test cards in the relevant insert must be observed.

SAMPLE PREPARATION

1. Blood sample should be collected by approved medical procedure.
2. Blood sample to be tested should be used as soon as possible after blood collection to reduce the risk of false-positive and false-negative results due to improper storage or contamination of the sample.
If a delay in testing occurs, sample should be stored at +2 to +8°C.
Blood drawn into EDTA should be tested within 7 days and sample treated with sodium citrate within 14 days after collection.
Blood bag / Donor Blood can be tested until the expiry date.

REAGENT PREPARATION

There is no preparation of the reagents required.
Take and use reagents directly from the vials.

PROCEDURE

Not provided material, additionally needed

Slide Method:	glass slide; Pasteur pipette; mixing stick; timer
Tube Centrifugation Method:	tubes (10 x 75 mm or 12 x 75 mm); microliter pipette; centrifuge; timer; isotonic saline (0,85 - 0,9% sodium chloride)
Card Method:	Cards: - BIO-RAD (DiaMed) „NaCl, Enzyme Test and cold Agglutinins“ - Grifols „DG Gel Neutral“; tubes; microliter pipette; centrifuge; timer; corresponding card centrifuge; card specific diluent; isotonic saline (0,85 - 0,9% sodium chloride)
Microplate Method:	microplate with 96 U-wells; microliter pipette; centrifuge; timer; microplate shaker; microplate centrifuge; isotonic saline (0,85 - 0,9% sodium chloride)

Test procedure

Slide Method

1. Use erythrocyte sediment only.
2. Place one drop (approximately 50 µL) of appropriate reagent on a marked glass slide.
3. Using a Pasteur pipette add one drop erythrocyte sediment (approximately 50 µL) to the drop of reagent on the glass slide.
4. Mix the erythrocytes with reagent well with a stick and spread to a circle with a diameter of approximately 2 cm.
5. By slightly rotating the slide, check for agglutination within 1 minute (reaction starts within seconds).
6. Document the result.

Tube Centrifugation Method

1. Prepare a 2% to 5% suspension of red blood cells in isotonic saline (red blood cells may be washed 1-3 times with isotonic saline).
2. At first put 100 µL of appropriate reagent in a marked tube subsequently add 100 µL of appropriate cell suspension in the tube.
Alternative one drop = approximately 50 µL cell suspension can be added to one drop = approximately 50 µL test serum.
3. Mix Erythrocytes-/Reagent mixture well by slightly shaking.
4. Incubate tube at room temperature for 1-15 min.
5. Centrifugation of tube for 1 minute at 2.000 rpm (approximately 800-1.000 x g).
6. Gently shake the red cells completely from the bottom of the tube and check macroscopically for agglutination within 3 minutes.
7. Document the result.

Card Method (manual method / valid for the cards:

- BIO-RAD (DiaMed) „NaCl, Enzyme Test and cold Agglutinins“
- Grifols „DG Gel Neutral“

1. Prepare 0,8 % suspension of red blood cells in card specific diluent (red blood cells may be washed 1-3 times with isotonic saline).
 2. Add 50 µL of appropriate cell suspension to a marked micro tube.
 3. Add 25 µL of appropriate reagent to the micro tube.
 4. Centrifuge the card in appropriate card centrifuge with the for this centrifuge unchangeable g-force within 30 minutes.
 5. Check macroscopically micro tube for agglutination within 30 minutes.
 6. Document the result.
- ### Microplate Method
1. Prepare 2-5% suspension of red blood cells in isotonic saline (red blood cells may be washed 1-3 times with isotonic saline).
 2. Add 50 µL of appropriate reagent to a marked well.
 3. Add 50 µL of appropriate cell suspension to the well.
 4. Shake the microplate for 30 seconds on microplate shaker with medium speed.
 5. Centrifuge the microplate in appropriate microplate centrifuge for 30 seconds at 400 x g.
 6. Shake microplate for 30 seconds on a microplate shaker with medium speed.
 7. Evaluate macroscopically on agglutination directly after shaking.
 8. Document the result.
 9. Tests with negative or doubtful results have to be incubated for 5 to 10 minutes at room temperature.
 10. Repeat steps 5 to 8 after incubation at room temperature

INTERPRETATION OF RESULTS

"Slightly rotating / shaking"

at Slide Method / at Tube Centrifugation Method / Microplate Method:

Positive results (+): visible agglutination of erythrocytes is a positive result and

indicates the presence of the corresponding antigen.

Negative results (-): No visible agglutination of erythrocytes is a negative result and

indicates the absence of the corresponding antigen.

Read and interpret the results of the card according to the instruction for use of the corresponding card.

LIMITATIONS OF THE PROCEDURE

1. Inaccuracy at compliance with instructions written under section "Procedures" and "Interpretation of results" may lead to incorrect results.
2. No valid conclusion concerning the test result can be reached, if controls with uncertain or false results occur.
3. Enzyme treated erythrocytes or addition of bovine albumin and/or other solutions containing protein may cause unspecific reactions.
4. Hemolyzed, turbid, contaminated or clotted samples should not be used.
5. With the slide method, unspecific reactions might appear due to drying of the reaction-formation or if the slide is heated.
6. With slide method most of weak D-antigens (D^{weak}) and categories are not recognized.
7. The reagents listed above react negatively with erythrocytes of category D^{VI}.
8. Due to variability of antigen expression, reactivity of these reagents against certain phenotypes may give weaker reactivity compared to control cells.
9. No one specific antiserum or technique can be guaranteed to detect all rare, weak or variant antigens.
10. Red blood cells coated with antibodies (cells that are positive in the direct antiglobulin-test (DAT)) may give false-positive results in card method. These cells react positive without testserum too.
11. Pay attention to all statements to limitations in the instruction for use of the cards.
12. Newborn blood is excluded for the clone P3x61 on the microtiter plate.
13. If products are taken out of 1 L flasks, the material of this bottle has to be used within 6 months.

PERFORMANCE CHARACTERISTICS

In compliance with Common Technical Specifications (CTS Commission Decision of 03. February 2009) a performance evaluation was conducted. Different samples (donor, patient, neonates-, panel blood) was used and compared with other reference methods / products.

Produkt	Positive	Negative	Falsch Positive	Falsch Negative
Anti- D (clone: Rum-1)	2449	592	0	0*
Anti- D (clone: P3x61)	2586	617	0	0*

The calculated values found for:

Produkt	Sensitivität	Spezifität
Anti- D (clone: Rum-1)	100%*	100%
Anti- D (clone: P3x61)	100%*	100%

* for Anti-D (clone RUM-1) in slide method 4 D^{weak} samples were tested false negative, so in this method sensitivity was 99,77% only and for Anti-D (clone P3x61) 8 D^{weak} samples were tested false negative in this method, so sensitivity was 99,56% only.

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
 Product Code	 Batch
 Store from - to	 Expiration Date
 In-Vitro Diagnostic	 EU CE symbol
 Manufacturer according to 98/79/EG	 Consult instructions for use
 Unique Device Identification	 Distributor

REF

01.011 - 10	Anti-D monoclonal, human IgM clone: RUM-1	10 ml
01.011 - 10.V	Anti-D monoclonal, human IgM clone: RUM-1	5 x 10 ml
01.011 - 10.X	Anti-D monoclonal, human IgM clone: RUM-1	10 x 10 ml
01.011 - 1K	Anti-D monoclonal, human IgM clone: RUM-1	1000 ml
01.167 - 10	Anti-D monoclonal, human IgM clone: P3X61	10 ml
01.167 - 10.V	Anti-D monoclonal, human IgM clone: P3X61	5 x 10 ml
01.167 - 10.X	Anti-D monoclonal, human IgM clone: P3X61	10 x 10 ml

730-15-0307 Version 007 / 01.09.2021

CE 0483

 ANTI TOXIN GmbH Industriestraße 88 69245 Bammental Germany
 +49 (0) 6223/ 8661-0  +49 (0) 6223/ 8661-13  qara@antitoxin-gmbh.de