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## INSTRUCTIONS FOR USE

**Anti-C (RH2), monoclonal IgM (human) OrthoClone**  
**Anti-E (RH3), monoclonal IgM (human) OrthoClone**  
**Anti-c (RH4), monoclonal IgM (human) OrthoClone**  
**Anti-e (RH5), monoclonal IgM (human) OrthoClone**

REF	690211A
REF	690321A
REF	690281A
REF	690361A

### INTENDED USE

Monoclonal Anti-C, -E, -c, -e reagents are produced from cell culture supernatants of heterohybridoma cell lines. The cells secrete an antibody of IgM-type that reacts specifically with the corresponding antigen. The antibody in each case is human protein. The reagents are used for In-Vitro-Diagnostic, to determine whether human red blood cells possess or lack the corresponding blood group antigens C, E, c and e. The procedures described are for the listed manual methods and application with Ortho BioVue® Reverse Diluent cassettes on the ORTHO VISION™ Analyzer and ORTHO AutoVue® Innova.

The test sera are intended to be used by qualified technical personnel only.

### PRINCIPLE OF PROCEDURE

The procedures used with these reagents are based on the principle of agglutination. Normal human erythrocytes, possessing the corresponding antigen, will agglutinate in the presence of the specific antibody directed towards the antigen.

### REAGENT

The listed reagent contains antibodies of the following cell clones:

Anti-C monoclonal, human IgM Clone: MS-273

Anti-c monoclonal, human IgM Clone: MS-35

Anti-E monoclonal, human IgM Clones: MS-12, MS-260

Anti-e monoclonal, human IgM Clones: MS-62, MS-69

These reagents contain <0.1% (w/v) sodium azide as preservative. Additionally these reagents are comprised of active antibody, sodium chloride, macromolecules and bovine albumin, which has been tested and certified by the US Veterinary service inspectors.

**CAUTION:** Please handle biological reagents with proper care: These test serum are prepared from supernatants of cell cultures. As biological products they should be looked upon as potentially infectious because of never complete exclusion of danger through excitants of disease. The test serum contains sodium azide that may be toxic and may react with lead or copper to form highly explosive salts. On disposal, flush with large quantities of water.

Found from the above reasons the reagents should be handled with proper care.

### STORAGE REQUIREMENTS

Store at +2 to +8 °C (unopened / opened), or shortly at room temperature while in use. If the devices ORTHO VISION™ Analyzer or ORTHO AutoVue® Innova are used, the reagents may remain a maximum of 5 days on the device (+15 bis +30 °C). Do not use reagents beyond the stated expiration date!

### REMARKS

1. It is recommended that each lot of reagents be tested with appropriate positive and negative controls.
2. Inappropriate storage impairs efficacy of the reagent.
3. Weak turbidity of the reagent does not affect its reactivity. Bacteria and chemical contamination of the product should be avoided. If a visible change is detected, the reagent should no longer be used, this sign may indicate a microbiological contamination.
4. Strength of positive reactions also depends on age of used blood.
5. Centrifugation outside the specified speed range may lead to false results. With the card method, the use of another card-specific centrifuge (each card centrifuge has its specified unchangeable g-force) may lead to false results due to the changed g-force.
6. The procedures described below are for manual methods and listed analyzers only. When using other automated or semi-automated instruments, follow the device manufacturer's procedures in the operator's manual. Laboratories must follow approved validation procedures.
7. For usage of these test sera all effective national laws, directives and guidelines in its current version must be followed. [In Germany especially: „Richtlinien zur Gewinnung von Blut und Blutbestandteilen und zur Anwendung von Blutprodukten (Hämotherapie)“<sup>1</sup>].

### SAMPLE PREPARATION

1. Blood sample should be collected by approved medical procedure.
2. Blood sample to be tested should be used as soon as possible after blood collection to reduce the risk of false-positive and false-negative results due to improper storage or contamination of the sample.

If a delay in testing occurs, sample should be stored at +2 to +8°C.

Blood drawn into EDTA should be tested within 7 days and sample treated with sodium citrate within 14 days after collection.

Blood bag / Donor Blood can be tested until the expiry date.

### REAGENT PREPARATION

There is no specific preparation of the reagents required. Use reagents directly from the vial.

### PROCEDURE

Material required but not provided:

#### Slide Method

1. Glass slide
2. Pasteur pipette
3. Mixing stick
4. Timer

#### Tube Centrifugation Method

1. Test tubes, 10 x 75 mm or 12 x 75 mm
2. Mikroliter pipette designed to deliver approximately 50µL/100 µL
3. Timer
4. Centrifuge
5. Isotonic saline (0,85 - 0,9% sodium chloride)
6. Disposable pipette tips

#### Microplate Method

1. Microplate with U-bottom, (optional: pre-treated)
2. Timer
3. Centrifuge, suitable for centrifugation of microplates
4. Microplate-carrier for centrifuge (option)
5. Microplate shaker
6. Microplatetest mirror (option)
7. Pipette for 50 µL
8. Isotonic saline; (0.85 – 0.9 % sodium chloride)
9. Disposable pipette tips
10. Test tubes, 10 x 75 mm or 12 x 75 mm

#### Application of Ortho BioVue® Cassettes

1. Ortho BioVue® System Centrifuge
2. Isotonic saline (0,85 – 0,9 % sodium chloride)
3. Mikroliter pipette for 10 µL, 40 µL and 50 µL
4. Disposable pipette tips
5. Ortho BioVue® System Work Rack
6. Ortho BioVue® System Reverse
7. Timer
8. Test tubes, 10 x 75 mm or 12 x 75 mm
9. Centrifuge

### Test procedure

#### Slide Method

1. Use erythrocyte sediment only.
2. Place one drop (approximately 50 µL) of appropriate reagent on a marked glass slide.
3. Using a Pasteur pipette add one drop erythrocyte sediment (approximately 50 µL) to the drop of reagent on the glass slide.
4. Mix the erythrocytes with reagent well with a stick and spread to a circle with a diameter of approximately 2 cm.
5. By slightly rotating the slide, check for agglutination within 1 minute (reaction starts within seconds).
6. Document the result.

#### Tube Centrifugation Method

1. Prepare a 2% to 5% suspension of red blood cells in isotonic saline (red blood cells may be washed 1-3 times with isotonic saline).
2. Add 100 µL (alternative: one drop, approximately 50 µL) of appropriate reagent in a marked tube.
3. Add 100 µL (alternative: one drop, approximately 50 µL) of appropriate cell suspension to the tube.
4. Mix well by slightly shaking.
5. Incubate tube at room temperature for 1-15 min.
6. Centrifugation of tube for 1 minute at 2.000 rpm (approximately 800-1.000 g).
7. Gently shake the red cells completely from the bottom of the tube and check macroscopically for agglutination within 3 minutes.
8. Document the result.

#### Microplate Method

1. Prepare a 2% to 5% suspension of red blood cells in isotonic saline (red blood cells may be washed 1-3 times with isotonic saline).
2. Add one drop (approximately 50 µL) of appropriate reagent into a marked well.
3. Add one drop (50 µL) of cell suspension into the well.
4. Mix microplate on a microplate shaker.

**Note:** Recommended: (1) To mix: 10 to 30 seconds at medium strength. (2) To resuspend: 10 to 30 seconds at medium strength respectively with duration and strength, allowing a full resuspension of sediment without destruction of positive reactions.

5. Centrifugation of microplate. Recommended duration of centrifugation: 30 seconds at approximately 400 x g, or as long as necessary to produce best positive reactions with antigen positive cells, with used centrifuge and microplate, and still have an easy resuspension of antigen negative cells.

**Note:** The used centrifugal force should be the minimum force to produce a sediment and a clear supernatant. Because of multitude offered centrifuges a general speed or duration for best results can not be recommended. Each laboratory has to test and validate its equipment for speed and duration to obtain best results.

6. Resuspend the erythrocyte knobs on the shaker and position microplate for reading.
7. Read microplate and document the result.
8. Tests with negative or doubtful results have to be incubated for 5 to 10 minutes at room temperature.
9. Repeat steps 5 to 7 after incubation at room temperature.

**Note:** All tests should be resuspended and evaluated directly after centrifugation.

## Manual Test Method in Ortho BioVue® Reverse Diluent Cassettes

The procedure described below is for manual Ortho BioVue® Cassettes testing only. Please refer to the instructions for use of the Ortho BioVue® Reverse Diluent Cassettes, sections STORAGE REQUIREMENT, SAMPLING AND PREPARATION, PREPARATION OF THE TEST SERIES, and INTERPRETATION OF RESULTS, as well as -where applicable- LIMITATIONS OF THE PROCEDURE.

- For preparation of the cassettes, please refer to the instructions for use for Ortho BioVue® Reverse Cassettes, chapter „procedure“ point 1 and 2.
  - Prepare a 3 - 5% erythrocyte suspension in isotonic sodium chloride solution. (red blood cells may be washed 1-3 times with isotonic saline)
  - Add 10 µL of a appropriate 3 to 5% red blood cell suspension to the marked reaction chamber of the cassette.
  - Add 40 µL appropriate reagent to the appropriate reaction chambers of the cassette.  
(The pipette tips shall not contact the rim of the reaction chambers. If this should adversely happen, please replace the tip with a new one, before you will fill the next reaction chamber.)
  - Centrifuge the cassette in the Ortho BioVue® System centrifuge with the for this centrifuge unchangeable g-force.
- Note:** Centrifugation must be performed within 30 minutes after addition of specimens to the reaction chambers.
- After completion of test, read each column from the cassette's front and rear sides for agglutination and/or hemolysis.
  - Protocol the final test result from the cassette side with the more clearly expressed positive reading.

## Test in Ortho BioVue® Reverse Diluent Cassettes on ORTHO VISION™ Analyzer or ORTHO AutoVue® Innova

Please follow the instructions in the user manual of the respective analyzer.

### INTERPRETATION OF RESULTS

"Slightly rotating / shaking" at Slide-, Tube- and Microplate Method

- Positive result (+):** visible agglutination of erythrocytes indicates the presence of the corresponding antigen.  
**Negative result (-):** no visible agglutination of erythrocytes indicates the absence of the corresponding antigens.

**Cassettes:** As for the manual method, for interpretation of results, please refer to the instructions for use of Ortho BioVue® Reverse Cassettes, chapter „INTERPRETATION OF RESULTS“. As for analyzers, reading and interpretation is effected automatically by each device.

### LIMITATIONS OF THE PROCEDURE

- The "Procedure" and "Interpretation of Results" sections must be followed closely to assure the accuracy of the test results.
- No valid conclusion concerning the test result can be reached, if controls with uncertain or false results occur.
- Enzyme treated erythrocytes or addition of bovine albumin and/or other solutions containing protein may cause unspecific reactions.
- Hemolyzed, turbid, contaminated or clotted samples should not be used.
- With the slide method, unspecific reactions might appear due to drying of the reaction-formation or if the slide is heated.
- Due to variability of antigen expression, reactivity of these reagents against certain phenotypes may produce a weaker reaction compared to control cells.
- Some red blood cells may have developed Rh variants that may not be recognized with monoclonal antibodies of a clone. Therefore the tests should / have to be carried out with at least two different clones / cloning mixtures as required in some countries. Deviating test results require further testing.
- No one specific antiserum or technique can be guaranteed to detect all rare, weak or variant antigens <sup>2</sup>.
- Red blood cells coated with alloantibodies or autoantibodies of the same or similar specificity as the reagent (i.e., cells that are positive in the direct antiglobulin-test [DAT]) may give false-positive results in card method. These cells react positive without testserum too.
- Red blood cells coated with antibodies (cells that are positive in the direct antiglobulin-test [DAT]) are not suitable for this test procedure.
- Pay attention to all statements to limitations in inserts of used Ortho BioVue® Reverse Diluent Cassettes and the user manual of respective analyzer.
- Negative test results with anti-e (clones: MS-62, MS-69) in microplate method have to be confirmed by another technique, as false negative results may occur.
- Do not use blood of newborn on microplates.
- With cord cells anti-c (clone: MS-33) may possess anti-I/i cold agglutinin activity and therefore show false positive results.

### PERFORMANCE

In compliance with Common Technical Specifications (CTS Commission Decision of 03. February 2009) a performance evaluation was conducted. Different samples (donor, patient, neonates-, panel blood) was used and compared with other reference methods / products.

Product	Positive	False negative	Sensitivity	Negative	False positive	Specificity
Anti-C monoclonal, human IgM Clone: MS-273	844	0	100%	391	0	100%
Anti-c monoclonal, human IgM Clone: MS-35	1005	0	100%	230	0	100%
Anti-E monoclonal, human IgM Clones: MS-12, MS-260	327	0	100%	908	0	100%
Anti-e monoclonal, human IgM Clones: MS-62, MS-69	1417	0*	100%*	65	0	100%

\* for Anti-e (clones MS-62, MS-69) in microplate method 42 samples were tested false negative, so in this method sensitivity was 97,12% only.  
(see "Limitations of the procedure" item 8)

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### LEGEND OF SYMBOLS

	Store from + to	<b>REF</b>	Product Code		Expiration Date		Manufacturer as to 98/79/EU
<b>LOT</b>	Lot	<b>CLON</b>	Clone(s)	<b>IVD</b>	In vitro diagnostic medical device		EU CE-symbol
<b>UDI</b>	Unique Device Identification		Consult Instrucion for use				

730-22-0518 Versione 018 / 01.07.2021

Antitoxin GmbH, Industriestraße 88, 69245 Bammental, Germany

+49 (0) 6223/ 8661-0 +49 (0) 6223/ 8661-13 gara@antitoxin-gmbh.de

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**Anti-C (RH2), monoclonal IgM (human) OrthoClone**  
**Anti-E (RH3), monoclonal IgM (human) OrthoClone**  
**Anti-c (RH4), monoclonal IgM (human) OrthoClone**  
**Anti-e (RH5), monoclonal IgM (human) OrthoClone**

REF	690211A
REF	690321A
REF	690281A
REF	690361A

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Anti-C monoclonale, human IgM Clone: MS-273

Anti-c monoclonale, human IgM Clone: MS-35

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Antitoxin GmbH, Industriestraße 88, 69245 Bammental, Niemcy