NUR ZUR IN-VITRO-DIAGNOSTIK

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**REF**

- 01.212-02** Anti-Jk<sup>a</sup> for card method monoclonal, human IgM clone: MS-15 2 ml  
**01.212-03** Anti-Jk<sup>a</sup> for card method monoclonal, human IgM clone: MS-15 3 ml  
**01.212-05** Anti-Jk<sup>a</sup> for card method monoclonal, human IgM clone: MS-15 5 ml  
**01.213-02** Anti-Jk<sup>b</sup> for card method monoclonal, human IgM clone: MS-8 2 ml  
**01.213-03** Anti-Jk<sup>b</sup> for card method monoclonal, human IgM clone: MS-8 3 ml  
**01.213-05** Anti-Jk<sup>b</sup> for card method monoclonal, human IgM clone: MS-8 5 ml

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FOR IN VITRO DIAGNOSTIC USE ONLY

**INTENDED USE**

Monoclonal agglutinating Anti-Jk<sup>a</sup> and Anti-Jk<sup>b</sup> reagents are produced from cell culture supernatants of heterohybridoma-cell lines. The cells are secreting an antibody of IgM-type that reacts specific with the corresponding antigen. The antibodies are human protein. The reagents are used to in-vitro determine qualitative whether red blood cells possess or lack the corresponding blood group antigens Jk<sup>a</sup> and Jk<sup>b</sup>. The reagents are intended to be used by qualified and technical personnel only.

**PRINCIPLE OF PROCEDURE**

The test methods used with these reagents are based on the principle of column agglutination. Normal human erythrocytes, possessing one of these antigens, will agglutinate in the presence of the specific antibody directed toward the antigen.

**REAGENTS**

The listed reagents contain antibodies of the following cell clones:

Anti-Jk<sup>a</sup> for card method monoclonal, human IgM clone: MS-15

Anti-Jk<sup>b</sup> for card method monoclonal, human IgM clone: MS-8

The reagents contain <0.1% (w/v) sodium azide as preservative. Additionally the reagents are prepared of active antibody, sodium chloride, macromolecules and bovine albumin, which has been tested and certified by the US Veterinary service inspectors.

**WARNING**

These reagents were prepared from supernatant of cell cultures. As biological products it should be looked upon as potentially infectious because of never complete exclusion of danger through excitants of disease. The reagents contain sodium azide that may be toxic and may react with lead or copper to form highly explosive salts. On disposal, flush with large quantities of water. For the reasons mentioned above the reagents should be handled with proper care.

**STORAGE REQUIREMENT**

Store opened and unopened products at +2 to +8°C. May be at room temperature while in use. In principle, store and use the reagents to declared expiry date only.

**REMARKS**

1. With each testing positive and negative controls should be performed.
2. Inappropriate storage impairs efficacy of the reagents.
3. Weak turbidity of the reagents does not affect its reactivity. Bacteria and chemical contamination of the product should be avoided. If a visible change is detected, the reagents should no longer be used, this sign may indicate a microbiological contamination.
4. Strength of positive reactions also depends on age of used blood.
5. Use the appropriate card centrifuge. The use of another card-specific centrifuge (each card centrifuge has its specified unchangeable g-force) may lead to false results due to the changed g-force.
6. The test methods identified below are for manual testing only and for the devices listed in this insert. When using other automated or semi-automated instruments, follow the procedures that are contained in the operator's manual provided by the device manufacturer. Laboratories must follow approved validation procedures to demonstrate compatibility of this product on such automated systems.
7. For usage of the reagents all effective national laws, directives and guidelines have to be observed, in Germany especially the „Richtlinien zur Gewinnung von Blut und Blutbestandteilen und zur Anwendung von Blutprodukten (Hämotherapie)“ in its actual form.
8. The information on the use of the test cards in the relevant insert must be observed.

**SAMPLE PREPARATION**

1. Blood sample should be collected by approved medical procedure.
2. Blood samples to be tested should be used as soon as possible after blood collection to reduce the risk of false-positive and false-negative results due to improper storage or contamination of the samples.  
If a delay in testing occurs, samples should be stored at +2 to +8°C. Blood drawn into EDTA should be tested within 7 days and samples treated with sodium citrate within 14 days after collection.  
Blood bag / Donor Blood can be tested until the expiry date.

**REAGENT PREPARATION**

There is no preparation of the reagents required. Take and use the reagents directly from the vials.

**PROCEDURE**

Not included but necessary materials:

**Card Method manuell:**

1. Cards: - Grifols „DG Gel Neutral“  
- Grifols „DG Gel Coombs“  
- Ortho BioVue® System Reverse cassettes
2. Microliter pipette
3. Glass tubes
4. Isotonic saline (0,85 - 0,9% sodium chloride)
5. Centrifuge
6. Card Centrifuge
7. Card specific diluent

**Card Method automatic:**

Vision® Analyzer

**Test procedure**

**Card Method** (manual method / valid for the cards:

- Grifols „DG Gel Neutral“
- Grifols „DG Gel Coombs“)

1. Prepare 0,8 % suspension of red blood cells in card specific diluent (red blood cells may be washed 1 - 3 times with isotonic saline).
2. Add 50 µL of appropriate cell suspension to each marked micro tube.
3. Add 25 µL of appropriate reagent to each micro tube.
4. Centrifuge the card in appropriate card centrifuge with the for this centrifuge unchangeable g-force.
5. Check macroscopically for agglutination within 30 minutes.
6. Document the result.

**Ortho BioVue® System** (manual method)

1. Prepare 3-5 % suspension of red blood cells in card specific diluent (red blood cells may be washed 1 - 3 times with isotonic saline)
2. Add 40 µL of appropriate reagent to the appropriate marked reaction chamber(s)
3. Add 10 µL of appropriate cell suspension to the appropriate reaction chamber(s)
4. Centrifuge the cassette in appropriate centrifuge with the for this centrifuge unchangeable g-force
5. The test results should be read directly after the end of centrifugation
6. Document the result

**Ortho BioVue® System** (Ortho Vision® Analyzer)

The use of the mentioned above reagents is validated on the Ortho Vision® Analyzer

The Ortho Vision® Analyzer offers the following options / settings:

- 50 µL 0,8% suspension of red blood cells in PBS buffer or
- 10 µL 3-5% suspension of red blood cells in PBS buffer

**INTERPRETATION OF RESULTS**

Positive results (+): visible agglutination of erythrocytes is a positive result and indicates the presence of the corresponding antigen.

Negative results (-): No visible agglutination of erythrocytes is a negative result and indicates the absence of the corresponding antigen.

Read and interpret the results of the card method according to the instruction for use of the corresponding card.

The interpretation of the results from Ortho Vision® Analyzer is performed automatically by the software.

**LIMITATIONS OF THE PROCEDURE**

1. Inaccuracy at compliance with instructions written under section "Procedures" and "Interpretation of results" may lead to incorrect results.
2. No valid conclusion concerning the test result can be reached, if controls with uncertain or false results occur.
3. Enzyme treated erythrocytes or addition of bovine albumin and/or other solutions containing protein may cause unspecific reactions.
4. Hemolyzed, turbid, contaminated or clotted blood samples must not be used in this test.
5. Due to variability of antigen expression on human red blood cells, reactivity of the reagents, mentioned above, against certain phenotypes may give weaker reactivity compared to control cells.
6. No one specific antiserum or technique can be guaranteed to detect all rare, weak or variant antigens.<sup>2</sup>
7. Red blood cells coated with alloantibodies or autoantibodies of the same or similar specificity, as that reagent used for the test (i.e. cells that are positive in the direct antiglobulintest (DAT)) are not suitable for this test procedure.
8. Red blood cells with a positive direct Coombs-test may cause false-positive reactions in card method.
9. It is described in the literature, that samples from patients treated with anti-CD38 monoclonal antibodies can cause false positive results in the Coombs test.<sup>5</sup>
10. Information on limits in the instructions for use of cards used must be observed.
11. Incubation at room temperature may increase potentially weak reactions on "Grifols DG Gel Neutral card".
12. For donors, a small pellet of non-agglutinated red blood cells in the bottom of the microtube may be observed in positive reactions (image of double population, DP). In patients, DPs could also be caused by prior transfusion.  
Additional information on patient history could be necessary for resolution.
13. When used on the Ortho Vision® Analyzer, the test serum may remain uninterruptedly on the device (+15 to +30 °C) for a maximum of 7 days. On the other hand, if the test serum is stored overnight at +2 to +8 °C for 12 hours, the period will be extended accordingly to 14 days.

**PERFORMANCE CHARACTERISTICS**

A performance evaluation for the products was conducted. Different samples material was used and compared with a reference method. When testing with the Grifols Neutral- and the Ortho Reverse card Blood bag with the additive solutions SAG-M und PAGGS-M, as well as EDTA- und Citrablood were used. When testing with the Grifols Coombs card Blood bag with the additive solution PAGGS-M, as well as EDTA- und Citrablood were used.

Anti-Jk<sup>a</sup> agreement of the results and performance data.

Method	Positive	Negative	Sensitivity	Specificity
Grifols Neutral Card	103 / 103	43 / 43	100%	100%
Ortho Reverse Card				
Grifols Coombs Card	54 / 54	20 / 20	100%	100%

Anti-Jk<sup>b</sup> agreement of the results and performance data:

Method	Positive	Negative	Sensitivity	Specificity
Grifols Neutral Card	106 / 106	36 / 36	100%	100%
Ortho Reverse Card				
Grifols Coombs Card	54 / 54	20 / 20	100%	100%

**LITERATURE**

1. Richtlinien zur Gewinnung von Blut und Blutbestandteilen und zur Anwendung von Blutprodukten (Hämotherapie)
2. CLSI, I/LA33-A Validatin of Automated System for Immunhematological Testing Before Implementation; Approved Guidline  
CLINICAL AND LABORATORY STANDARDS INSTITUTE December 2009
3. Peter D. Issit, David J. Anstee Applied Blood Group Serology, 4. edition, Montgomery Scientific Publications 1998
4. Christian Mueller-Eckhardt, Volker Kiefel Transfusionsmedizin, 3.edition, Springer-Verlag 2004.
5. Blood Transfusion Management for Patients Treated with Anti-CD38 monoclonal Antibodies.  
Frontiers in Immunology November 2018 / Article 2616

<b>REF</b> Product Code	<b>LOT</b> Batch
 Store from - to	 Expiration Date
<b>IVD</b> In-Vitro Diagnostic	<b>CE</b> EU CE symbol
 Manufacturer according to 98/79/EG	 Consult instructions for use

**REF**

<b>01.212-02</b>	Anti-Jk <sup>a</sup> for card method monoclonal, human IgM	clone: MS-15	2 ml
<b>01.212-03</b>	Anti-Jk <sup>b</sup> for card method monoclonal, human IgM	clone: MS-15	3 ml
<b>01.212-05</b>	Anti-Jk <sup>a</sup> for card method monoclonal, human IgM	clone: MS-15	5 ml
<b>01.213-02</b>	Anti-Jk <sup>b</sup> for card method monoclonal, human IgM	clone: MS-8	2 ml
<b>01.213-03</b>	Anti-Jk <sup>a</sup> for card method monoclonal, human IgM	clone: MS-8	3 ml
<b>01.213-05</b>	Anti-Jk <sup>b</sup> for card method monoclonal, human IgM	clone: MS-8	5 ml

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**REF**

<b>01.212-02</b>	Anti-Jk <sup>a</sup> for card method monoclonal, human IgM	clone: MS-15	2 ml
<b>01.212-03</b>	Anti-Jk <sup>b</sup> for card method monoclonal, human IgM	clone: MS-15	3 ml
<b>01.212-05</b>	Anti-Jk <sup>2</sup> for card method monoclonal, human IgM	clone: MS-15	5 ml
<b>01.213-02</b>	Anti-Jk <sup>a</sup> for card method monoclonal, human IgM	clone: MS-8	2 ml
<b>01.213-03</b>	Anti-Jk <sup>b</sup> for card method monoclonal, human IgM	clone: MS-8	3 ml
<b>01.213-05</b>	Anti-Jk <sup>2</sup> for card method monoclonal, human IgM	clone: MS-8	5 ml

**CE 0483**

 ANTITOXIN GmbH Industriestrasse 88 69245 Bammental Allemagne  
730-15-5307 Version 007 / 20.05.2021

**Anti-Jk<sup>a</sup>** for card method monoclonal, human IgM

Italiano

**Anti-Jk<sup>b</sup>** for card method monoclonal, human IgM

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**REF**

<b>01.212-02</b>	Anti-Jk <sup>a</sup> for card method monoclonal, human IgM	clone: MS-15	2 ml
<b>01.212-03</b>	Anti-Jk <sup>a</sup> for card method monoclonal, human IgM	clone: MS-15	3 ml
<b>01.212-05</b>	Anti-Jk <sup>a</sup> for card method monoclonal, human IgM	clone: MS-15	5 ml
<b>01.213-02</b>	Anti-Jk <sup>b</sup> for card method monoclonal, human IgM	clone: MS-8	2 ml
<b>01.213-03</b>	Anti-Jk <sup>b</sup> for card method monoclonal, human IgM	clone: MS-8	3 ml
<b>01.213-05</b>	Anti-Jk <sup>b</sup> for card method monoclonal, human IgM	clone: MS-8	5 ml

**CE 0483**

 ANTITOXIN GmbH Industriestrasse 88 69245 Bammental Germania  
730-15-5307 Versione 007 / 20.05.2021

# Anti-Jk<sup>a</sup> for card method monoclonal, human IgM

# Anti-Jk<sup>b</sup> for card method monoclonal, human IgM

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**REF**

- 01.212-02 Anti-Jk<sup>a</sup> for card method monoclonal, human IgM clone: MS-15 2 ml  
 01.212-03 Anti-Jk<sup>a</sup> for card method monoclonal, human IgM clone: MS-15 3 ml  
 01.212-05 Anti-Jk<sup>a</sup> for card method monoclonal, human IgM clone: MS-15 5 ml  
 01.213-02 Anti-Jk<sup>b</sup> for card method monoclonal, human IgM clone: MS-8 2 ml  
 01.213-03 Anti-Jk<sup>b</sup> for card method monoclonal, human IgM clone: MS-8 3 ml  
 01.213-05 Anti-Jk<sup>b</sup> for card method monoclonal, human IgM clone: MS-8 5 ml

CE 0483

 ANTITOXIN GmbH Industriestrasse 88 69245 Bammental Alemania  
 730-15-5307 Versión 007 / 20.05.2021

# Anti-Jk<sup>a</sup> for card method monoclonal, human IgM Anti-Jk<sup>b</sup> for card method monoclonal, human IgM

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ANTITOXIN GmbH Industriestrasse 88 69245 Bammental Almanya