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Anti-E (RH3), monoclonal IgM (human)

REF 690291A

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REF 690291A

INTENDED USE

Monoclonal agglutinating Anti-E reagent is produced from cell culture supernatants of hetero-hybridoma cell lines. The cells secrete an antibody of IgM-type that reacts specifically with the corresponding antigen. The antibody is human protein. The reagent is used for In-Vitro-Diagnostic, to determine whether human red blood cells possess or lack the corresponding blood group antigen E.
The reagent is intended to be used by qualified and technical personnel only.

PRINCIPLE OF PROCEDURE

The procedures used with this reagent is based on the principle of agglutination. Normal human erythrocytes, possessing the corresponding antigen, will agglutinate in the presence of the specific antibody directed towards the antigen.

REAGENT

The listed reagent contains antibodies of the following cell clone:

Anti-E monoclonal, human IgM clones: MS-258, 906

This reagent contains <0.1% (w/v) sodium azide as preservative. Additionally the reagent is comprised of active antibody, sodium chloride, macromolecules and bovine albumin, which has been tested and certified by the US Veterinary service inspectors.

CAUTION: Please handle biological reagent with proper care: This test serum was prepared from supernatants of cell cultures. As biological product it should be looked upon as potentially infectious because of never complete exclusion of danger through excitants of disease. The test serum contains sodium azide that may be toxic and may react with lead or copper to form highly explosive salts. On disposal, flush with large quantities of water. Found from the above reasons test serum should be handled with proper care.

STORAGE REQUIREMENTS

Store at +2 to +8 °C (unopened / opened), or at room temperature while in use. Do not use reagent beyond its labelled expiration date!

REMARKS

1. It is recommended that each lot of reagent be tested with appropriate positive and negative controls.
2. Inappropriate storage impairs efficacy of the reagent.
3. Weak turbidity of the reagent does not affect its reactivity. Bacteria and chemical contamination of the product should be avoided. If a visible change is detected, the reagent should no longer be used, this sign may indicate a microbiological contamination.
4. Strength of positive reactions also depends on age of used blood.
5. Centrifugation outside the specified speed range may lead to false results.
6. The test methods identified below are for manual testing only. When using automated or semi-automated instruments, follow the procedures that are contained in the operator's manual provided by the device manufacturer. Laboratories must follow approved validation procedures
7. For usage of this reagent all effective national laws, directives and guidelines must be followed.
[In Germany especially: „Richtlinien zur Gewinnung von Blut und Blutbestandteilen und zur Anwendung von Blutprodukten (Hämotherapie)“.¹ in its actual form.]

SAMPLE PREPARATION

1. Blood sample should be collected by approved medical procedure.
2. Blood sample to be tested should be used as soon as possible after blood collection to reduce the risk of false-positive and false-negative results due to improper storage or contamination of the sample.
If a delay in testing occurs, sample should be stored at +2 to +8°C.
Blood drawn into EDTA should be tested within 7 days and sample treated with sodium citrate within 14 days after collection.
Blood bag / Donor Blood can be tested until the expiry date.

REAGENT PREPARATION

There is no specific preparation of the reagent required. Use reagent directly from the vials.

PROCEDURE

Material required but not provided:

Slide Method

1. Glass slide
2. Pasteur pipette
3. Mixing stick
4. Timer

Tube Centrifugation Method

1. Test tubes, 10 x 75 mm or 12 x 75 mm
2. Pipettes designed to deliver 50 µL/100 µL
3. Timer
4. Centrifuge
5. Isotonic saline (0,85 - 0,9% Natriumchlorid)
6. Disposable pipette tips

Microplate Method

1. Microplate with U-bottom, (optional: pre-treated)
2. Centrifuge, suitable for centrifugation of microplates
3. Microplate-carrier for centrifuge (option)
4. Microplate shaker
5. Microplatetest mirror (option)
6. Pipette for 50 µL
7. Isotonic saline (0.85 – 0.9 % sodium chloride)
8. Disposable pipette tips

Test procedure

Slide Method

1. Use erythrocyte sediment only.
2. Place one drop (approximately 50 µL) of the reagent on a marked glass slide.
3. Using a Pasteur pipette add one drop erythrocyte sediment (approximately 50 µL) to the drop of reagent on the glass slide.
4. Mix the erythrocytes with reagent well with a stick and spread to a circle with a diameter of approximately 2 cm.
5. By slightly rotating the slide, check for agglutination within 1 minute (reaction starts within seconds).
6. Document the result.

Tube Centrifugation Method

1. Prepare a 2% to 5% suspension of red blood cells in isotonic saline (red blood cells may be washed 1-3 times with isotonic saline).
2. At first put 100 µL of the reagent in a marked tube subsequently add 100 µL of appropriate cell suspension in the tube.
Alternative one drop = approximately 50 µL cell suspension can be added to one drop = approximately 50 µL test serum.
3. Mix Erythrocytes-/Reagent mixture well by slightly shaking.
4. Incubate tube at room temperature for 1-15 min.
5. Centrifugation of tube for 1 minute at 2.000 rpm (approximately 800-1.000 x g).
6. Gently shake the red cells completely from the bottom of the tube and check macroscopically for agglutination within 3 minutes.
7. Document the result.

Microplate Method

1. Prepare 2-5% suspension of red blood cells in isotonic saline (red blood cells may be washed 1-3 times with isotonic saline).
2. Add one drop (50 µL) of reagent into each well.
3. Add one drop (50 µL) of cell suspension into each well.
4. Mix microplate on a microplate shaker.

Note: Recommended : (1) To mix: 10 to 30 seconds at medium strength. (2) To resuspend: 10 to 30 seconds at medium strength respectively with duration and strength, allowing a full resuspension of sediment without destruction of positive reactions.

5. Centrifugation of microplate. Recommended duration of centrifugation: 30 seconds at approximately 400 x g, or as long as necessary to produce best positive reactions with antigen positive cells, with used centrifuge and microplate, and still have an easy resuspension of antigen negative cells.

Note: The used centrifugal force should be the minimum force to produce a sediment and a clear supernatant. Because of multitude offered centrifuges a general speed or duration for best results can not be recommended. Each laboratory has to test and validate its equipment for speed and duration to obtain best results.

6. Loosen the erythrocyte knobs on the shaker and position microplate for reading.
7. Read microplate and document the result.
8. Tests with negative or doubtful results have to be incubated for 5 to 10 minutes at room temperature.
9. Repeat steps 5 to 7 after incubation at room temperature.

Note: All tests should be resuspended and evaluated directly after centrifugation.

INTERPRETATION OF RESULTS

" Slightly rotating / shaking " at Slide Method, Tube Centrifugation Method and Microplate Method:

Positive result (+): visible agglutination of erythrocytes indicates the presence of the corresponding antigen.

Negative result (-): no visible agglutination of erythrocytes indicates the absence of the corresponding antigen.

LIMITATIONS OF THE PROCEDURE

1. The "Procedure" and "Interpretation of Results" sections must be followed closely to assure the accuracy of the test results.
2. No valid conclusion concerning the test result can be reached, if controls with uncertain or false results occur.
3. Enzyme treated erythrocytes or addition of bovine albumin and/or other solutions containing protein may cause unspecific reactions.
4. Hemolyzed, turbid, contaminated or clotted samples should not be used.
5. With the slide method, unspecific reactions might appear due to drying of the reaction-formation or if the slide is heated.
6. Due to variability of antigen expression, reactivity of this reagent against certain phenotypes may produce a weaker reaction compared to control cells.
7. No one specific antiserum or technique can be guaranteed to detect all rare, weak or variant antigens².
8. Some red blood cells may have developed Rh variants that may not be recognized with monoclonal antibodies of a clone.
Therefore the tests should / have to be carried out with at least two different clones / cloning mixtures as required in some countries.
Deviating test results require further testing.
9. Do not use blood of newborn on microplates.

PERFORMANCE

In compliance with Common Technical Specifications (CTS Commission Decision of 03. February 2009) a performance evaluation was conducted. Different samples (donor, patient, panel blood) was used and compared with other products.

Product	Positive	False negative	Sensitivity	Negative	False positive	Specificity
Anti-E monoclonal, human IgM clones: MS-258, 906	327	0	100%	908	0	100%

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SYMBOL - LEGENDE

 Store from - to	 Product Code	 Expiration Date	 Manufacturer as to 98/79/EU
 Lot	 Clone(s)	 In vitro diagnostic medical device	 EU CE-symbol
 Unique Device Identification	 Consult Instrucion for use		

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Antitoxin GmbH, Industriestraße 88, 69245 Bammental, Njemačka