Bromelin solution

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
FOR IN VITRO DIAGNOSTIC USE ONLY

INTENDED USE

Bromelin is a preserved stable enzyme reagent whose active ingredient is produced from pineapple. After appropriate processing of the enzyme protein, the ready to use Bromelin solution is formed.

The reagent is intended to be used by qualified and technical personnel only.

PRINCIPLE OF PROCEDURE

Bromelin solution is used in blood group serology to determination very weak Isoagglutinates, for the compatibility test as well as to enhance the response of immune antibodies. The reagent is therefore used in the routine search for antibodies in sera of patient and probands.

The effect is due to their ability as a proteolytic enzyme to change the surface of erythrocytes in such a way that a better accessibility of some antigens for the corresponding immunoglobulins is achieved.

However, a number of blood group antigens are destroyed by bromelin, such as MNSs, Duffy, Xg, Pr and T.

REAGENTS

The listed reagent is offered in the following form:

Bromelin solution ready for use (1-drop-technique)

The reagent contains 0.01% (w/v) neomycinsulfat as preservative.

Beside of the part of active enzym component, the reagent contains sodium chloride, macromolecules and buffer solutions.

WARNING

The reagent is prepared from the vegetable extract of pineapple.

As biological product it should be looked upon as potentially infectious because of never complete exclusion of danger through excitants of disease. The reagent contains neomycinsulfat that may be toxic.

On disposal, flush with large quantities of water.

For the reasons mentioned above the reagent should be handled with proper care.

STORAGE REQUIREMENT

Store opened and unopened products at +2 to +8°C.

May be at room temperature while in use.

In principle, store and use the reagent to declared expiry date only.

REMARKS

1. With each testing positive and negative controls should be performed.
2. Inappropriate storage impairs efficacy of the reagent
3. Weak turbidity of the reagent does not affect its reactivity.
Bacteria and chemical contamination of the reagent should be avoided.
If a visible change is detected, the reagent should no longer be used, this sign may indicate a microbiological contamination.
4. Strength of positive reactions also depends on age of used blood.
5. Centrifugation highly different from appointed relative centrifugal force may lead to false results.
6. The test method identified below is for manual testing only. When using automated or semi-automated instruments, follow the procedures that are contained in the operator's manual provided by the device manufacturer. Laboratories must follow approved validation procedures.
7. For usage of this reagent all effective national laws, directives and guidelines have to be observed in its current valid version, in Germany especially the „Richtlinien zur Gewinnung von Blut und Blutbestandteilen und zur Anwendung von Blutprodukten (Hämotherapie)“.

SAMPLE PREPARATION

1. Blood sample should be collected by approved medical procedure.
2. Blood samples to be tested should be used as soon as possible after blood collection to reduce the risk of false-positive and false-negative results due to improper storage or contamination of the reagent.
If a delay in testing occurs, samples should be stored at +2 to +8°C.
Blood drawn into EDTA should be tested within 7 days and samples treated with sodium citrate within 14 days after collection.
Blood bag / Donor Blood can be tested until the expiry date.

REAGENT PREPARATION

There is no preparation of the reagent required.

Take and use the reagent directly from the vials.

PROCEDURE

Not included but necessary materials:

- Tube Centrifugation Method:
1. test tubes (10 x 75 mm or 12 x 75 mm)
 2. microliter pipette
 3. centrifuge
 4. timer
 5. isotonic saline (0,85 - 0,9% sodium chloride)

Test procedure

Tube Centrifugation Method

1. a) Prepare 2% to 5% suspensions of red blood cells in isotonic saline (red blood cells may be washed 1-3 times with isotonic saline) or
b) use testerythrocytes of a commercially acquired cell panel.
2. In a labeled test tube
a) give 100µl of an equivalent test serum, or
b) give 100µl serum/plasma of a patient/proband to be tested
3. add to the test tube 100µl
a) of the appropriate cell suspension or
b) testerythrocytes of a coomerccially acquired cell panel.
4. Then add one drop of Bromelin solution
5. Mix Erythrocytes / Reagentmixture well by slightly shaking.
6. Incubate tube at room temperature for 5-15 min.
7. Centrifugation of the tube for 1 minute at 1.000 rpm (approximately 180-270 x g).
8. Gently shake the red cells completely from the bottom of the tube and check macroscopically for agglutination within 3 minutes.
9. Document the result.

INTERPRETATION OF RESULTS

"Slightly shaking" at Tube Centrifugation Method:

- Positive results (+): Visible agglutination of erythrocytes is a positive result and indicates
- a) the presence of the corresponding antigen
 - b) the presence of one or several antibodies which can be identified with the help of the antigram of the cell panel.
- Negative results (-): No visible agglutination of erythrocytes is a negative result and indicates
- a) the absence of the corresponding antigen
 - b) an irregular antibody does not exist.

LIMITATIONS OF THE PROCEDURE

1. Inaccuracy at compliance with instructions written under section "Procedures" and "Interpretation of results" may lead to incorrect results.
2. No valid conclusion concerning the test result can be reached, if controls with uncertain or false results occur.
3. Hemolyzed, turbid, contaminated or clotted blood samples must not be used in this test.
4. No one specific antiserum or technique can be guaranteed to detect all rare, weak or variant of antigens.²
5. Red blood cells coated with alloantibodies or autoantibodies (i.e., cells that are positive in the direct antiglobulin-test (DAT)) are unsuitable for testing.
6. Some antigens as e.g. MNSs, Duffy, Xg, Pr or T are destroyed by Bromelin.
7. The concentration, activity and volume used of the enzyme are very important. Due to incorrect or too long storage of the enzyme solution (inactivation of the enzyme), but also the use of too much enzyme solution (destruction of antigens), can lead to false negative results. In case of a negative test result, it is recommended to check the enzyme activity with "Bromelin control solution". Bromelin control solution only reacts positively if bromelin has been added to the test approach with sufficient activity.

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REF	Product Code	LOT	Lot
	Store from - to		Expiration Date
IVD	In-Vitro Diagnostic	EU CE symbol	
	Manufacturer as to 98/79/EG		Observe instruction for use
UDI	Unique Device Identification		Distributor

REF

BG1613 Bromelin solution 10 ml

730-13-1906 Version 006 / 01.04.2022

 ANTITOXIN GmbH Industriestraße 88 69245 Bammmental Germany
 +49 (0) 6223/ 8661-0  +49 (0) 6223/ 8661-13  gara@antitoxin-gmbh.de

sifin

sifin diagnostics gmbh
Berliner Allee 317-321
13088 Berlin / Germany
 +49 30 700 144-0