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FOR IN VITRO DIAGNOSTIC USE ONLY

INTRODUCTION

Monoclonal agglutinating Anti-Le^a and -Le^b reagents are produced from cell culture supernatants of hybridoma-cell lines, which secreting antibodies that reacts specifically with the corresponding antigen. The reagents are used to in-vitro-determine qualitative whether human red blood cells possess or lack the corresponding blood group antigens Le^a and Le^b. The reagents are intended to be used by qualified and technical personnel only.

PRINCIPLE OF PROCEDURE

The test method used with these reagents is based on the principle of column agglutination. Normal human erythrocytes, possessing one of these antigens, will agglutinate in the presence of the specific antibody directed toward the antigen.

REAGENTS

The listed reagents contain antibodies of the following clones:

Anti-Le^a for card method monoclonal, mouse clone: GA2

Anti-Le^b for card method monoclonal, mouse clone: LEB1

Both reagents contain <0.1% (w/v) sodium azide as preservative. Beside the parts active antibody the reagents contain sodium chloride, macromolecules and bovine albumin, which has been tested and certified by the US Veterinary service inspectors.

WARNING

These reagents were prepared from supernatants of cell cultures. Nevertheless as biological products it should be looked upon as potentially infectious because of never complete exclusion of danger through excitants of disease. The reagents contain sodium azide that may be toxic and may react with lead or copper to form highly explosive salts.

On disposal, flush with large quantities of water.

For the reasons mentioned above the reagents should be handled with proper care.

STORAGE REQUIREMENT

Store opened and unopened products at 2 to 8 °C. May be at room temperature while in use. In principle store and use the reagents only to declared expiry date.

REMARKS

- With each testing positive and negative controls should be performed.
- Inappropriate storage impairs efficacy of the reagents.
- Weak turbidity of the reagents does not affect its reactivity. Bacterial and chemical contamination of the reagents should be avoided. If a visible change is detected, the reagents should no longer be used this sign may indicate a microbiological contamination.
- Strength of positive reactions also depends on age of used blood.
- Use the Grifols card centrifuge. The use of another card-specific centrifuge (each card centrifuge has its specified unchangeable g-force) may lead to false results due to the changed g-force.
- The test method identified below is for manual testing only. When using automated or semi-automated instruments, follow the procedures that are contained in the operator's manual provided by the device manufacturer. Laboratories must follow approved validation procedures.
- For usage of these reagents all effective national laws, directives and guidelines have to be observed, in Germany especially the „Richtlinien zur Gewinnung von Blut und Blutbestandteilen und zur Anwendung von Blutprodukten (Hämotherapie)“¹ in its actual form.
- The information on the use of the Grifols test cards in the relevant insert must be observed.

SAMPLE PREPARATION

- Blood sample should be collected by approved medical procedure.
- Blood samples to be tested should be used as soon as possible after blood collection to reduce the risk of false-positive and false-negative results due to improper storage or contamination of the sample. If a delay in testing occurs, samples should be stored at 2 to 8°C. Blood drawn into EDTA should be tested within 7 days and samples treated with sodium citrate within 14 days after collection. Blood bag / Donor Blood can be tested until the expiry date.

REAGENT PREPARATION

There is no preparation of the reagents required. Use reagents directly from the vials.

PROCEDURE

Not provided material additionally needed.

Card Method:

- Grifols "DG Gel Neutral"
- Microliter pipette
- Tubes
- Grifols Card centrifuge
- DG Gel Sol "card specific diluent"

Test procedure

Card Method

- Prepare 0,8 % suspensions of red blood cells in DG Gel Sol (card specific diluent).
- Add 50 µL of appropriate cell suspension to each labeled micro tube.
- Add 25 µL of appropriate reagent to each micro tube.
- Only for the reagent **Anti-Le^a**: Incubate the card at room temperature for 0-15 minutes.
- Centrifugation of card in appropriate centrifuge with the, for the respective centrifuge, unchangeable g-force.
- Check macroscopically for agglutination within 30 minutes.
- Document the result.

INTERPRETATION OF RESULTS

Positive results (+): visible agglutination of erythrocytes is a positive result and indicates the presence of the corresponding antigen.

Negative results (-): No visible agglutination of erythrocytes is a negative result and indicates the absence of the corresponding antigen.

Read and interpret the results of the card method according to the Grifols card instruction for use.

LIMITATIONS OF THE PROCEDURE

- Inaccuracy at compliance with instructions written under section "Procedures" and "Interpretation of results" may produce incorrect results.
- No valid conclusion concerning the test result can be reached, if controls with uncertain or false results occur.
- Enzyme treated erythrocytes or addition of bovine and/or other solutions containing protein may cause unspecific reactions.
- Hemolyzed, turbid, contaminated or clotted blood samples must not be used in this test.
- Due to variability of antigen expression on human red cells, reactivity of the reagents, mentioned above, against certain phenotypes may give weaker reactivity compared to control cells.
- No one specific antiserum or technique can be guaranteed to detect all rare, weak or variant antigens.²
- Red blood cells coated with alloantibodies or autoantibodies of the same or similar specificity as the reagent (i.e., cells that are positive in the direct antiglobulin-test (DAT)) are not suitable for this test procedure.
- Red blood cells with a positive direct Coombs test may cause false-positive reactions in card method.
- Pay attention to all statements to limitations in inserts of Grifols card.
- In some cases, Le^b positive cells of blood group A, A₁B and A₂B may experience attenuated or absent responses with the reagent Anti-Le^b.
- Incubation at room temperature could increase reaction of Anti-Le^b.
- Positive weak reactions need to be confirmed with extensive testing. Anti-Le^a could weakly react with patients blood group A₂ or O with the antigen Le(a-b+).
- For donors, a small pellet of non-agglutinated red blood cells in the bottom of the microtube may be observed in positive reactions (image of double population, DP). In patients, DPs could also be caused by prior transfusion. Additional information on patient history could be necessary for resolution.


LITERATURE

- Richtlinien zur Gewinnung von Blut und Blutbestandteilen und zur Anwendung von Blutprodukten (Hämotherapie)
- CLSI, I/LA33-A Validation of Automated System for Immunhematological Testing Before Implementation: Approved Guideline CLINICAL AND LABORATORY STANDARDS INSTITUTE December 2009
- Peter D. Issitt, David J. Anstee Applied Blood Group Serology, fourth edition, Montgomery Scientific Publications 1998
- Christian Mueller-Eckhardt, Volker Kiefel Transfusionsmedizin, 3. Auflage, Springer-Verlag 2004.

PRESENTATION

213217 Anti-Le^a for card method monoclonal, mouse clone: GA2 5 ml
213219 Anti-Le^b for card method monoclonal, mouse clone: LEB1 5 ml

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Antitoxin GmbH - Industriestraße 88 - 69245 Bammental