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730-13-2706 Version 006 / 15. Januar 2021

Antitoxin GmbH - Industriestraße 88 - 69245 Bammental

For Slide and Tube Test
FOR IN VITRO DIAGNOSTIC USE ONLY

INTENDED USE

Monoclonal agglutinating Anti-P₁ reagent is produced from cell culture supernatants of hybridoma-cell line. The cells are secreting an antibody of IgM-type that reacts specific with the corresponding antigen. The reagent is used to in-vitro determine qualitative whether red blood cells possess or lack the corresponding blood group antigen. The reagent is intended to be used by qualified and technical personnel only.

PRINCIPLE OF PROCEDURE

The test methods used with this reagent are based on the principle of agglutination. Normal human erythrocytes, possessing the antigen, will agglutinate in the presence of the corresponding antibody.

REAGENTS

The listed reagent contains antibodies of the following cell clone:

Anti-P₁ monoclonal, mouse clone: 650

The reagent contains <0.1% (w/v) sodium azide as preservative. Beside the part active antibody, the reagent contains sodium chloride, macromolecules and bovine albumin, which has been tested and certified by the US Veterinary service inspectors.

WARNING

This reagent is prepared from supernatants of cell cultures. As biological product it should be looked upon as potentially infectious because of never complete exclusion of danger through excitants of disease. The reagent contains sodium azide, that may be toxic and may react with lead or copper to form highly explosive salts. On disposal, flush with large quantities of water.

For the reasons mentioned above the reagent should be handled with proper care.

STORAGE REQUIREMENT

Store opened and unopened products at 2 to 8°C.

May be at room temperature while in use.

In principle, store and use the reagent to declared expiry date only.

REMARKS

1. With each testing positive and negative controls should be performed.
2. Inappropriate storage impairs efficacy of the reagent.
3. Weak turbidity of the reagent does not affect its reactivity. Bacteria and chemical contamination of the product should be avoided. If a visible change is detected, the reagent should no longer be used, this sign may indicate a microbiological contamination.
4. Strength of positive reactions also depends on age of used blood.
5. Centrifugation highly different from appointed relative centrifugal force may lead to false results.
6. The test method identified below is for manual testing only. When using automated or semi-automated instruments, follow the procedures that are contained in the operator's manual provided by the device manufacturer. Laboratories must follow approved validation procedures.
7. For usage of this reagent all effective national laws, directives and guidelines have to be observed, in Germany especially the „Richtlinien zur Gewinnung von Blut und Blutbestandteilen und zur Anwendung von Blutprodukten (Hämotherapie)“²¹ in its actual form.

SAMPLE PREPARATION

1. Blood sample should be collected by approved medical procedure.
2. Blood samples to be tested should be used as soon as possible after blood collection to reduce the risk of false-positive and false-negative results due to improper storage or contamination of the sample. If a delay in testing occurs, samples should be stored at 2 to 8°C. Blood drawn into EDTA should be tested within 7 days and samples treated with sodium citrate within 14 days after collection. Blood bag / Donor Blood can be tested until the expiry date.

REAGENT PREPARATION

There is no preparation of the reagent required. Take and use the reagent directly from the vials.

PROCEDURE

Not included but necessary materials:

Slide Method

1. Glass slide
2. Pasteur pipette
3. Mixing stick
4. Timer

Tube Centrifugation Method

1. Test tubes, 10 x 75 mm or 12 x 75 mm
2. Microliter pipette
3. Timer
4. Centrifuge
5. Isotonic saline (0,85 - 0,9% sodium chloride)

Test procedure

Slide Method

1. Use erythrocyte sediment only.
2. Place one drop (approximately 50 µl) of appropriate reagent on a glass slide.
3. Using a Pasteur pipette add one drop erythrocyte sediment (approximately 50 µl) to the drop of reagent on the glass slide.
4. Mix the erythrocytes with reagent well with a stick and spread to a circle with a diameter of approximately 2 cm.
5. By slightly rotating the slide, check for agglutination within 1 minute (reaction starts within seconds).

Tube Centrifugation Method

1. Prepare 2-5 % suspensions of red blood cells in Isotonic saline (red blood cells may be washed one time or up to three times with isotonic saline)
2. At first put 100 µl of appropriate reagent in each marked tube subsequently add 100 µl of appropriate cell suspension in each tube. Alternative one drop = 50 µl cell suspension can be added to one drop = 50 µl reagent.
3. Mix Erythrocytes-/Reagent-mixture well by slightly shaking.
4. Incubate tube at room temperature for 15 min.
5. Centrifugation of the tubes for 1 minute at 2.000 rpm (approximately 800-1.000 x g).
6. Gently shake the red cells from the bottom of the tube and check macroscopically for agglutination within 3 minutes
7. Document the results.

INTERPRETATION OF RESULTS

"Slightly rotating/shaking" at Slide Method and at Tube Centrifugation Method

Positive results (+): visible agglutination of erythrocytes is a positive result and indicates the presence of the corresponding antigen.

Negative results (-): No visible agglutination of erythrocytes is a negative result and indicates the absence of the corresponding antigen.

LIMITATIONS OF THE PROCEDURE

1. Inaccuracy at compliance with instructions written under section "Procedures" and "Interpretation of results" may lead to incorrect results.
2. No valid conclusion concerning the test result can be reached, if controls with uncertain or false results occur.
3. Enzyme treated erythrocytes or addition of bovine albumin and/or other solutions containing protein may cause unspecific reactions.
4. Hemolyzed, turbid, contaminated or clotted blood samples must not be used in this test.
5. Due to variability of antigen expression on human red blood cells, reactivity of the reagent, mentioned above, against certain phenotypes may give weaker reactivity compared to control cells.
6. With the slide method unspecific reactions might appear due to drying of the reaction-formation or if the slide is heated.
7. No one specific antiserum or technique can be guaranteed to detect all rare, weak or variant antigens.²

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213071 Anti-P₁ monoclonal, mouse clone: 650 5 ml

730-13-2706 Version 006 / 15. January 2021

	In vitro diagnostics device
	Batch Code
	Use by
	Temperature limitation
	Consult instructions for use
	Catalogue Number
	Manufacturer

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Antitoxin GmbH - Industriestrasse 88 - 69245 Bammental