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CE 0483

Antitoxin GmbH - Industriestraße 88 - 69245 Bammental

For Slide, Tube, Microplate and Card Method
FOR IN VITRO DIAGNOSTIC USE ONLY

INTENDED USE

Monoclonal agglutinating Anti-K-reagents are produced from cell culture supernatants of hetero-hybridoma-cell lines. The cells are secreting antibodies of IgM-type that reacts specific with the corresponding antigen. The antibody is human protein. The reagents are used to in-vitro-determine qualitative whether human red blood cells possess or lack the corresponding blood group antigen K. The reagents are intended to be used by qualified and trained technical personnel only.

PRINCIPLE OF PROCEDURE

The procedures used with these reagents are based on the principle of agglutination. Normal human erythrocytes, possessing the corresponding antigen, will agglutinate in the presence of the specific antibody directed toward the antigen.

REAGENTS

The listed reagents contain antibodies of the following cell clones:

Anti-K (KEL1) monoclonal, human IgM clone: MS-56

Anti-K (KEL1) monoclonal, human IgM clone: AEK-4

The reagents contain <0.1% (w/v) sodium azide as preservative. Additionally the reagents are prepared of active antibody, sodium chloride, macromolecules and bovine albumin, which has been tested and certified by the US Veterinary service inspectors.

WARNING

These reagents are prepared from supernatants of cell cultures. As biological products they should be looked upon as potentially infectious because of never complete exclusion of danger through excitants of disease. The reagents contain sodium azide that may be toxic and may react with lead or copper to form highly explosive salts.

On disposal, flush with large quantities of water.

For the reasons mentioned above reagents should be handled with proper care.

STORAGE REQUIREMENT

Store opened and unopened products at 2 to 8°C. May be at room temperature while in use. In principle, store and use the reagents to declared expiry date only.

REMARKS

- With each testing positive and negative controls should be performed.
- Inappropriate storage impairs efficacy of the reagents.
- Weak turbidity of the reagents does not affect its reactivity. Bacteria and chemical contamination of the product should be avoided. If a visible change is detected, it may indicate a microbiological contamination, the reagent should no longer be used.
- Strength of positive reactions also depends on age of used blood.
- Centrifugation outside the specified speed range may lead to false results. In the card method use the Grifols card centrifuge. The use of another card-specific centrifuge (each card centrifuge has its specified unchangeable g-force) may lead to false results due to the changed g-force.
- The test method identified below is for manual testing only. When using automated or semi-automated instruments, follow the procedures that are contained in the operator's manual provided by the device manufacturer. Laboratories must follow approved validation procedures.
- For usage of these reagents all effective national laws, directives and guidelines have to be observed, in Germany especially the „Richtlinien zur Gewinnung von Blut und Blutbestandteilen und zur Anwendung von Blut-produkten Hämotherapie“¹ in its actual form.
- The information on the use of the Grifols test cards in the relevant insert must be observed.

SAMPLE PREPARATION

- Blood sample should be collected by approved medical procedure.
- Blood samples to be tested should be used as soon as possible after blood collection to reduce the risk of false-positive and false-negative results due to improper storage or contamination of the sample.

If a delay in testing occurs, samples should be stored at 2 to 8°C.

Blood drawn into EDTA should be tested within 7 days and samples treated with sodium citrate within 14 days after collection.

Blood bag / Donor Blood can be tested until the expiry date.

REAGENT PREPARATION

There is no preparation of the reagents required. Use reagents directly from the vials.

PROCEDURE

Not provided material, additionally needed:

at Slide Method: Glass slide; Pasteur pipette; Mixing stick; timer

at Tube Centrifugation Method: Test tubes, 10 x 75 mm or 12 x 75 mm; microliter pipette; centrifuge; timer; isotonic saline (0,85 - 0,9% sodium chloride)

at Card Method: Cards: Grifols "DG Gel Neutral"; Test tubes, 10 x 75 mm or 12 x 75 mm; microliter pipette; Grifols card centrifuge; card specific diluent (Grifols "DG Gel Sol");

at Microplate Method: Microplate with 96 U-wells; Test tubes, 10 x 75 mm or 12 x 75 mm; microliter pipette; timer; microplate-centrifuge; microplate-shaker; Isotonic saline (0,85 - 0,9% sodium chloride)

Test procedure

Slide Method (use only with clone: MS-56)

- Use erythrocyte sediment only.
- Place one drop (approximately 50 µL) of appropriate reagent on a glass slide.
- Using a Pasteur pipette add one drop erythrocyte sediment (approximately 50 µL) to the glass slide.
- Mix the erythrocytes with reagent well with a stick and spread to a circle with a diameter of approximately 2 cm.
- By slightly rotating the slide, check for agglutination within 1 minute (reaction starts within seconds). Unspecific reactions might appear due to drying of the reaction-formation or if the slide is heated.
- Document the result.

Tube Centrifugation Method

- Prepare 2% to 5% suspensions of red blood cells in isotonic saline only (cells washed one time or up to three times with isotonic saline).
- At first put 100 µL of appropriate reagent in a marked tube subsequently add 100 µL of appropriate cell suspension in the tube. Alternative one drop = approximately 50 µL cell suspension can be added to one drop = approximately 50 µL serum
- Mix well by slightly shaking.
- Incubate tube at room temperature for 15 min.
- Centrifugation of tube for 1 minute at 2.000 rpm (approximately 800-1.000 x g).
- Gently resuspend the red cells and check macroscopically for agglutination within 3 minutes.
- Document the result.

Card Method

- Prepare 0,8 % suspensions of red blood cells in DG Gel Sol (card specific diluent).
- Add 50 µL of appropriate cell suspension to each marked micro tube.
- Add 25 µL of the appropriate reagent to each micro tube.
- Centrifuge the card in a Grifols card centrifuge with the for this centrifuge unchangeable g-force within 30 minutes.
- Check macroscopically for agglutination within 30 minutes.
- Document the result.

Microplate Method

- Prepare 2% to 5% suspensions of red blood cells in isotonic saline (cells may be washed one time or up to three times with isotonic saline).
- Add 50 µL of appropriate reagent to a marked well.
- Add 50 µL of appropriate cell suspension to the well.
- To mix both shake microplate for 30 seconds on microplate-shaker with medium speed.
- Centrifuge microplate in appropriate microplate-centrifuge for 30 seconds at 400 x g.
- Shake microplate for 30 seconds on a microplate shaker with medium speed.
- Evaluated macroscopically for agglutination directly after centrifugation.
- Document the result.
- Tests with negative or doubtful results have to be incubated for 5 to 10 Minutes at room temperature.
- Repeat steps 5 to 8 after incubation at room temperature

INTERPRETATION OF RESULTS

"Slightly rotating/shaking" at.

Slide Method, Tube Centrifugation Method; Microplate Method:

Positive results (+): visible agglutination of erythrocytes is a positive result and indicates the presence of the corresponding antigen.

Negative results (-): No visible agglutination of erythrocytes is a negative result and indicates the absence of the corresponding antigen.

Read and interpret the results of the card method according to the Grifols card instruction for use.

LIMITATIONS OF THE PROCEDURE

- Inaccuracy at compliance with instructions written under section "Procedures" and "Interpretation of results" may lead to incorrect results.
- No valid conclusion concerning the test result can be reached, if controls with uncertain or false results occur.
- Enzyme treated erythrocytes or addition of bovine albumin and/or other solutions containing protein may cause unspecific reactions.
- Hemolyzed, turbid, contaminated or clotted blood samples must not be used in this test.
- With the slide method, unspecific reactions might appear due to drying of the reaction-formation or if the slide is heated.
- Due to variability of antigen expression on human red blood cells, reactivity of the reagents, mentioned above, against certain phenotypes may give weaker reactivity compared to control cells.
- No one specific antiserum or technique can be guaranteed to detect all rare, weak or variant antigens. ²
- Red blood cells coated with alloantibodies or autoantibodies of the same or similar specificity as that reagent used for the test (i.e. cells that are positive in the direct antiglobulin-test (DAT)) are not suitable for this test procedure.
- Red blood cells with a positive direct Coombs-test may cause false-positive reactions in card method. These cells react positive without testserum too.
- Pay attention to all statements to limitations in inserts of Grifols card.

PERFORMANCE CHARACTERISTICS

In compliance with Common Technical Specifications (CTS Commission Decision of 03. February 2009) a performance evaluation was conducted. Different samples (donor, patient, neonates-, panel blood) was used and compared with other reference methods / products.

Product	Positive	Negative	False Positive	False Negative
Anti-K (clone: MS-56)	614	1302	0	0
Anti-K (clone: AEK-4)	569	546	0	0

The calculated values found for:

Product	Sensitivity	Specificity
Anti-K (clone: MS-56)	100%	100%
Anti-K (clone: AEK-4)	100%	100%

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
Anti-K (KEL1) monoclonal, IgM human clone: MS-56
Anti-K (KEL1) monoclonal, IgM human clone: AEK-4

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Antitoxin GmbH - Industriestrasse 88 - 69245 Bammental