# Control solution

for monoclonal reagents

# Deutsch

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**REF**

**33110** Control solution for monoclonal reagents 10ml  
**33110.X** Control solution for monoclonal reagents 10 x 10ml

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for monoclonal reagents

English

for Slide-, Tube-, Microplate and Card-Method  
FOR IN VITRO DIAGNOSTIC USE ONLY

## INTENDED USE

A doubtful result with a monoclonal reagent may occur by autoimmunoantibodies, protein anomalies, rouleaux formation or a positive DCT. With each testing with monoclonal reagents a control should be done with Control solution to demonstrate that a positive result depends on the specific antibody.

Control solution for monoclonal reagents can be used in slide-, tube-, microplate and card method. The control reagent is intended to be used by qualified and technical personnel only.

## PRINCIPLE OF PROCEDURE

The test method used with this control reagent is based on the principle of agglutination. Control-solution for monoclonal reagents is a control reagent containing all ingredients of monoclonal testsera with the exception of a specific antibody. Is Control solution used instead of a monoclonal reagent the result has to be negative. A positive result with Control solution with erythrocytes tested in recommended techniques demands a retesting of specific reaction of used monoclonal reagent.

## REAGENTS

The listed reagent is available in following form:

Control solution for monoclonal reagent

The reagent contains <0.1% (w/v) sodium azide as preservative. Additional the reagent contains sodium chloride, macromolecules and bovine albumin, which has been tested and certified by the US Veterinary service inspectors.

## WARNING

The reagent is prepared from biological products.

Nevertheless, as biological product it should be looked upon as potentially infectious because of never complete exclusion of danger through excitants of disease. The reagent contains sodium azide that may be toxic and may react with lead or copper to form highly explosive salts. On disposal, flush with large quantities of water.

For the reasons mentioned above, the reagent should be handled with proper care.

## STORAGE REQUIREMENT

Store opened and unopened products at +2 to +8°C.

May be at room temperature while in use.

In principle, store and use the reagent to indicated expiry date only.

## REMARKS

1. Inappropriate storage impairs efficacy of the reagent.
2. Weak turbidity of the reagent does not affect its reactivity. Bacteria and chemical contamination of the product should be avoided. If a visible change is detected, it may indicate a microbiological contamination, the reagent should no longer be used.
3. Centrifugation outside the specified speed range may lead to false results.
4. The test methods identified below are for manual testing only. When using automated or semi-automated instruments, follow the procedures that are contained in the operator's manual provided by the device manufacturer. Laboratories must follow approved validation procedures.
5. For usage of the reagent all effective national laws, directives and guidelines have to be observed in its current valid version, in Germany especially the „Richtlinien zur Gewinnung von Blut und Blutbestandteilen und zur Anwendung von Blutprodukten (Hämotherapie)“<sup>1</sup>.

## SAMPLE PREPARATION

1. Blood sample should be collected by approved medical procedure.
2. Blood sample to be tested should be used as soon as possible after blood collection to reduce the risk of false-positive and false-negative results due to improper storage or contamination of the sample.  
If a delay in testing occurs, sample should be stored at +2 to +8°C.  
Blood drawn into EDTA should be tested within 7 days and sample treated with sodium citrate within 14 days after collection.  
Blood bag / Donor Blood can be tested until the expiry date.

## REAGENT PREPARATION

There is no preparation of the reagent required.

Take and use reagent directly from the vials.

## PROCEDURE

Not provided material, additionally needed

- |                             |                                                                                                                                                                                                                                                                                                                                                                                                                               |
|-----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Slide Method:               | <ol style="list-style-type: none"><li>1. glass slide</li><li>2. pasteur pipette</li><li>3. mixing stick</li><li>4. timer</li></ol>                                                                                                                                                                                                                                                                                            |
| Tube Centrifugation Method: | <ol style="list-style-type: none"><li>1. tubes (10 x 75 mm or 12 x 75 mm)</li><li>2. microliter pipette</li><li>3. timer</li><li>4. centrifuge</li><li>5. isotonic saline (0,85 - 0,9% sodium chloride)</li></ol>                                                                                                                                                                                                             |
| Microplate Method:          | <ol style="list-style-type: none"><li>1. microplate with 96 U-wells</li><li>2. microliter pipette</li><li>3. centrifuge</li><li>4. timer</li><li>5. microplate shaker</li><li>6. microplate centrifuge</li><li>7. isotonic saline (0,85 - 0,9% sodium chloride)</li></ol>                                                                                                                                                     |
| Card Method:                | <ol style="list-style-type: none"><li>1. Cards: - BIO-RAD (<i>DiaMed</i>) „NaCl, Enzyme Test and cold Agglutinins“<br/>- Grifols „DG Gel Neutral“;<br/>- Ortho BioVue® System Reverse</li><li>2. test tubes</li><li>3. microliter pipette</li><li>4. timer</li><li>5. centrifuge</li><li>6. corresponding card centrifuge</li><li>7. card specific diluent</li><li>8. isotonic saline (0,85 - 0,9% sodium chloride)</li></ol> |

## Test procedure

### Slide Method

1. Use erythrocyte sediment only.
2. Place one drop (approximately 50 µL) of the Control solution for monoclonal reagents on a marked glass slide.
3. Using a Pasteur pipette add one drop erythrocyte sediment (approximately 50 µL) to the drop of reagent on the glass slide.
4. Mix the erythrocytes with the control reagent well with a stick and spread to a circle with a diameter of approximately 2 cm.
5. By slightly rotating the slide, check for agglutination within 1 minute (reaction starts within seconds).
6. Document the result.

### Tube Centrifugation Method

1. Prepare 2% to 5% suspensions of red blood cells in isotonic saline (red blood cells may be washed 1-3 times with isotonic saline).
2. At first put 100 µL of the Control solution for monoclonal reagents in a marked tube, subsequently add 100 µL of appropriate cell suspension.  
Alternative one drop = approximately 50 µL cell suspension can be added to one drop = approximately 50 µL Control solution for monoclonal reagents.
3. Mix Erythrocytes-/Control solution mixture well by slightly shaking.
4. Incubate tube at room temperature for 5-15 min.
5. Centrifuge tube for 1 min at 2.000 rpm (approximately 800-1.000 x g).
6. Gently shake the red cells from the bottom of the tube and check macroscopically for agglutination within 3 minutes.
7. Document the result.

### Microplate Method

1. Prepare 2-5% suspension of red blood cells in isotonic saline (red blood cells may be washed 1-3 times with isotonic saline).
2. Add 50 µL of appropriate Control solution for monoclonal reagents to a marked well.
3. Add 50 µL of appropriate cell suspension to the well.
4. To mix both shake microplate for 30 seconds on microplate-shaker with medium speed.
5. Centrifuge microplate in appropriate microplate-centrifuge for 30 seconds at 400 x g.
6. Shake microplate for 30 seconds on a microplate shaker with medium speed.
7. Evaluated macroscopically for agglutination directly after centrifugation.
8. Document the result.
9. Tests with negative or doubtful results have to be incubated for 5 to 10 minutes at room temperature.
10. Repeat steps 5 to 8 after incubation at room temperature.

### Card Method (manual method / valid for the cards:

- BIO-RAD (*DiaMed*) „NaCl, Enzyme Test and cold Agglutinins“
- Grifols „DG Gel Neutral“

1. Prepare 0,8 % suspension of red blood cells in card specific diluent (red blood cells may be washed 1-3 times with isotonic saline).
2. Add 50 µL of appropriate cell suspension to a marked micro tube.
3. Add 25 µL of the Control solution for monoclonal reagents to the micro tube.
4. Centrifuge the card after latest 30 minutes in appropriate card centrifuge with the for this centrifuge unchangeable g-force.
5. Check macroscopically micro tube for agglutination within 30 minutes.
6. Document the result.

### Ortho BioVue® System

1. Prepare 3-5 % suspension of red blood cells in isotonic saline (red blood cells may be washed 1-3 times with isotonic saline).
2. Add 40 µL of the Control solution for monoclonal reagents to a marked reaction chamber.
3. Add 10 µL of appropriate cell suspension to the reaction chamber.
4. Centrifuge the cassette after latest 15 minutes in appropriate card centrifuge with the for this centrifuge unchangeable g-force.
5. The test results should be read macroscopically for agglutination directly after the end of centrifugation.
6. Document the result.

## INTERPRETATION OF RESULTS

"Slightly rotating" at Slide Method

"Slightly shaking" at Tube Centrifugation Method, .

Positive results (+): visible agglutination of erythrocytes is a positive result and indicates here a doubtful result of the used monoclonal reagent.

Negative results (-): No visible agglutination of erythrocytes is a negative result and indicates the reactions of the monoclonal testsera can be used.  
There is no incorrect determination.

Read and interpret the results of the card method according to the instruction for use of the corresponding card.

## LIMITATIONS OF THE PROCEDURE

1. Inaccuracy at compliance with instructions written under section "Procedures" and "Interpretation of results" may lead to incorrect results.
2. No valid conclusion concerning the test result can be reached, if controls with uncertain or false results occur.
3. Enzyme treated erythrocytes or addition of bovine albumin and/or other solutions containing protein may cause unspecific reactions.
4. Hemolyzed, turbid, contaminated or clotted samples should not be used.
5. With the slide method, unspecific reactions might appear due to drying of the reaction-formation or if the slide is heated.
6. Red blood cells coated with alloantibodies or autoantibodies of the same or similar specificity as that reagent used for the test (i.e., cells that are positive in the direct antiglobulin-test (DAT)) are not suitable for this test procedure.

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CE-IMMUNDIAGNOSTIKA GmbH  
Karl-Landsteiner-Str. 6 D-69151 Neckargemuend

|                                                                                   |                              |                                                                                   |                             |
|-----------------------------------------------------------------------------------|------------------------------|-----------------------------------------------------------------------------------|-----------------------------|
| <b>REF</b>                                                                        | Product Code                 | <b>LOT</b>                                                                        | Lot                         |
|  | Store from - to              |  | Expiration Date             |
| <b>IVD</b>                                                                        | In-Vitro Diagnostic          | <b>CE</b>                                                                         | EU CE symbol                |
|  | Manufacturer as to 98/79/EG  |  | Observe instruction for use |
| <b>UDI</b>                                                                        | Unique Device Identification |  | Distributor                 |

|                |                                                    |
|----------------|----------------------------------------------------|
| <b>REF</b>     |                                                    |
| <b>33110</b>   | Control solution for monoclonal reagents 10ml      |
| <b>33110.X</b> | Control solution for monoclonal reagents 10 x 10ml |

730-13-3012 Version 012a \_\_ / 19.11.2025

 Antitoxin GmbH Industriestraße 88 69245 Bammental Germany  
 +49 (0) 6223/ 8661-0  +49 (0) 6223/ 8661-13  [qara@antitoxin-gmbh.de](mailto:qara@antitoxin-gmbh.de)