AHG serum blend

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 +49 (0) 6223/ 8009400  +49 (0) 6223/ 8009499  info@ce.immundiagnostika.com

AHG serum blend

Coombs serum, polyspecific,
polyclonal, rabbit Anti-IgG / monoclonal, mouse Anti-C3 clone: Bric-8

English

FOR IN VITRO DIAGNOSTIC USE ONLY

INTENDED USE

The product itself has no dedicated intended use, but can only be used in combination with other diagnostic products for diagnostic procedures.

USER

The use of this AHG serum blend is for qualified and trained personnel only.

PRINCIPLE OF THE PROCEDURE

This polyspecific AHG serum blend contains antibodies to human immunoglobulins as well as complement factors. For the procedure and method, please refer to the instructions for use of the respective test serum used in combination with this reagent.

COMPOSITION

This AHG serum blend contains polyclonal reagent from immunised rabbit plasma (Anti-IgG) and monoclonal reagent obtained from supernatants of cell cultures of mouse hybridoma cell line Bric-8 (Anti-C3).

The reagent is coloured green.

The reagent contains < 0.1% (w/v) sodium azide as preservative.

Additionally, the reagent contains sodium chloride, macromolecules and bovine albumin (BSA).

The BSA used in this formulation is derived from US-sourced animals from USDA and APHIS-approved facilities. It complies with Regulations (EC) No 1069/2009 and (EU) No 142/2011 for use in vitro diagnostic reagents and has been certified free from both Vesicular Stomatitis Virus and Bluetongue Virus.

WARNING

The AHG serum blend is prepared from animal plasma and supernatants of cell cultures. As a biological product, this reagent should be regarded as potentially infectious, since the risk of pathogen transmission can never be completely excluded. The reagent contains sodium azide, a toxic substance that may react with lead or copper to form highly explosive salts. Upon disposal, flush with large quantities of water.

For the reasons mentioned above, the reagent should be handled with proper care.

STORAGE REQUIREMENT

Store opened and unopened products at +2 to +8 °C. It may be kept at room temperature while in use. Use only until the declared expiry date, which is given in the format year-month-day (YYYY-MM-DD).

REMARKS

Follow the instructions for use for the respective test serum used in combination with this product.

SAMPLE PREPARATION

Follow the instructions for use for the respective test serum used in combination with this product.

REAGENT PREPARATION

There is no preparation of the reagent required. Take and use the reagent directly from the vials.

TEST PROCEDURE

Follow the instructions for use for the respective test serum used in combination with this product.

INTERPRETATION OF RESULTS

Follow the instructions for use for the respective test serum used in combination with this product.

LIMITATIONS OF THE PROCEDURE


Follow the instructions for use for the respective test serum used in combination with this product.

LITERATURE

Follow the instructions for use for the respective test serum used in combination with this product.

GLOSSARY OF SYMBOLS

The following symbols may be used in the labelling of the device:

| | | | |
|---|--------------------------------------|---|-----------------------------|
|  REF | Product code |  LOT | Batch |
|  | Store from - to |  | Expiry Date |
|  IVD | In-Vitro Diagnostic |  | EC CE symbol |
|  | Manufacturer according to (EC) 98/79 |  | Consult instruction for use |
|  UDI | Unique Device Identification |  | Distributor |

REF

| | | |
|---------|-----------------|------------|
| 51310 | AHG serum blend | 10 ml |
| 51310.X | AHG serum blend | 10x 10 ml |
| 51310.L | AHG serum blend | 50 x 10 ml |

ANTITOXIN GmbH Industriestraße 88 69245 Bammental Germany

+49 (0) 6223/ 8661-0

+49 (0) 6223/ 8661-13

gara@antitoxin-gmbh.de

730-13-1709/ Version 012/ 2026-06-10

CE-IMMUNDIAGNOSTIKA GmbH
Karl-Landsteiner-Str. 6 D-69151 Neckargemuend