Anti-H Lectin (Laburnum alpinum)

Deutsch

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NUR ZUR IN-VITRO-DIAGNOSTIK

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FOR IN VITRO DIAGNOSTIC USE ONLY

INTENDED USE

Anti-H (A₂) reagent is obtained from seeds of the Laburnum alpinum. The lectin reacts specifically with H-substance. The reagent is used for qualitative in-vitro detection of the presence or absence of free H-substance on human red blood cells. It is used to differentiate between A₁ and A₂ blood. In A₁A₁ bleeding, the H substance is completely blocked, in A₁B, A₁O and especially BO types, small residues of H substance may remain and weak reactions may occur. Strong reactions occur with A₂ types or even weaker A types. As defined, strongest reaction occurs with O-bloods. The reagent is intended to be used by qualified and technical personnel only.

PRINCIPLE OF PROCEDURE

The test method used with this reagent is based on the principle of agglutination. Normal human erythrocytes, which have free H-substance, will agglutinate in the presence of the corresponding antibody.

REAGENTS

The listed reagent is offered in the following form:

Anti-H Lectin (Laburnum alpinum)

This reagent contains <0.1% (w/v) sodium azide as preservative. Beside the part of active lectin extract, the reagent contains sodium chloride, macromolecules and bovine albumin, which has been tested and certified by the US Veterinary service inspectors.

WARNING

This reagent is prepared from the vegetable extract of Laburnum alpinum. As biological product it should be looked upon as potentially infectious because of never complete exclusion of danger through excitants of disease. The reagent contains sodium azide that may be toxic and may react with lead or copper to form highly explosive salts. On disposal, flush with large quantities of water. For the reasons mentioned above the reagent should be handled with proper care.

STORAGE REQUIREMENT

Store opened and unopened products at +2 to +8°C. May be at room temperature while in use. In principle, store and use the reagent to declared expiry date only.

REMARKS

1. With each testing positive and negative controls should be performed.
2. Inappropriate storage impairs efficacy of the reagent
3. Weak turbidity of the reagent does not affect its reactivity. Bacteria and chemical contamination of the reagent should be avoided. If a visible change is detected, the reagent should no longer be used, this sign may indicate a microbiological contamination.
4. Strength of positive reactions also depends on age of used blood.
5. Centrifugation highly different from appointed relative centrifugal force may lead to false results.
6. The test method identified below are for manual testing only. When using automated or semi-automated instruments, follow the procedures that are contained in the operator's manual provided by the device manufacturer. Laboratories must follow approved validation procedures.
7. For usage of this reagent all effective national laws, directives and guidelines have to be observed in its current valid version, in Germany especially the „Richtlinien zur Gewinnung von Blut und Blutbestandteilen und zur Anwendung von Blutprodukten (Hämotherapie)“.

SAMPLE PREPARATION

1. Blood sample should be collected by approved medical procedure.
2. Blood samples to be tested should be used as soon as possible after blood collection to reduce the risk of false-positive and false-negative results due to improper storage or contamination of the reagent. If a delay in testing occurs, samples should be stored at +2 to +8°C. Blood drawn into EDTA should be tested within 7 days and samples treated with sodium citrate within 14 days after collection. Blood bag / Donor Blood can be tested until the expiry date.

REAGENT PREPARATION

There is no preparation of the reagent required. Take and use the reagent directly from the vials.

PROCEDURE

Not included but necessary materials:

- Tube Centrifugation Method:
1. test tubes (10 x 75 mm or 12 x 75 mm)
 2. microliter pipette
 3. centrifuge
 4. timer
 5. isotonic saline (0,85 - 0,9% sodium chloride)

Test procedure

Tube Centrifugation Method

1. Prepare 2% to 5% suspensions of red blood cells in isotonic saline (red blood cells may be washed 1-3 times with isotonic saline).
2. At first put 100 µl of the reagent in a marked tube, subsequently add 100 µl of appropriate cell suspension. Alternative one drop = 50 µl cell suspension can be added to one drop = 50 µl reagent.
3. Mix Erythrocytes- / Reagent mixture well by slightly shaking.
4. Incubate tube at room temperature for 15 min.
5. Centrifugation of the tube for 1 minute at 1.000 rpm (approximately 180-270 x g).
6. Gently shake the red cells completely from the bottom of the tube and check macroscopically for agglutination within 3 minutes.
7. Document the result.

INTERPRETATION OF RESULTS

"Slightly shaking" at Tube Centrifugation Method:


- Positive results (+): Visible agglutination of erythrocytes is a positive result and indicates the presence of the corresponding antigen.
- Negative results (-): No visible agglutination of erythrocytes is a negative result and indicates the absence of the corresponding antigen.

LIMITATIONS OF THE PROCEDURE

1. Inaccuracy at compliance with instructions written under section "Procedures" and "Interpretation of results" may lead to incorrect results.
2. No valid conclusion concerning the test result can be reached, if controls with uncertain or false results occur.
3. Enzyme treated erythrocytes or addition of bovine albumin and/or other solutions containing protein may cause unspecific reactions.
4. Hemolyzed, turbid, contaminated or clotted blood samples must not be used in this test.
5. No one specific antiserum or technique can be guaranteed to detect all rare, weak or variant of antigens.²
6. Due to the different expression of the free H substance on human red blood cells, the reactivity of the reagent against certain bloods may give weaker reaction than with control cells.
7. Red blood cells coated with alloantibodies or autoantibodies of the same or similar specificity as that reagent used for the test (i.e., cells that are positive in the direct antiglobulin-test (DAT)) are unsuitable.

LITERATURE

1. Richtlinien zur Gewinnung von Blut und Blutbestandteilen und zur Anwendung von Blutprodukten (Hämotherapie)
2. CLSI, I/LA33-A Validatin of Automated System for Immunhematological Testing Before Implementation; Approved Guidline CLINICAL AND LABORATORY STANDARDS INSTITUTE December 2009
3. Peter D. Issit, David J. Anstee Applied Blood Group Serology, 4. Edition, Montgomery Scientific Publications 1998
4. Christian Mueller-Eckhardt, Volker Kiefel Transfusionsmedizin, 3. edition, Springer-Verlag 2004.

 REF	Product Code	 LOT	Lot
	Store from - to		Expiration Date
 IVD	In-Vitro Diagnostic		EU CE symbol
	Manufacturer as to 98/79/EG		Observe instruction for use
 UDI	Unique Device Identification		Distributor

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ANTITOXIN GmbH Industriestraße 88 69245 Bammental Germany