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- 50 µL 0,8% Erythrozytensuspension in PBS Puffer oder

- 10 µL 3-5% Erythrozytensuspension in PBS Puffer

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## REF

- 02.182-02** Anti-Fy<sup>a</sup> for card method Coombs-reactive, polyclonal, human 2 ml  
**02.182-03** Anti-Fy<sup>a</sup> for card method Coombs-reactive, polyclonal, human 3 ml  
**02.182-05** Anti-Fy<sup>a</sup> for card method Coombs-reactive, polyclonal, human 5 ml  
**02.183-02** Anti-Fy<sup>b</sup> for card method Coombs-reactive, polyclonal, human 2 ml  
**02.183-03** Anti-Fy<sup>b</sup> for card method Coombs-reactive, polyclonal, human 3 ml  
**02.183-05** Anti-Fy<sup>b</sup> for card method Coombs-reactive, polyclonal, human 5 ml

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730-15-5106 Version 006 / 15. Januar 2021

ImuMed

# Anti-Fy<sup>a</sup> for card method Coombs-reactive, polyclonal, human Anti-Fy<sup>b</sup> for card method Coombs-reactive, polyclonal, human

English

FOR IN VITRO DIAGNOSTIC USE ONLY

## INTRODUCTION

Polyclonal Coombs-reactive Anti-Fy<sup>a</sup> and Anti-Fy<sup>b</sup> reagents are produced from human plasma that contains antibody of IgG-type, which reacts specifically with the corresponding antigen. The reagents are used to in-vitro determine qualitative whether red blood cells possess or lack the blood group antigens Fy<sup>a</sup> and Fy<sup>b</sup>. The reagents are intended to be used by qualified and technical personnel only.

## PRINCIPLE OF PROCEDURE

The test methods used with these reagents are based on the principle of column agglutination. Normal human erythrocytes, possessing one of these antigens, will be recognized and coated by the corresponding specific antibody and then the cells will be agglutinated by a secondary antibody that reacts with human IgG-molecules.

## REAGENTS

The listed reagents are available in a formulation, with differences in specificity for the antibody:

Anti-Fy<sup>a</sup> for card method Coombs-reactive, polyclonal, human

Anti-Fy<sup>b</sup> for card method Coombs-reactive, polyclonal, human

Both reagents contain <0.1% (w/v) sodium azide as preservative. Beside the parts active antibody and human serum the reagents contain sodium chloride, macromolecules and bovine albumin.

## WARNING

These reagents were prepared from human plasma. The raw material for these products are tested for HBsAg, HIV 1/2- and HCV-antibodies and found to be negative. Nevertheless as biological products it should be looked upon as potentially infectious because of never complete exclusion of danger through excitants of disease. The reagents contain sodium azide that may be toxic and may react with lead or copper to form highly explosive salts. On disposal, flush with large quantities of water. For the reasons mentioned above the reagent should be handled with proper care.

## STORAGE REQUIREMENT

Store opened and unopened products at +2 to +8°C. May be at room temperature (+15 to +30°C) while in use. In principle store and use the reagents only to declared expiry date.

## REMARKS

- Strength of positive reactions also depends on age of used blood.
- With each testing positive and negative controls should be executed.
- Weak turbidity of the reagents does not affect its reactivity.
- Inappropriate storage impairs efficacy of the reagents.
- The use of another card-specific centrifuge (each card centrifuge has its specified unchangeable g-force) may lead to false results due to the changed g-force.
- Blood samples to be tested should be used as soon as possible. If a delay in testing occurs, samples should be stored at +2 to +8°C. Blood drawn into EDTA should be tested within 7 days and samples treated with sodium citrate within 14 days after collection.
- The procedures identified below are for manual testing only and for the devices listed in this insert. When using other automated or semi-automated instruments, follow the procedures that are contained in the operator's manual provided by the device manufacturer. Laboratories must follow approved validation procedures to demonstrate compatibility of this product on such automated systems.
- For usage of these reagents all effective national laws, directives and guidelines have to be observed, in Germany especially the „Richtlinien zur Gewinnung von Blut und Blutbestandteilen und zur Anwendung von Blutprodukten (Hämotherapie)“.
- The information on the use of the test cards in the relevant insert must be observed.

## REAGENT PREPARATION

There is no preparation of the reagents required. Use reagents directly from the vials.

## PROCEDURE

Not provided material additionally needed:

### Card Method:

- Cards:
  - BIO-RAD (*DiaMed*) ID-Card "Liss/Coombs
  - Diagnostic Grifols „DG Gel Coombs“
  - Ortho BioVue® System AHG polyspecific cassettes or
  - Ortho BioVue® System Anti-Humanglobulin Anti-IgG cassettes
- Microliter pipette
- Glass tubes
- Isotonic saline (0,85 - 0,9% sodium chloride)
- Centrifuge
- Card Centrifuge
- Card incubator
- Card specific diluent

### Test procedure

**Card Method** (manual method / valid for the cards:

- BIO-RAD (*DiaMed*) ID-Card "Liss/Coombs
- Diagnostic Grifols „DG Gel Coombs“)

- Use 0,8 % suspension of red blood cells in card specific diluent (red blood cells may be washed 1 - 3 times with isotonic saline).
- Add 50 µL of appropriate cell suspension to each micro tube.
- Add 25 µL of appropriate reagent to each micro tube.
- Observe that the content of the reaction chamber(s) are combined. If necessary tap gently on side of card.
- Incubate card at +37°C for 15 minutes in an adequate incubator.
- Centrifuge the card in appropriate card centrifuge with the for this centrifuge unchangeable g-force.
- Check macroscopically for agglutination within 30 minutes.
- Document the result.

## Ortho BioVue® System (manual method)

- Use 3 - 5 % suspension of red blood cells in card specific diluent (red blood cells may be washed 1 - 3 times with isotonic saline).
- Add 40 µL of appropriate reagent to the appropriate reaction chamber(s)
- Add 10 µL of appropriate cell suspension to the appropriate reaction chamber(s).
- Check that the contents of the reaction chambers have been mixed. If necessary tap gently on the side of the card.
- Incubate card at +37°C for 15 minutes in an adequate incubator
- Centrifugation of cassette in appropriate centrifuge after latest 15 minutes with the, for the respective centrifuge, unchangeable g-force.
- The test results should be read directly after the end of centrifugation from front and back side of each micro tube.
- Document the result.

## Ortho BioVue® System (Ortho Vision® Analyzer)

The use of both reagents is validated on the Ortho Vision® Analyzer

The Ortho Vision® Analyzer offers the following options / settings:

- 50 µL 0,8% suspension of red blood cells in PBS buffer or
- 10 µL 3-5% suspension of red blood cells in PBS buffer

## INTERPRETATION OF RESULTS

Positive results (+): visible agglutination of erythrocytes is a positive result and indicates the presence of the corresponding antigen.

Negative results (-): No visible agglutination of erythrocytes is a negative result and indicates the absence of the corresponding antigen.

Carry out the reading and interpretation of the results according to the instruction for use of the corresponding card.

## LIMITATIONS OF THE PROCEDURE

- Inaccuracy at compliance with instructions written under section "Procedures" and "Interpretation of results" may lead to incorrect results.
- No valid conclusion concerning the test result can be reached, if controls with uncertain or false results occur.
- Enzyme treated erythrocytes or addition of bovine albumin and/or other solutions containing protein may cause unspecific reactions.
- Hemolyzed, turbid, contaminated or clotted blood samples must not be used in this test.
- Due to variability of antigen expression on human red blood cells, reactivity of the reagents, mentioned above, against certain phenotypes may give weaker reactivity compared to control cells.
- Red blood cells coated with alloantibodies or autoantibodies of the same or similar specificity as the reagent (i.e. cells that are positive in the direct antiglobulin-test (DAT)) are not suitable for this test procedure.
- Red blood cells with a positive direct Coombstest may cause false-positive reactions in card method.
- Information on limits in the instructions for use of cards used must be observed
- When used on the Ortho Vision® Analyzer, the test serum may remain uninterrupted on the device (+15 to +30 °C) for a maximum of 7 days. On the other hand, if the test serum is stored overnight at +2 to +8 °C for 12 hours, the period will be extended accordingly to 14 days.

## PERFORMANCE CHARACTERISTICS

A performance evaluation was conducted. Different samples (donor, patient, "panel blood") was used and compared with other products.

Product	Positive	Negative	False Positive	False Negative
Anti-Fy <sup>a</sup>	129	57	1	0
Anti-Fy <sup>b</sup>	142	43	2	0

The calculated values found for:

Product	Sensitivity	Specificity
Anti-Fy <sup>a</sup>	100%	98,28%
Anti-Fy <sup>b</sup>	100%	95,56%

## LITERATUR

- Peter D. Issitt, David J. Anstee Applied Blood Group Serology, fourth edition, Montgomery Scientific Publications 1998
- Christian Mueller-Eckhardt, Volker Kiefel Transfusionsmedizin, 3. Auflage, Springer-Verlag 2004

<b>REF</b> Product Code	<b>LOT</b> Batch
 Store from - to	 Expiration Date
<b>IVD</b> In-Vitro Diagnostic	<b>CE</b> EU CE symbol
 Manufacturer according to 98/79/EG	 Consult instructions for use

## REF

- |                  |   |      |
|------------------|---|------|
| <b>02.182-02</b> | Anti-Fy <sup>a</sup> for card method Coombs-reactive, polyclonal, human | 2 ml |
| <b>02.182-03</b> | Anti-Fy <sup>a</sup> for card method Coombs-reactive, polyclonal, human | 3 ml |
| <b>02.182-05</b> | Anti-Fy <sup>a</sup> for card method Coombs-reactive, polyclonal, human | 5 ml |
| <b>02.183-02</b> | Anti-Fy <sup>b</sup> for card method Coombs-reactive, polyclonal, human | 2 ml |
| <b>02.183-03</b> | Anti-Fy <sup>b</sup> for card method Coombs-reactive, polyclonal, human | 3 ml |
| <b>02.183-05</b> | Anti-Fy <sup>b</sup> for card method Coombs-reactive, polyclonal, human | 5 ml |

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**REF**

- 02.182-02** Anti-Fy<sup>a</sup> for card method Coombs-reactive, polyclonal, human 2 ml
- 02.182-03** Anti-Fy<sup>a</sup> for card method Coombs-reactive, polyclonal, human 3 ml
- 02.182-05** Anti-Fy<sup>a</sup> for card method Coombs-reactive, polyclonal, human 5 ml
- 02.183-02** Anti-Fy<sup>b</sup> for card method Coombs-reactive, polyclonal, human 2 ml
- 02.183-03** Anti-Fy<sup>b</sup> for card method Coombs-reactive, polyclonal, human 3 ml
- 02.183-05** Anti-Fy<sup>b</sup> for card method Coombs-reactive, polyclonal, human 5 ml

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- REF**
- 02.182-02 Anti-Fy<sup>a</sup> for card method Coombs-reactive, polyclonal, human 2 ml
  - 02.182-03 Anti-Fy<sup>a</sup> for card method Coombs-reactive, polyclonal, human 3 ml
  - 02.182-05 Anti-Fy<sup>b</sup> for card method Coombs-reactive, polyclonal, human 5 ml
  - 02.183-02 Anti-Fy<sup>b</sup> for card method Coombs-reactive, polyclonal, human 2 ml
  - 02.183-03 Anti-Fy<sup>b</sup> for card method Coombs-reactive, polyclonal, human 3 ml
  - 02.183-05 Anti-Fy<sup>b</sup> for card method Coombs-reactive, polyclonal, human 5 ml

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## REF

- 02.182-02 Anti-Fy<sup>a</sup> for card method Coombs-reactive, polyclonal, human 2 ml
- 02.182-03 Anti-Fy<sup>a</sup> for card method Coombs-reactive, polyclonal, human 3 ml
- 02.182-05 Anti-Fy<sup>a</sup> for card method Coombs-reactive, polyclonal, human 5 ml
- 02.183-02 Anti-Fy<sup>b</sup> for card method Coombs-reactive, polyclonal, human 2 ml
- 02.183-03 Anti-Fy<sup>b</sup> for card method Coombs-reactive, polyclonal, human 3 ml
- 02.183-05 Anti-Fy<sup>b</sup> for card method Coombs-reactive, polyclonal, human 5 ml

CE 0483

ANTITOXIN GmbH Industriestrasse 88 69245 Bammental Alemanha

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ImuMed

**Anti-Fy<sup>a</sup> for card method Coombs-reactive, polyclonal, human**  
**Anti-Fy<sup>b</sup> for card method Coombs-reactive, polyclonal, human**  
**PARA USO EXCLUSIVO EN DIAGNÓSTICO IN VITRO**

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### REF

02.182-02	Anti-Fy <sup>a</sup> for card method Coombs-reactive, polyclonal, human	2 ml
02.182-03	Anti-Fy <sup>a</sup> for card method Coombs-reactive, polyclonal, human	3 ml
02.182-05	Anti-Fy <sup>a</sup> for card method Coombs-reactive, polyclonal, human	5 ml
02.183-02	Anti-Fy <sup>b</sup> for card method Coombs-reactive, polyclonal, human	2 ml
02.183-03	Anti-Fy <sup>b</sup> for card method Coombs-reactive, polyclonal, human	3 ml
02.183-05	Anti-Fy <sup>b</sup> for card method Coombs-reactive, polyclonal, human	5 ml

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#### REF

<b>02.182-02</b>	Anti-Fy <sup>a</sup> for card method Coombs-reactive, polyclonal, human	2 ml
<b>02.182-03</b>	Anti-Fy <sup>a</sup> for card method Coombs-reactive, polyclonal, human	3 ml
<b>02.182-05</b>	Anti-Fy <sup>a</sup> for card method Coombs-reactive, polyclonal, human	5 ml
<b>02.183-02</b>	Anti-Fy <sup>b</sup> for card method Coombs-reactive, polyclonal, human	2 ml
<b>02.183-03</b>	Anti-Fy <sup>b</sup> for card method Coombs-reactive, polyclonal, human	3 ml
<b>02.183-05</b>	Anti-Fy <sup>b</sup> for card method Coombs-reactive, polyclonal, human	5 ml

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ANTITOXIN GmbH Industriestrasse 88 69245 Bammental Almanya