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REF

BG1431 Anti-Le^a monoclonal, mouse clone: LEA2 2 ml

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Anti-Le^a monoclonal

For Tube Test
FOR IN VITRO DIAGNOSTIC USE ONLY

English

INTENDED USE

The reagent is used to in-vitro determine qualitative whether human red blood cells possess or lack the corresponding blood group antigen Lea.
The reagent is intended to be used by qualified and technical personnel only to perform immunohematology screening tests as part of the practice of transfusion medicine in the general population.
The test method used with this reagent is based on the principle of agglutination, performed not automated.
Normal human erythrocytes, possessing the corresponding antigen, will agglutinate in the presence of the specific antibody directed toward the antigen.

INDICATION / CONTRA-INDICATION

The Anti-Lea monoclonal blood grouping reagent is used to test patient or donor red cells for the presence of the Lea antigen. Typing of donor cells facilitates the selection of suitable antigen-negative units for transfusion to patients with this antibody. Cell typing also serves as final verification of the identification of Anti-Lea in patient or donor sera.
There is no contra-indication to perform the in-vitro-test on blood samples.
The product was validated with sample collected in Europe from patients of unknown ethnic background.

REAGENTS

The listed reagents contains antibodies of the following clones:
Anti-Le^a monoclonal, mouse clone: LEA2

The reagents contains <0.1% (w/v) sodium azide as preservative. Beside the part active antibody, the reagents contains sodium chloride, macromolecules and bovine albumin, which has been tested and certified by the US Veterinary service inspectors.

WARNING

The reagents are prepared from supernatants of cell cultures. Nevertheless as biological products it should be looked upon as potentially infectious because of never complete exclusion of danger through excitants of disease. The reagents contains sodium azide that may be toxic and may react with lead or copper to form highly explosive salts.
On disposal, flush with large quantities of water.
For the reasons mentioned above, reagents should be handled with proper care.

STORAGE REQUIREMENT

Store opened and unopened products at +2 to +8°C.
May be at room temperature while in use.
In principle, store and use the reagents to indicated expiry date only.

REMARKS

1. With each testing positive and negative controls should be performed.
2. Inappropriate storage impairs efficacy of the reagents.
3. Weak turbidity of the reagents does not affect its reactivity. Bacteria and chemical contamination of the product should be avoided. If a visible change is detected, it may indicate a microbiological contamination, the reagent should no longer be used.
4. Strength of positive reactions also depends on age of used blood
5. Centrifugation outside the specified speed range may lead to false results.
6. The test method identified below is for manual testing only. When using automated or semi-automated instruments, follow the procedures that are contained in the operator's manual provided by the device manufacturer. Laboratories must follow approved validation procedures.
7. For usage of the reagents all effective national laws, directives and guidelines have to be observed in its current valid version, in Germany especially the „Richtlinien zur Gewinnung von Blut und Blutbestandteilen und zur Anwendung von Blutprodukten (Hämotherapie)“¹.

SAMPLE PREPARATION

1. Blood sample should be collected by approved medical procedure.
2. Blood sample to be tested should be used as soon as possible after blood collection to reduce the risk of false-positive and false-negative results due to improper storage or contamination of the sample.
If a delay in testing occurs, sample should be stored at +2 to +8°C.
Blood drawn into EDTA should be tested within 7 days and sample treated with sodium citrate within 14 days after collection.
Blood bag / Donor Blood can be tested until the expiry date.

REAGENT PREPARATION

There is no preparation of the reagents required.
Take and use reagents directly from the vials.

PROCEDURE

Not provided material, additionally needed

- Tube Centrifugation Method:
1. tubes (10 x 75 mm or 12 x 75 mm)
 2. microliter pipette
 3. timer
 4. centrifuge
 5. isotonic saline (0,85 - 0,9% sodium chloride)

Test procedure

Tube Centrifugation Method

1. Prepare 2% to 5% suspensions of red blood cells in isotonic saline (red blood cells may be washed 1-3 times with isotonic saline).
2. At first put 100 µL of appropriate reagent in a marked tube, subsequently add 100 µL of appropriate cell suspension.
Alternative one drop = approximately 50 µL cell suspension can be added to one drop = approximately 50 µL reagent.
3. Mix Erythrocytes-/Reagent-mixture well by slightly shaking.
4. Incubate tube at room temperature for 15 min.
5. Centrifuge tube for 1 min at 1.000 rpm (approximately 180-270 x g).
6. Gently shake the red cells from the bottom of the tube and check macroscopically for agglutination within 3 minutes.
7. Document the result.

INTERPRETATION OF RESULTS

"Slightly shaking" at Tube Centrifugation Method.

- Positive results (+): visible agglutination of erythrocytes is a positive result and indicates the presence of the corresponding antigen.
Negative results (-): No visible agglutination of erythrocytes is a negative result and indicates the absence of the corresponding antigen.

LIMITATIONS OF THE PROCEDURE

1. Inaccuracy at compliance with instructions written under section "Procedures" and "Interpretation of results" may lead to incorrect results.
2. No valid conclusion concerning the test result can be reached, if controls with uncertain or false results occur.
3. Enzyme treated erythrocytes or addition of bovine albumin and/or other solutions containing protein may cause unspecific reactions.
4. Hemolyzed, turbid, contaminated or clotted samples should not be used.
5. Due to variability of antigen expression, reactivity of the reagents against certain phenotypes may give weaker reactivity compared to control cells.
6. No one specific antiserum or technique can be guaranteed to detect all rare, weak or variant Antigens.²
7. Red blood cells coated with alloantibodies or autoantibodies of the same or similar specificity as that reagent used for the test (i.e., cells that are positive in the direct antiglobulin-test (DAT)) are not suitable for this test procedure.
8. There may be weak positive reactions with Le^a or Le^b negative cells. This phenomenon occurs especially with commercial test cells and may be reduced by washing the cells.

INCIDENTS RELATED TO THE DEVICE

Any serious incident that has occurred in relation to the device must be reported to the manufacturer and the competent authority of the Member State where the user and/or patient is established.

PERFORMANCE CHARACTERISTICS

A performance evaluation for the products was carried out.
The required samples were used and compared with other reference methods / products.

Technique	Product	Performance Characteristics			
		Positive Bleeds n	Sensitivity	Negative Bleeds n	Specificity
Tube methode		17/17	100 %	63/63	100 %

Diagnostic Sensitivity: The probability that the device gives a positive result in the presence of the target marker.

Diagnostic Specificity: The probability that the device gives a negative result in the absence of the target marker.

The product is equivalent and does not differ in quality from comparable reagents available on the market.

DIFFERENCES BETWEEN BATCHES

Validation between three batches over the entire shelf life showed no differences.

INTERFERENCE STUDY

The interference studies showed no impairment for the qualitative test when using the following interfering substances:

Heparin 720 U/dl, Albumin 15000 mg/dl, Triglycerides 1500 mg/dl, Bilirubin 40mg/dl, Ethanol 620 mg/dl, Glucose 1000 mg/dl.

For the anticoagulants (EDTA, Sodium Citrate, ACD, CPD-A, PAGGS-M) three times the recommended concentration was tested.

SUMMARY OF SAFETY AND PERFORMANCE

The Summary of Safety and Performance of this reagent is available via ANTITOXIN (www.antitoxin-gmbh.de) and can be accessed via the EUDAMED database.

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REF	Product Code	LOT	Lot
	Store from - to		Expiration Date
IVD	In-Vitro Diagnostic	EU CE symbol	
	Manufacturer as to 98/79/EG		Observe instruction for use
UDI	Unique Device Identification		Distributor

REF
BG1431 Anti-Le^a monoclonal, mouse clone: LEA2 2 ml

Version R002 / 16.08.2024

CE 0483

 ANTITOXIN GmbH Industriestraße 88 69245 Bammental Germany
 +49 (0) 6223/ 8661-0  +49 (0) 6223/ 8661-13  qara@antitoxin-gmbh.de

sifin

sifin diagnostics gmbh
 Berliner Allee 317-321
 13088 Berlin / Germany
 +49 30 700 144-0