Bromelin control solution

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Bromelin control solution

For Tube Method
FOR IN VITRO DIAGNOSTIC USE ONLY

English

INTENDED USE

Bromelin control solution serves as an indicator of whether in a bromelin solution test approach with a negative test result, the added Bromelin solution has optimal activity and has been added in the correct amount. The Bromelin control solution is a mixture of various incomplete antibodies of the Rhesus-system. This mixture leads to agglutination in the presence of Bromelin solution with all human red cells.

The reagent is intended to be used by qualified and technical personnel only.

PRINCIPLE OF PROCEDURE

The test method use with this reagent is based on the principle of agglutination. Normal human erythrocytes, are agglutinated by the presence of antibodies in the presence of active Bromelin solution. If the enzyme solution is missing in the reaction approach, the amount is not sufficient or no longer active, no agglutination takes place.

REAGENTS

The reagent is available in following formulation:

Bromelin control solution incomplete, polyclonal, human

The reagent contains <0.1% (w/v) sodium azide as preservative. Beside the parts active antibodies and human serum, the reagent contains sodium chloride, macromolecules and bovine albumin, which has been tested and certified by the US Veterinary service inspectors.

WARNING

The reagent is prepared from human plasma.

The raw materials for this reagent have been tested for HBsAg, HIV1/2- and HCV-antibodies and found to be negative. Nevertheless, as biological product it should be looked upon as potentially infectious because of never complete exclusion of danger through excitants of disease. The reagent contains sodium azide that may be toxic and may react with lead or copper to form highly explosive salts.

On disposal, flush with large quantities of water.

For the reasons mentioned above, the reagent should be handled with proper care.

STORAGE REQUIREMENT

Store opened and unopened products at +2 to +8°C.

May be at room temperature while in use.

In principle, store and use the reagent to indicated expiry date only.

REMARKS

1. With each testing positive and negative controls should be performed.
2. Inappropriate storage impairs efficacy of the reagent.
3. A weak turbidity of the reagent does not affect its reactivity. Bacteria and chemical contamination of the product should be avoided. If a visible change is detected, it may indicate a microbiological contamination, the reagent should no longer be used.
4. Strength of positive reactions also depends on age of used blood
5. Centrifugation outside the specified speed range may lead to false results.
6. The test method identified below is for manual testing only. When using automated or semi-automated instruments, follow the procedures that are contained in the operator's manual provided by the device manufacturer. Laboratories must follow approved validation procedures.
7. For usage of this reagent all effective national laws, directives and guidelines have to be observed in its current valid version, in Germany especially the „Richtlinien zur Gewinnung von Blut und Blutbestandteilen und zur Anwendung von Blutprodukten (Hämotherapie)“.

SAMPLE PREPARATION

1. Blood sample should be collected by approved medical procedure.
2. Blood sample to be tested should be used as soon as possible after blood collection to reduce the risk of false-positive and false-negative results due to improper storage or contamination of the sample.
If a delay in testing occurs, sample should be stored at +2 to +8°C.
Blood drawn into EDTA should be tested within 7 days and sample treated with sodium citrate within 14 days after collection.
Blood bag / Donor Blood can be tested until the expiry date.

REAGENT PREPARATION

There is no preparation of the reagent required.
Take and use reagent directly from the vials.

PROCEDURE

Not provided material, additionally needed

Tube Centrifugation Method: 1. tubes (10 x 75 mm or 12 x 75 mm)
2. microliter pipette
3. timer
4. centrifuge
5. isotonic saline (0,85 - 0,9% sodium chloride)

Test procedure

Tube Centrifugation Method

Bromelin Test

1. a) Prepare 2% to 5% suspensions of red blood cells in isotonic saline (red blood cells may be washed 1-3 times with isotonic saline) or
b) use testerythrocytes of a commercially acquired cell panel.
2. In a labeled test tube
a) give 100µl of an equivalent test serum, or
b) give 100µl serum/plasma of a patient/proband to be tested
3. add to the test tube 100µl
a) of the appropriate cell suspension or
b) testerythrocytes of a commercially acquired cell panel
4. Then add one drop of Bromelin solution
5. Mix Erythrocytes / Reagentmixture well by slightly shaking.
6. Incubate tube at room temperature for 5-15 min.
7. Centrifugation of tube for 1 minute at 1.000 rpm (approximately 180-270 x g).
8. Gently shake the red cells completely from the bottom of the tube and check macroscopically for agglutination within 3 minutes.
9. Document the result.

Bromelin Control Test

1. Then add 100 µl Bromelin control solution to all test tubes with a negative test result.
2. Mix Erythrocytes / Reagentmixture well by slightly shaking.
3. Centrifugation of tube for 1 minute at 1.000 rpm (approximately 180-270 x g).
4. Gently shake the red cells completely from the bottom of the tube and check macroscopically for agglutination within 3 minutes.
5. Document the result. The test tube must now show agglutinations.

INTERPRETATION OF RESULTS

"Slightly shake up" of the erythrocyte sediment.

Bromelin Test:

Positive results (+): Visible agglutination of erythrocytes is a positive result and indicates
a) the presence of the corresponding antigen
b) the presence of one or several antibodies which can be identified with the help of the antigam of the cell panel.

Negative results (-): No visible agglutination of erythrocytes is a negative result and indicates
a) the absence of the corresponding antigen
b) an irregular antibody does not exist.

Bromelin Control Test:


A **negative** result after the addition of Bromelin control solution indicates, that the Bromelin solution has not been added in the correct amount or it has no longer optimal activity and a completely new approach of the test series with fresh Bromelin solution required!

LIMITATIONS OF THE PROCEDURE

1. Inaccuracy at compliance with instructions written under section "Procedures" and "Interpretation of results" may lead to incorrect results.
2. No valid conclusion concerning the test result can be reached, if controls with uncertain or false results occur.
3. Hemolyzed, turbid, contaminated or clotted samples should not be used.
4. Due to variability of antigen expression, reactivity of this reagent against certain phenotypes may give weaker reactivity compared to control cells.
5. Red blood cells with a positive direct Coombs-test may cause false-positive reactions. The cells react positively even without a test serum.

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 REF	Product Code	 LOT	Lot
	Store from - to		Expiration Date
 IVD	In-Vitro Diagnostic		EU CE symbol
	Manufacturer as to 98/79/EG		Observe instruction for use
 UDI	Unique Device Identification		Distributor

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ANTITOXIN GmbH Industriestraße 88 69245 Bammental Germany