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REF

- 01.198-03 Anti-H for card method monoclonal, mouse clone: 10934C11 3 ml
- 01.198-05 Anti-H for card method monoclonal, mouse clone: 10934C11 5 ml

FOR IN VITRO DIAGNOSTIC USE ONLY

INTENDED USE

Monoclonal agglutinating Anti-H (A₂) reagent is produced from cell culture supernatants of hybridoma-cell line. The cells are secreting an antibody of IgM-type. The reagent is used to in-vitro determine qualitative whether red blood cells possess or lack free H-substance on human erythrocytes. It is used to differentiate between A₁ and A₂ blood. There is no free H-substance on A₁ A₁-blood, but with A₁B-, A₁O- and especially with BO-types there may still be free rests of H-substance and show weak reactions. Strong reactions are visible with A₂-blood or weaker A-variants. As defined strongest reaction are seen with O-bloods. The reagent is intended to be used by qualified and technical personnel only.

PRINCIPLE OF PROCEDURE

The test methods used with this reagent is based on the principle of column agglutination. Normal human erythrocytes which have free H-substance are agglutinated by the corresponding antibody.

REAGENTS

The listed reagent contains antibodies of the following cell clone:

Anti-H for card method monoclonal, mouse clone: 10934C11

The reagent contains <0.1% (w/v) sodium azide as preservative. Additionally the reagent is prepared of active antibody, sodium chloride, macromolecules and bovine albumin, which has been tested and certified by the US Veterinary service inspectors.

WARNING

This reagent is prepared from supernatant of cell cultures. As biological product it should be looked upon as potentially infectious because of never complete exclusion of danger through excitants of disease. The reagent contains sodium azide that may be toxic and may react with lead or copper to form highly explosive salts. On disposal, flush with large quantities of water. For the reasons mentioned above the reagent should be handled with proper care.

STORAGE REQUIREMENT

Store opened and unopened products at +2 to +8°C. May be at room temperature while in use. In principle, store and use the reagents to declared expiry date only.

REMARKS

1. With each testing positive and negative controls should be performed.
2. Inappropriate storage impairs efficacy of the reagent.
3. Weak turbidity of the reagent does not affect its reactivity. Bacteria and chemical contamination of the product should be avoided. If a visible change is detected, the reagent should no longer be used, this sign may indicate a microbiological contamination.
4. Strength of positive reactions also depends on age of used blood.
5. Use the appropriate card centrifuge. The use of another card-specific centrifuge (each card centrifuge has its specified unchangeable g-force) may lead to false results due to the changed g-force.
6. The test methods identified below are for manual testing only and for the devices listed in this insert. When using other automated or semi-automated instruments, follow the procedures that are contained in the operator's manual provided by the device manufacturer. Laboratories must follow approved validation procedures.
7. For usage of the reagent all effective national laws, directives and guidelines have to be observed, in Germany especially the „Richtlinien zur Gewinnung von Blut und Blutbestandteilen und zur Anwendung von Blutprodukten (Hämotherapie)“¹ in its actual form.
8. The information on the use of the test cards in the relevant insert must be observed.

SAMPLE PREPARATION

1. Blood sample should be collected by approved medical procedure.
2. Blood samples to be tested should be used as soon as possible after blood collection to reduce the risk of false-positive and false-negative results due to improper storage or contamination of the samples. If a delay in testing occurs, samples should be stored at +2 to +8°C. Blood drawn into EDTA should be tested within 7 days and samples treated with sodium citrate within 14 days after collection. Blood bag / Donor Blood can be tested until the expiry date.

REAGENT PREPARATION

There is no preparation of the reagent required. Take and use the reagent directly from the vials.

PROCEDURE

Not included but necessary materials:

Card Method manuell:

1. Cards: - Grifols „DG Gel Neutral“
- Ortho BioVue® System Reverse
2. Microliter pipette
3. Tubes
4. Isotonic saline (0,85 - 0,9% sodium chloride)
5. Centrifuge
6. Card Centrifuge
7. Card specific diluent

Card Method automatic:

Vision® Analyzer

Test procedure

Card Method (manual method / valid for the cards:

- Grifols „DG Gel Neutral“)

1. Prepare 0,8 % suspension of red blood cells in card specific diluent (red blood cells may be washed 1 - 3 times with isotonic saline).
2. Add 50 µL of appropriate cell suspension to each marked micro tube.
3. Add 25 µL of the reagent to each micro tube.
4. Centrifuge the card in appropriate card centrifuge with the for this centrifuge unchangeable g-force.
5. Check macroscopically for agglutination within 30 minutes.
6. Document the result.

Ortho BioVue® System (manual method)

1. Prepare 3-5 % suspensions of red blood cells in card specific diluent (red blood cells may be washed 1 - 3 times with isotonic saline).
2. Add 40 µL of the reagent to the appropriate marked reaction chambers.
3. Add 10 µL of appropriate cell suspension to the appropriate reaction chambers.
4. Centrifuge the cassette in appropriate centrifuge with the for this centrifuge unchangeable g-force within 15 minutes.
5. The test results should be read directly after the end of centrifugation.
6. Document the result.

Ortho BioVue® System (Ortho Vision® Analyzer)

The use of the mentioned above reagent is validated on the Ortho Vision® Analyzer

The Ortho Vision® Analyzer offers the following options / settings:

- 50 µL 0,8% suspension of red blood cells in PBS buffer or
- 10 µL 3-5% suspension of red blood cells in PBS buffer

INTERPRETATION OF RESULTS

Positive results (+): visible agglutination of erythrocytes is a positive result and indicates the presence of the corresponding antigen.

Negative results (-): No visible agglutination of erythrocytes is a negative result and indicates the absence of the corresponding antigen.

Read and interpret the results of the card method according to the instruction for use of the corresponding card.

The interpretation of the results from Ortho Vision® Analyzer is performed automatically by the software.

LIMITATIONS OF THE PROCEDURE

1. Inaccuracy at compliance with instructions written under section "Procedures" and "Interpretation of results" may lead to incorrect results.
2. No valid conclusion concerning the test result can be reached, if controls with uncertain or false results occur.
3. Enzyme treated erythrocytes or addition of bovine albumin and/or other solutions containing protein may cause unspecific reactions.
4. Hemolyzed, turbid, contaminated or clotted blood samples must not be used in this test.
5. Due to variability of free H-substance on human red blood cells, reactivity of the reagent, mentioned above, against certain phenotypes may give weaker reactivity compared to control cells.
6. No one specific antiserum or technique can be guaranteed to detect all rare, weak or variant of "antigens".²
7. Red blood cells with a positive direct Coombs-test may cause false-positive reactions in card method.
8. Information on limits in the instructions for use of cards used must be observed.
9. When used on the Ortho Vision® Analyzer, the test serum may remain uninterruptedly on the device (+15 to +30 °C) for a maximum of 7 days. On the other hand, if the test serum is stored overnight at +2 to +8 °C for 12 hours, the period will be extended accordingly to 14 days.

LITERATURE

1. Richtlinien zur Gewinnung von Blut und Blutbestandteilen und zur Anwendung von Blutprodukten (Hämotherapie)
2. CLSI, I/LA33-A Validatin of Automated System for Immunhematological Testing Before Implementation: Approved Guideline CLINICAL AND LABORATORY STANDARDS INSTITUTE December 2009
3. Peter D. Issit, David J. Anstee Applied Blood Group Serology, 4. edition, Montgomery Scientific Publications 1998
4. Christian Mueller-Eckhardt, Volker Kiefel Transfusionsmedizin, 3.edition, Springer-Verlag 2004.

 Product Code	 Batch
 Store from - to	 Expiration Date
 In-Vitro Diagnostic	 EU CE symbol
 Manufacturer according to 98/79/EG	 Consult instructions for use

REF

01.198-03 Anti-H for card method monoclonal, mouse clone: 10934C11 3 ml

01.198-05 Anti-H for card method monoclonal, mouse clone: 10934C11 5 ml

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ANTITOXIN GmbH Industriestrasse 88 69245 Bammental Nemčija