ABO-trol[®] forte plus

NUR ZUR IN-VITRO-DIAGNOSTIK

CE 0483

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REF

01.901-30 d-d-e ABO-trol[®] forte plus ready to use for verification of ABO-compatibility and D confirmation prior to blood transfusion

01.901-30-2 ABO-trol[®] forte plus 15x4 Teststreifen / Test Strips ready to use for verification of ABO-compatibility and D confirmation prior to blood transfusion

ABO-trol[®] forte plus

FOR IN-VITRO-DIAGNOSTIC USE ONLY CE 0483

Purpose

ABO-trol[®] forte plus is a ready for use bedside-test card for the verification of ABO- and D-compatibility prior to blood transfusion. Predropped fluid antisera, with approximately 70 µl of each, anti-A (mouse), anti-B (mouse) and anti-D (human), sufficient for one determination, are produced from cell culture supernatants of hybridoma cell lines that secrete antibodies of the IgM type.

The test cards are intended to be used by qualified and technical personnel only.

Prepared test sera, per card:

ABO-trol[®] forte plus

2x ca. 70 µl Anti-A monoclonal, IgM clone: Birma-1

2x ca. 70 µl Anti-B monoclonal, IgM clone: LB-2

2x ca. 70 µl Anti-D monoclonal, human IgM clone: RUM-1

One pack contains 30 cards.

ABO-trol[®] forte plus

15x4 Teststreifen / Test Strips

4x ca. 70 µl Anti-A monoclonal, IgM clone: Birma-1

4x ca. 70 µl Anti-B monoclonal, IgM clone: LB-2

4x ca. 70 µl Anti-D monoclonal, human IgM clone: RUM-1

One pack contains 15 cards.

Principle

The test sera are used to confirm the presence or absence of the blood grouping antigens A and B, respectively on human erythrocytes. They are characterized by quick and strong reactions with possibly present and corresponding antigens through an antigen-antibody reaction. The used test serum anti-D reacts with strong D antigens only within 1 minute in this technique.

Storage and Stability

ABO-trol[®] forte plus is to be stored at +2 to +8 °C.

Under adequate storage the antisera keep their reactivity up to shelf life printed on the card.

Warning

These reagents were not tested for murine viruses.

Remarks

As preservative the test sera in ABO-trol[®] forte plus contain sodium azide <1 g/l that may have a toxic effect and form explosive salts with lead or copper.

Basics

„Right before a transfusion of erythrocyte concentrates, native blood, granulocyte concentrates and plasma exchange therapy, the transfusing physician himself or, under his direct supervision, the ABO-identity test (bedside-test) should be performed directly on the recipient (e.g. on test cards). The test material is disposed of. It serves to confirm the previously determined ABO-blood group characteristics of the recipient. The result must be documented in writing. In the case of discrepancies, the laboratory or the transfusion medical institution must be informed immediately.“ (Source: *“Guidelines for the Preparation of Blood and Blood Components and for the Use of Blood Products (Hemotherapy)”*).

Product

ABO-trol[®] forte plus consists of a plastic card with 6 trough-shaped cavities

ABO-trol[®] forte plus 15x4 Teststreifen / Test Strips consists of a plastic card with 12 trough-shaped cavities serving as reaction chambers.

Each of these chambers is filled with approximately 70 µl of antiserum, blue coloured anti-A on the left, yellow coloured anti-B in the middle and uncoloured Anti-D on the right.

The reaction chambers filled with antisera are sealed with a transparent film. The chambers are intended for confirmation testing of patient erythrocytes and donor erythrocytes relating to ABO compatibility and D-antigen.

Test Implementation

For identity testing, erythrocyte concentrate, native blood and blood with anticoagulatives (citrate, EDTA) can be used.

- Note data of patient and/or donor on adhesive label.
- Collect blood from the patient.
- Penetrate the transparent film of chambers „anti-A“ with the cannula that is attached to the syringe and add 1 single small drop of patient erythrocytes, free-falling, into the chamber.
- With chamber “Anti-B” and “Anti-D”, please proceed as described under 3.
- Please treat donor blood as under 3. and 4. with empty chambers.
- Slightly shake the bedside testcard up and down, mix.
- After 1 minute, check on agglutination (= positive reaction; reactions with Anti-A and Anti-B visible after 10 seconds).
The maximum agglutination is shown after 3 minutes.
- Note result, date and signature on the adhesive label. Remove the label and place it into the medical file.

Evaluation

The bedside test is an ABO identity test and not a complete blood group determination. The following chart shows the compatibility of uneven, low in plasma erythrocyte concentrates.

Anti-A	Anti-B	Confirmation blood group	Compatible donor erythrocyte concentrate
pos	neg	A	A or O
neg	pos	B	B or O
pos	pos	AB	AB, A, B or O
neg	neg	O	O

If the blood group of the tested patient's erythrocytes is identical with the blood group marked on the blood bag for the transfusion, the ABO compatibility is verified. In case the agglutination of the patient's erythrocytes on Anti-A or Anti-B is not identical with what is marked on the blood bag for the transfusion, the ABO compatibility is not verified. In this case, a clarification is required.

For the result with the anti-D this applies only limited, since only strong D antigens are indicated. With a different agglutination picture, special attention must therefore be paid to the already existing laboratory results!

Performance Data

The ABO-trol[®] forte plus bedside cards are based on the antisera and were tested on their specificity/sensitivity with blood samples.

The test results are shown in the following chart.

Test serum	positive	negative	specificity	sensitivity
Anti-A	864	798	100%	100%
Anti-B	344	1318	100%	100%
Anti-D	574	190	100%	100%

Remarks and Limitations of the Procedure

- False positive reactions may occur by:
 - bacterial contamination of serum and/or of the blood samples
 - Wharton's Jelly
 - antibodies diverted against colour ingredients or preservatives
 - medicinal therapy or pathological discrepancies, dextrans or other plasma expanders
 - erythrocytes strongly loaded with antibodies (e.g. autoantibodies)
- False negative results may occur by:
 - drying up of test serum in reaction chambers (e.g. storage temperature too high)
 - wrong concentration of erythrocytes (dilution of blood sample)
 - chambers overfilled with blood (more than 1 single small drop)
- The test serum Anti-B does not react with “acquired B”.
- Newborn blood is excluded from the test.
- False positive or false negative results lead to discrepancies between bedside test and previous test results. In this case ABO-compatibility is not confirmed and further test have to be carried out with different test systems.

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REF Product Code	LOT Lot
Store from - to	Expiration Date
IVD In-Vitro Diagnostic	CE EU CE symbol
Manufacturer as to 98/79/EG	Observe instructions for use
UDI Unique Device Identification	

REF

01.901-30 d-d-e ABO-trol[®] forte plus ready to use for verification of ABO-compatibility and D confirmation prior to blood transfusion

01.901-30-2 ABO-trol[®] forte plus 15x4 Teststreifen / Test Strips ready to use for verification of ABO-compatibility and D confirmation prior to blood transfusion

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